



Technical Report
RAL-TR-2001-029



The Application of the NeXus Data Format to ISIS Muon Data

D Flannery S P Cottrell and P J C King

7th August 2001



© Council for the Central Laboratory of the Research Councils 2001

Enquiries about copyright, reproduction and requests for additional copies of this report should be addressed to:

The Central Laboratory of the Research Councils
Library and Information Services
Rutherford Appleton Laboratory
Chilton
Didcot
Oxfordshire
OX11 0QX
Tel: 01235 445384 Fax: 01235 446403
E-mail library@rl.ac.uk

ISSN 1358-6254

Neither the Council nor the Laboratory accept any responsibility for loss or damage arising from the use of information contained in any of their reports or in any communication about their tests or investigations.

The Application of the NeXus Data Format to ISIS Muon Data

D. Flannery, S.P. Cottrell, P.J.C. King

ISIS Muon Facility
Rutherford Appleton Laboratory

Contents

Contents	3
1. INTRODUCTION	7
1.1 Experiment Data Acquisition	8
1.2 Why choose a different format?	8
1.3 Why use NeXus?	9
2. ADAPTING THE NeXus FORMAT FOR μ SR DATA	11
2.1 ISIS Muon NeXus Instrument Definition	11
2.2 Abstract Model of Instrument Definition	14
2.3 Data Storage using the ISIS Muon NeXus Instrument Defintion	16
3. CONVERT_NEXUS USER MANUAL	17
3.1 Introduction – Welcome to Convert_NeXus!	17
3.1.1 Features	17
3.1.2 Using this manual	18
3.1.3 Why use the NeXus format?	19
3.2 Installation Instructions	20
3.3 Getting Started	21
3.3.1 Choosing a beam-line (location of files)	22
3.3.2 Converting Files from Beam-lines 1-4	23
3.3.3 Converting Files from Beam-line 5 (Current Directory)	24
3.3.4 Diagnostics	24
3.4 Advanced	26
3.4.1 Selecting Modes	26
3.4.2 Selecting files from beam-lines 1-4	27
3.4.3 Selecting files from Beam-line #5	28
3.4.4 Specifying an Instrument Specific File	28
3.4.5 Including a User Information File	29
3.4.6 Including Collimator Information	30
3.4.7 Entering T0 Information	31
3.4.8 Specifying a Dead-Time File	31
3.4.9 Specifying a grouping file from the command line	32
3.5 Creating a User Information File	33
3.5.1 Example User Information File	35
3.6 Trouble Shooting	36
3.6.1 System Messages	36

3.6.1.1	Error Messages	36
3.6.1.2	Warning Messages	36
3.6.1.3	Diagnostic Messages	38
3.6.2	Hints and Tips	38
4.	CONVERT_NEXUS DESIGN DOCUMENT	39
4.1	General Structure	39
4.1.1	Convert_NeXus Pseudo Code	40
4.2	Source Handler	41
4.2.1	Source Handler Pseudo Code	41
4.2.2	Notes	42
4.3	MACQ Handler	43
4.3.1	MACQ Handler Pseudo Code	43
4.3.2	Notes	44
4.4	TLOG_Handler	45
4.4.1	TLOG Handler Pseudo Code	45
4.4.2	Notes	47
4.5	UIF Handler	48
4.5.1	UIF Handler Pseudo Code	48
4.5.2	Format of UIF	48
4.5.3	Implementation	49
4.5.4	Notes	50
4.6	ISF Handler	51
4.6.1	Instrument Specific File	51
4.6.1.1	Instrument Specific File Pseudo Code	52
4.6.2	Deadtimes	52
4.6.2.1	Dead-Times Pseudo Code	53
4.6.3	Groupings	54
4.6.3.1	Grouping Pseudo Code	54
4.6.4	Time Zero	55
4.6.4.1	Time Zero Pseudo Code	55
4.7	Command Line Handler	57
4.7.1	Command Line Handler Pseudo Code	57
4.7.2	Note	58
4.7.3	Collimator	58
4.8	Histogram Handler	59
4.8.1	Get Data Groups Pseudo Code	59
4.8.2	Split Data Groups Pseudo Code	59

4.8.3	Notes	59
4.9	ISotime	61
4.10	TESTING PROCEDURE FOR CONVERT_NEXUS	63
4.10.1	MCS File	63
4.10.2	Temperature_logs.....	63
4.10.3	MACQ files.....	63
4.10.4	Grouping file	64
4.10.5	Deadtime file.....	64
4.10.6	User Information File	64
4.10.7	Instrument Specific File	65
4.10.8	Collimator Information	65
4.10.9	Time Zero Information.....	65
5.	NeXus_READER.....	67
5.1	Introduction.....	67
5.1.1	Common features	67
5.1.2	ISIS Muon NeXus Instrument Definition	67
5.2	C Version	67
5.2.2	Nexus_reader.c.....	70
5.2.3	Test_nexus_reader.c.....	73
5.2.4	Compilation Instructions	74
5.2.5	Limitations	74
5.3	F77 Reader	75
5.3.1	MUON_DEF.....	75
5.3.2	NEXUS_READER.FOR.....	77
5.3.3	TEST_NEXUS_READER.F77.....	79
5.3.4	Compilation Instructions	80
5.3.5	Limitations	80
5.4	Error Messages.....	81
6.	TMOGGER.....	83
6.1	Tmogger User Manual	83
6.1.1	How To Compile Tmogger	83
6.1.2	How To Set-Up Tmogger.....	84
6.1.3	Getting Started	85
6.2	Tmogger Design.....	89
6.2.1	Tmogger Module.....	89
6.2.2	Get_application_choice ()	90
6.2.3	Plot_macq_data ()	90

6.2.4	Plot_dat ().....	90
6.2.5	Plot_tlog_and_macq_data ()	91
7.	NEXUS_UDA.....	93
7.1	Introduction.....	93
7.2	How To Set Up UDA.....	93
7.3	How to use UDA.....	94
8.	CONCLUSIONS AND PLANS FOR THE FUTURE	97
	FURTHER READING.....	99
	ACKNOWLEDGEMENTS	101
	APPENDIX A: Nexus_reader.h.....	103
	APPENDIX B: Nexus_reader.c	105
	APPENDIX C: Test_nexus_reader.c	117
	APPENDIX D: MUON_DEF_F77.INC	121
	APPENDIX E: MUON_DEF_F90.INC.....	125
	APPENDIX F: NEXUS_READER.FOR.....	129
	APPENDIX G: TEST_NEXUS_READER.FOR	139

1. INTRODUCTION

ISIS is the world's most intense pulsed source of muons and neutrons available for condensed matter research. The MuSR, EMU and ARGUS spectrometers, together with the DEVA development beam, make up the ISIS muon facility; the MuSR spectrometer is shown in Figure 1.1. These instruments are used by researchers from around the world to perform muon spin rotation, relaxation and resonance (μ SR) experiments.




Figure 1.1 – Diagram of the MuSR Spectrometer at the Muon Facility, ISIS

ISIS is one of a number of facilities around the world that are involved in μ SR experiments. The other principle sources are located at Switzerland (PSI), Canada (TRIUMF) and Japan (KEK). A desire for a common format for μ SR data was expressed at a workshop held during the 1999 International μ SR Conference, as it was felt that it would be of great benefit to researchers who are mobile between facilities and would also assist the exchange of ideas and software.

Despite having originally been designed for the use with neutron and X-ray experiments, the flexibility of the NeXus data format makes it equally suitable for storing μ SR data and opens up the possibility of sharing software between these user communities. Indeed, since NeXus uses the freely available Hierarchical Data Format (HDF) (a portable, binary, extensible format that is self-describing) many general tools are already in existence to manipulate NeXus files.

The work presented in this report explores the use of the NeXus format to store μ SR data collected at ISIS. Described are a suite of programs to convert existing ISIS muon binary files into the new format, carry out analysis using these NeXus files and plot the temperature and data logs contained within the new file format. In addition, subroutines to read the NeXus files are described.

1.1 Experiment Data Acquisition

The μ SR technique involves the implantation of spin polarised positive muons into the sample and the subsequent monitoring of the time dependence of the decay positron asymmetry (muons are unstable particles with a lifetime of $2.2\mu\text{s}$). Detected counts are accumulated during the run (which typically lasts 1-2 hours), and are sorted by the Data Acquisition Electronics (DAE) according to the time and detector in which they were measured. At the end of the run the data is written as a binary data file and automatically copied across to the central data archive.

MCS is the software currently used to control the collection of data and sample environment equipment for the ISIS muon spectrometers. Prior to run-time, MCS gathers information about the measurement, including the sample, magnetic field and temperature settings applied to a particular run. Experiments typically involve carrying out many such runs with varying physical parameters.

The ISIS muon facility typically carries out 100 different experiments each year, each lasting about 3-4 days. This activity results in approximately 15000 data files being produced annually, occupying over 3750 MB of disk storage. In addition, log files are written for each run to record the variation in time of a number of physical parameters associated with the experiment. At present this includes the temperature (TLOG) and the number of detected counts over fixed time intervals (MACQLOG).

1.2 Why choose a different format?

The data files written by MCS currently use a binary format that is unstructured and unique to ISIS. While adequate in the past, this method of data storage presents a number of obvious problems for continued use. These include the difficulty of exchanging scientific data between computing institutions, the problem of adding additional information into the file as the μ SR technique expands and the portability of the data files between operating systems. In addition, the logged data recorded during an experiment is stored separately to the corresponding run data file and in a different file format (ASCII). Such an isolation between related data sets is undesirable.

Specifically a new data format should:

- Bring together all information relevant to an experiment, including temperature, data and count logs.
- Be portable and platform independent.

- Be compatible with other non- μ SR techniques to facilitate the sharing of information and existing programs.
- Be based on an industry standard format to take advantage of existing software.
- Be self describing, so that any application that can read the format will be able to read all sections of the data.
- Be structured.
- Be extensible so that additional information can be added to data files without affecting the ability of earlier versions of the read routines to read the data files.

1.3 Why use NeXus?

NeXus is a data format for the exchange of neutron and X-ray scattering data between facilities and user institutions. It has been developed by an international team of scientists and computer programmers from neutron and X-ray facilities around the world. The NeXus format uses the Hierarchical Data Format (HDF) which is portable, binary, extensible and self-describing.

NeXus defines the structure and contents of these HDF files in order to facilitate the visualisation and analysis of neutron and X-ray data. In addition, an Application Program Interface (API) has been produced to simplify the reading and writing of NeXus files.

The NeXus format should offer the μ SR community significant benefits because:

- It provides the opportunity to include a comprehensive description of the experiment within a single data file. Instrument specific information (e.g. the detector arrangement), temperature and count log data can all be added.
- It is a portable, platform independent data format and therefore data can easily be moved between different operating systems encouraging users to visit the facility most suited to their experiment.
- It is compatible with other non- μ SR techniques, encouraging the sharing of information and existing programs.
- It reduces the need for local expertise.
- It reduces the total number of conversion utilities required.

- It reduces redundant software development and enhances co-operation since effort can easily be shared between institutions.
- It simplifies the development of sophisticated visualisation software since existing HDF tools can be adapted.
- It is a binary format enabling efficient storage of the data and allowing the use of the standard data compression routines built into the format.
- It is a format that is structured, self-describing and extensible; data can be added to the files as need arises without affecting the functionality of existing read routines.

Evolving a common μ SR data format is clearly a desirable objective for the μ SR community. The NeXus 2001 workshop (held at PSI, Switzerland, March 2001) demonstrated that there is considerable support for the NeXus format from the neutron and X-ray communities and the upgraded data acquisition system for the ISIS neutron facility will write NeXus files directly during the experiments. For these reasons we have decided to develop the NeXus format for storing muon data and as a first step in this process we have created an application to translate the existing ISIS muon binary data files into the new NeXus format.

This document considers the steps necessary to adapt the NeXus format for use with μ SR data and presents an Instrument Definition (essentially a structure for the NeXus file) suitable for the ISIS muon instruments. The design and use of the conversion utility (Convert_NeXus) is then presented (chapters 3 and 4). Subroutines have been written in both Fortran 77 and C to read the muon NeXus data format and a full description of these routines and their use is given (chapter 5). Finally, example applications to analyse data stored in the muon NeXus format (UDA_NeXus) and also to plot logged temperature values and counts accumulated (Tmogger) are discussed (chapters 6 and 7). Chapters 3 onwards are self contained and are designed to be read independently for information on a particular aspect of the work.

2. ADAPTING THE NeXus FORMAT FOR μ SR DATA

The initial stage in the development of NeXus format for μ SR data required a complete description of the structure and content of the NeXus file (the Instrument Definition). It would clearly be very difficult for a common data format to define all possible data items for every type of muon instrument used in the world. To overcome this, standard Instrument Definitions for a limited number of generic instruments can be devised; the user need then only know the type of definition being used to enable many common characteristics to be read.

In September 2000 a draft Instrument Definition was produced, primarily designed for the muon instruments at ISIS. A NeXus workshop was held at PSI, Switzerland, in March 2001 and this opportunity was taken to discuss this Instrument Definition with other members of the μ SR community. The meeting concluded that it would be difficult to devise a single Instrument Definition that would be appropriate to all muon instruments, but instead defined a subset of essential element names, units, coordinate systems and data types; it was recognised that laboratories would extend these as required. It was also hoped that certain elements, unnecessary for every type of μ SR experiment but likely to be common in a particular field of μ SR, could be standardised (e.g. histogram resolution in muon Time Differential experiments). Following the workshop the ISIS muon Instrument Definition was revised to incorporate suggestions made at the meeting and is now published as version 1.

2.1 ISIS Muon NeXus Instrument Definition

The ISIS muon Instrument Definition is given in the following table. The table shows the structure and content of the ISIS muon NeXus files and specifies the units, coordinate system and data types (including the rank of the data types) for each data item. The final column denotes those data items considered essential elements in a general Instrument Definition for the μ SR community and read routines designed to read only essential items should be able to read a NeXus μ SR data file from any laboratory. It is recognised, however, that each laboratory will want to extend this set of essential data items (this has clearly been done for the ISIS muon Instrument Definition) and produce a correspondingly enhanced set of read routines. It is important to realise that our Instrument Definition will evolve with time, with additional data items added according to need. The read routines published here will, however, continue to read a subset of the information contained in later versions of the ISIS muon NeXus files.

Item	Description	Type	E
Nexus File			
NeXus_version	attribute describing version of NeXus API used to create file	Integer	*
user	attribute denoting scientist who performed experiment	String	*
Nxentry run			
IDF_version	version of IDF that NeXus file conforms to	Integer	*
program_name	name of creating program - "MCS"	String	*
version=	attribute - version of creating program	String	*
number	run number	Integer	*
title	string containing sample, temperature and field	String	*
notes	comment from MCS file	String	*
analysis	type of muon experiment - "muonTD" (muon, time differential)	String	*
lab	origin of experiment - "ISIS" (collected at the ISIS facility)	String	*
beamline	particular beamline used for experiment	String	*
start_time	start time and date of measurement	String	*
end_time	stop time and date of measurement	String	*
duration	calculated duration	String	*
switching_states	"1" - Normal data collection, "2" - Red/Green mode	Integer	*
NXuser user			
name	Scientist(s) name	String	*
experiment_number	RB number	String	*
	Additional information may be added it here using the UIF facility		
NXsample sample			
name	sample name	String	*
temperature	temperature setting	Float	*
units=	attribute, "Kelvin"	String	*
magnetic_field	magnetic Field setting	Float	*
units=	attribute, "Gauss"	String	*
shape	sample orientation	String	*
magnetic_field_state	mode, e.g. "TF"	String	*
magnetic_field_vector	vector describing magnetic field orientation	Float []	*
coordinate_system=	attribute, "cartesian"	String	*
units=	attribute, "gauss"	String	*
environment	rig, e.g. "CCR"	String	*
temperature_log_1	link to log of temperature values obtain from '.TLOG' file	NXlink	
Nxinstrument instrument			
name	instrument name	String	*
NXdetector detector			
number	number of detectors	Integer	*
orientation	detector arrangement, 'Longitudinal ("l")' or 'Transverse ("t")'	Character	
angles	2D array defining detector positions (<i>see Note 1</i>)	Float [] []	
coordinate_system=	attribute, "spherical"	String	*
available=	attribute "1" if angles are available, otherwise "0"	Integer	*
deadtimes	1D array of detector deadtime values (<i>see Note 2</i>)	Float []	
units=	attribute, "microseconds"	String	*
available=	attribute, "1" if deadtime values are available, otherwise "0"	Integer	*
NXcollimator collimator			
type	e.g. "Slits" (not defined in MCS file)	String	*
aperture	e.g. slit setting (not defined in MCS file)	String	
NXbeam beam			
event_log_1	link to log of events obtained from '.MACQLOG' file	NXlink	
total_counts	total number of counts	Float	
units=	attribute, "Mev"	String	*
daereads	number of readouts from DAE	Integer	*
frames	number of ISIS frames collected	Integer	*
NXdata histogram_data_1			
counts	(<i>see Note 3</i>) 2D array of counts: (detector number*switching_states) vs. bin	Integer [] []	*
units=	"counts"; attribute to describe data units	String	*
signal=	"1"; attribute to indicate signal to be plotted	Integer	*
number=	attribute, number of histograms in NXdata group (<i>see Note 4</i>)	Integer	*
length=	attribute, length of histogram	Integer	*
t0_bin=	attribute, t0 bin value for histograms	Integer	*
first_good_bin=	attribute, first good bin values for histograms	Integer	*
last_good_bin=	attribute, last good bin values for histograms	Integer	*
offset=	attribute giving offset to centre of 1 st bin - 0.5*histogram_resolution	Float	*
units=	"picoseconds"; attribute to describe histogram offset units	String	*
histogram_resolution	histogram resolution, set to zero if not applicable	Integer	*
units=	"picoseconds"; attribute to describe resolution units	String	*
time_zero	time zero for muon measurements (<i>see Note 5</i>)	Integer	*
units=	attribute, "microseconds"	String	*
raw_time	scale for time axis (raw time) in microseconds	Float []	*
axis=	"1"; fastest varying index	Integer	*
primary=	"1"; raw time is the default	Integer	*
units=	attribute, "microseconds"	String	*
corrected_time	scale for time axis (corrected for 'histogram_timezero')	Float []	*
axis=	"1"; fastest varying index (2 nd axis)	Integer	*
units=	attribute, "microseconds"	String	*

grouping	1D array defining grouping for histograms in NXdata (see Note 6)	Integer []	*
available=	attribute, number of groups ("0" if information not available)	Integer	*
alpha	2D array, alpha, for pairs of groups defined in 'grouping' (see Note 7)	Float [] []	*
available=	attribute, number of pairs for which alpha is defined ("0" if not in use)	Integer	*
NXlog temperature_log_1	(see Note 8)		
name	name of log, "sample temperature"	String	
available=	attribute, number of temperature values	Integer	
values	log of temperature values obtained from '.TLOG' file	Float []	
units=	attribute, units of temperature logged "Kelvin"	String	
time	time stamp, seconds from start of run	Float []	
units=	attribute, "seconds"	String	
NXlog event_log_1			
name	name of log, "ISIS beam"	String	
available=	attribute, number of event values	Integer	
values	log of events obtained from '.MACQLOG' file	Float []	
units=	attribute, units of events logged "counts"	String	
time	time stamp, seconds from start of run	Float []	
units=	attribute, "seconds"	String	

* denotes elements that must be contained within any μ SR NeXus file

Coordinate systems: Coordinate axes are taken with the z-axis pointing along the beam direction and the x-axis defined as vertically upwards, the origin is taken to be at the centre of the spectrometer (the sample position). Cartesian, cylindrical and spherical polar coordinate systems may be used with angles θ (rotation about the beam axis) and ϕ (rotation about an axis perpendicular to the beam direction) increasing for clockwise rotations viewed from the origin of the coordinate system.

Note 1: 2D array describing detector positions and solid angles. The index of the 1st rank represents the detector number and the first three elements of the 2nd rank define the position in the specified coordinate system. The fourth element of the 2nd rank defines the detector solid angle.

Note 2: The deadtime values for each detector. The values, read from the default deadtime file 'DT(E)PAR.DAT', are stored as a 1D array where the index represents the detector number. The attribute 'deadtimes_available' should be set to "1" if this information is read successfully, otherwise "0".

Note 3: MCS can generate datasets where individual histograms can have different resolutions, lengths, t0 bins, first good bins and last good bins! In practice this feature is rarely used. For efficiency, the muon NeXus instrument definition specifies that sets of histograms where these five parameters are common are grouped into a single NXdata group. Most ISIS muon NeXus files will therefore contain a single NXdata group.

Note 4: The value of the attribute 'histogram_number' should be set equal to the total number of histograms in the NXdata group. The value should reflect any histogram doubling resulting from the Red/Green data collection mode.

Note 5: Not defined in the MCS file, the value should be read from the file 'BASETIME.UDA' if this is available in the current directory. The attribute 'available' should be set accordingly.

Note 6: If there is a grouping file ('long.uda' for longitudinal data and 'trans.uda' for transverse data) in the directory where the NeXus file is being created then this should be used. If this is not available, the default grouping file corresponding to the current instrument orientation should be used. A check should be made that the total number of histograms listed in the grouping file corresponds to the value of 'histogram_number' in the NeXus file, with a mismatch indicating an error. The attribute 'grouping_available' should be set equal to the total number of groups unless grouping information is not available or there is an error condition, when it should be set to "0". The grouping is represented as a 1D array with the index representing the histogram number and the value corresponding to the group into which it should be placed.

Note 7: The balance parameter, alpha, is given for pairs of groups defined in 'grouping'. Values are stored in a 2D array together with the numbers for the forward and backward groups. The attribute 'alpha_available' should be set to the number of alpha values defined or zero if not in use.

Note 8: NXLOG groups can be added as required for any variables logged during an experiment. Logged groups should contain a values array and units and a corresponding time array and units.

2.2 Abstract Model of Instrument Definition

The ISIS Muon Instrument definition is a complete representation of the structure and content of the ISIS muon NeXus files. Summarised below and in Figure 2.2 are the various external elements stored in the muon NeXus file:

- Count log files: the counts accumulated over fixed intervals are stored in a MACQLOG file named according to the run number. This file is converted to a suitable format and included within the NeXus file.
- User Information Files: Convert_NeXus introduces a mechanism for users to include additional user supplied information in the converted NeXus file.
- Temperature log files: temperature readings are taken at fixed time intervals and stored in a TLOG file named according to the run number. However, for historic reasons, this information is frequently spread over several VMS file versions. These versions are

collated together by the Convert_NeXus application and included in the NeXus file, along with the time each temperature reading was taken.

- Instrument Specific File: at present, details of the instrument and apparatus used to perform an experiment is not stored in the ISIS muon data files. Provision has been made to include this information within the ISIS muon NeXus files by developing an Instrument Specific File (ISF) for each of the ISIS muon instruments.
- Collimator information: this enables the horizontal slit setting to be stored. At present, however, there is no way of recording this information directly from the experiment and the information will be included in the ISIS muon NeXus file during the conversion.
- Grouping information: the summation of histograms most suited to the type and method of data collection.
- Dead-time information: after a detector counts a positron there is a very short period when further positrons are missed – the dead-time. In this short interval of time a number of possible muon decays could have occurred and been missed. Calibration measurements are regularly carried out to quantify the dead-time for each detector, these values can then be used to correct data during experiments. The dead-time information for appropriate spectrometer is incorporated within the ISIS muon NeXus file.




Figure 2.2 – A Software Engineering Context diagram, denoting the interactions between the terminators and the Convert_NeXus application.

2.3 Data Storage using the ISIS Muon NeXus Instrument Definition

An example ISIS muon NeXus binary data file viewed using the Fortner HDF Browser is shown in Figure 2.3. The hierarchical nature of the data entries can clearly be seen (left most window pane, labelled '30000'), with the entries following the scheme given by the ISIS muon Instrument Definition. The remaining window panes show the contents of various data items, e.g. the array of histogram counts, the array of detector deadtimes and the run note.




Figure 2.3 – Muon NeXus file viewed using the Fortner HDF browser.

3. CONVERT_NEXUS USER MANUAL

3.1 Introduction – Welcome to Convert_NeXus!

Convert_NeXus is a Fortran 90 application that converts the current ISIS muon data files from an inextensible binary format into the structured HDF-based, binary NeXus format. This application has been specially designed to run on the VMS operating system at ISIS, and conforms to the ISIS muon Instrument Definition version 1 (described in Chapter 2). An Instrument Definition is a comprehensive description of the structure and contents of a NeXus file.

3.1.1 Features

Beyond the clear advantages of adopting a common data format, there are a number of other reasons why it is beneficial to use Convert_NeXus. Of particular value is the ability to include elements which have previously resided separate to the run data files. A brief description of the numerous features of Convert_NeXus are listed below and detailed information on the various options can be found in section 3.4:

- **MACQ Log files** – Created separately to the run data files are log files which store the flow of events from the data acquisition electronics. This information is automatically selected by Convert_NeXus and included within the NeXus file, along with a calculated time stamp for each entry in the log file.
- **Temperature Log files** – During a run, temperature readings are taken at fixed time intervals and stored in temperature log files. Due to legacy reasons this information is frequently spread over several versions. These versions are collated together and included in the NeXus file, along with the calculations of the time at which each temperature reading was taken.
- **Instrument Specific File** – Convert_NeXus allows incorporation of information concerning the Instrument used to perform the experiment. This information principally consists of detector geometries of the spectrometer.
- **User Information Files** – Convert_NeXus introduces a mechanism for users to include additional information. This would not have been possible with the original binary data files.

- **Collimator Information** – At present the collimator type and aperture information are not electronically recorded during an experiment. However it is possible to include this information within the NeXus file during conversion.
- **Grouping Information** – Histograms are usually grouped for analysis. Grouping information is extracted from an appropriate grouping file that is dependent on the orientation of the spectrometers detector array. Values for alpha are also obtained from this grouping file.
- **Dead-time Information** – Each detector can have a dead-time associated with it, which is a missed count amendment factor. These dead-times are calibrated periodically and stored in a dead-time file categorised by instrument.
- **T0 (Time zero)** – Time Zero signifies the exact time of the centre of the muon pulse. This Time Zero information is taken periodically and stored in a basetime file, which is where `Convert_NeXus`, by default, searches for this value. It is also possible for the user to include their own T0 value.

3.1.2 Using this manual

This is a ‘feature oriented’ manual. Each subsection covers a major `Convert_NeXus` topic. It is intended to be both as a tutorial introduction to `Convert_NeXus` and as a reference manual. The remainder of this section describes why we have chosen to adopt this format.

Section 3.2 describes how to install `Convert_NeXus` from scratch. This involves compiling source files and linking to the NeXus API and HDF libraries. This is for the involved programmer. Read at your peril!

Section 3.3 demonstrates a step-by-step tutorial on how to convert a few example muon data files into the NeXus format.

Section 3.4 introduces some of the more advanced concepts when using `Convert_NeXus`. This chapter is divided into sections, each carefully explaining the available features of `Convert_NeXus` and how to utilise them.

Section 3.5 demonstrates how to create a User Information File using `UIF_creator`, and illustrates the constituent parts of a User Information File.

Section 3.6 is a trouble shooting guide, which includes hints and tips on the operation of `Convert_NeXus`, and describes the system messages that `Convert_NeXus` produces.

3.1.3 Why use the NeXus format?

As instrumentation becomes more complex and data visualisation becomes more challenging, individual scientists, or even institutions, have found it difficult to keep up with new developments. The introduction of the NeXus format has been designed to alleviate these difficulties.

There are many advantages of a common format:

- Reduce the need for local expertise – Because different institutions will inevitably store data using a unique format, the Scientist/Programmer must learn how to obtain meaningful data from a range of existing data files formats. Producing a standardised format with standardised notations will help avoid this.
- Reduce the number of conversion utilities – To electronically compare and use data from different institutions requires the conversion of the data files into some common format. Writing individual tools is not a problem, but with the growing number of research groups performing experiments we are faced with an explosion of the number of converters.
- Reduce redundant software development – With a common data format, data visualisation tools could easily be moved between institution reducing the amount of duplicate software. Establishments can easily share effort to produce increasingly sophisticated software.
- Increase cooperation in software development – Many techniques of data manipulation are common in different scientific communities. These communities will be able to share software solutions to common problems.
- Increased functionality of generic software – NeXus is developed on top of the HDF format. This means that we can immediately take advantage of the tools already developed for this format.

Despite having been originally designed for the use with neutron and X-ray data, the flexibility of the NeXus data format makes it equally suitable for storing μ SR data and opens up the possibility of sharing software between these user communities.

3.2 Installation Instructions

Convert_NeXus has been specifically designed to run on the Open VMS operating system at ISIS. To recompile the source code the following files will be required along with a standard Fortran 90 compiler:

- | | |
|----------------------------|--|
| • NXmodule.f90 | NeXus API |
| • NXUmodule.f90 | NeXus utility API |
| • ISOtime.f90 | Exclusively written date and time module |
| • Command_line_handler.f90 | Module to handle command line arguments |
| • Source_handler.f90 | Module to handle MCS file |
| • ISF_handler.f90 | Module to handle Instrument Specific Files |
| • UIF_handler.f90 | Module to handle User Information Files |
| • TLOG_handler.f90 | Module to handle TLOG files |
| • MACQ_handler.f90 | Module to handle MACQ files |
| • Histogram_handler.f90 | Module to handle histogram properties |
| • Convert_nexus.f90 | Main Application |

In these instructions it is assumed that the HDF libraries and NeXus API's are currently installed. If these are not yet installed please refer to the NeXus Data Format Homepage (http://www.neutron.anl.gov/nexus/NeXus_API.html#Install) for detailed instructions.

Step 1 – Compile above files in sequence

Step 2 – Link Convert_Nexus, ISOtime, Tlog_handler, Macq_handler, Uif_handler, Isf_handler, Command_line_handler, Source_handler, Histogram_handler, NeXus API (where Nexus Library exists)

See Figure 3.1 overleaf for DEC Fortran 90 compilation instructions on the Open VMS operating system at ISIS.

```

$set verify
$define hdf_root applib$disk:[4.1r3_alphavms7.1.] ! where you installed HDF, trailing "." is needed
$define nexus_root user$disk:[df01.] ! where the nexus library exists
$F90 NXmodule.f90
$F90 NXUmodule.f90
$F90 ISOtime.f90
$F90 command_line_handler
$F90 source_handler
$F90 isf_handler
$F90 uif_handler
$F90 tlog_handler
$F90 macq_handler
$F90 histogram_handler
$F90 convert_nexus
$LINK convert_nexus,isotime,tlog_handler,macq_handler,error_handler,uif_handler,isf_handler,command_line_handler,source_handler,
histogram_handler,nexus_root:[nexus]NEXUS/LIB,HDF_ROOT:[LIB]MFHDF/LIB,DF/LIB,LIB2/LIB,LIBJPEG/LIB

```

Figure 3.1 – DEC Fortran 90 Compilation Instructions Screenshot

3.3 Getting Started

This chapter introduces the basic concepts of Convert_NeXus, enabling the user to run the application using minimal effort. More advanced topics are covered in subsequent chapters.

To run the application the user must first initialise the application. This can be achieved by typing: “@ *muon\$disk:[muonmgr.software.utilities.NeXus]NeXus_setup*” at the VMS command prompt (See Figure 3.2).

The above command will execute a command file which will install a suite of Muon NeXus applications in the user’s current directory. The applications installed are:

- *Convert_NeXus* – Application which converts ISIS Muon MCS data files into NeXus format.
- *NXbrowse* – Open source NeXus file browser.
- *NeXus_UDA* - μ SR Data Analysis program which can interpret ISIS Muon NeXus files.
- *F77_NeXus_Reader* – ISIS Muon NeXus read routine, written in F77 (allows users to read NeXus data into their own programs).
- *C_NeXus_Reader* - ISIS Muon NeXus read routine written in C (allows users to read NeXus data into their own programs).
- *Tmogger* – Application which plots temperature and event logs stored in ISIS Muon NeXus files.
- *UIF_creator* – Application which produces UIF files for Convert_NeXus.

It should also be noted that the command file copies across any additional information needed by these applications including a selection of test data, and will also initialise any logical variables needed.

We will only concern ourselves with the ‘Convert_NeXus’ application here. For more information on the other available software please refer to the Muon NeXus homepage (www.isis.rl.ac.uk/muons/nexus).

To execute Convert_NeXus just type “*CONVERT_NEXUS*”

```

DECterm 1
File Edit Commands Options Print Help
$ @MUON$DISK:[MUONMGR.SOFTWARE.UTILITIES.NEXUS]NEXUS_SETUP

Alpha versions of software being set up . . .

Copying test data to current directory...

Commands available:
  CONVERT_NEXUS - Convert Raw data files to Nexus format
  NXBROWSE - Browse converted Nexus files
  NEXUS_UDA - Data analysis program that can also read Nexus files
  F77_NEXUS_READER - Read Nexus files and output data to screen
  C_NEXUS_READER - Read Nexus files and output data to screen
  TMOGGER - Plot temperature and event log information from Nexus files
  UIF_CREATOR - Create User Information Files for CONVERT_NEXUS application

$

```

Figure 3.2 – Muon / NeXus Utilities Screenshot

3.3.1 Choosing a beam-line (location of files)

Once the application is underway, a welcome message appears along with a short menu. A prompt appears asking the user to choose a beam-line. The location of the binary Muon Data Files are categorised by their beam-line.

<u>Beam-line</u>	<u>Location</u>
1. EMU	“emu\$disk0:[data.emu]”
2. MUSR (old data sequence)	“musr\$disk0:[data.musr]”
3. MUSR (new data sequence)	“musr\$disk0:[data.musr]”
4. MUT	“mut\$disk0:[data.mut]”
5. TEST_DATA	[<i>Current Directory</i>]

For example, data recorded from experiments performed on the EMU spectrometer are located under beam-line 1, data recorded from the MUSR spectrometer will be stored under beam-lines

2 and 3, data recorded from the DEVA spectrometer are stored under beam-line 4 (see Figure 3.3).

Beam-line 5 allows you to specify MCS files from your current directory where we have placed specially selected muon data files with differing characteristics to demonstrate the capability of Convert_NeXus.

The files included within this directory are:

- EMU30000.RAW – which represents a run with red/green data collection.
- R35000.RAL – which represents a conventional non red/green run.
- MUT00477.RAW – which represents a run containing differing histogram properties.





Figure 3.3 – Beam-line Menu Screenshot

3.3.2 Converting Files from Beam-lines 1-4

After selecting a beam-line between 1 and 4 (here option 1 has been chosen, the EMU beam-line) a prompt appears asking the user for the range of files to be converted. Each time an experiment is performed it is assigned a unique 5-digit run number. On this screen the user must firstly enter the run number of the first file to be converted. If this is a valid 5-digit number the user will then be able to enter the last file to be converted. This must also be a 5-digit number and must be greater than or equal to the run number of the first file.

Convert_NeXus then commences to convert the specified range of files. If it encounters any problems reading an MCS file it discards this file and continues to process the remaining files in that range.




Tip:

If the user only wishes to convert a single file then they must make the first file equal to the last file.

Figure 3.4 – File Range Entry Screenshot

3.3.3 Converting Files from Beam-line 5 (Current Directory)

If you select beam-line 5, a slightly different prompt appears, asking the user for the full filename of the MCS file they wish to convert. Convert_NeXus then searches the current directory for the specified file. It should be noted that by choosing this beam-line you can only specify one MCS file at a time.



Note:

Convert_NeXus also looks in the current directory for TLOG and MACQLOG files when beam-line 5 is chosen


Figure 3.5 – File Selection Screenshot

3.3.4 Diagnostics

Once the user has selected the appropriate file(s) the conversion process will begin. While processing each file Convert_NeXus will display selected data read from each MCS file processed. Diagnostic information will also be shown for each major process undertaken by Convert_NeXus.

Such processes include:

- Reading and manipulation of Temperature log information.
- Reading and manipulation of Event log information.
- Processing of User Information Files.
- Inclusion of Dead-time, Time-Zero and Grouping information.
- Inclusion of Detector geometries read from appropriate Instrument Specific File.

The image is a screenshot of a terminal window titled "DECterm 2". The window has a menu bar with "File", "Edit", "Commands", "Options", "Print", and "Help". The main text area displays diagnostic messages from the Convert_NeXus software. The messages indicate that an MCS file 'emu\$disk0:[data.emu]r35000.ral' is open and read successfully. It then displays the contents of the MCS file as a list of parameters and their values. After a separator line, it shows that a Nexus file ':35000.nxs' is open and that a TLOG file 'emu\$disk0:[data.emu.tlog]r35000.tlog:1' is also open and read successfully.

```
DECterm 2
File Edit Commands Options Print Help

DIAGNOSTIC : MCS file 'emu$disk0:[data.emu]r35000.ral' is open
DIAGNOSTIC : reading MCS file
DIAGNOSTIC : MCS file read succesfully
DIAGNOSTIC : closing MCS file

DIAGNOSTIC : contents of MCS file ...

t0bins      -->      19
fgbins      -->      26
lgbins      -->     1000
his len (1) -->     1000
resolution (1) -->    16000
hisnum      -->      32
reads       -->      14
frames      -->     5000
det no      -->      32
run no      -->     35000
mcs version -->     RD
comment     -->  slits=10mm, zf
grp name    -->  PAVIA
rb num      -->   11293
duration    -->  00:23:39
analysis    -->  muonTO
shape       -->   n/a
start_date  -->  2000-07-14 11:14:50
finish date -->  2000-07-14 11:38:29
progrnam    -->  MCS
title       -->
Temperature = 10 K, Magnetic Field = 0 Gauss, Sample = LACUO
det orientation --> 1
field_state  -->  n/a
instrument name --> EMU
rig          -->  n/a

DIAGNOSTIC : Nexus file name is :35000.nxs
DIAGNOSTIC : Nexus file open

-----

DIAGNOSTIC : TLOG filename is 'emu$disk0:[data.emu.tlog]r35000.tlog:1
DIAGNOSTIC : TLOG file open
DIAGNOSTIC : TLOG file is within NeXus run time
DIAGNOSTIC : TLOG file read succesfully
DIAGNOSTIC : closing TLOG file
```

Figure 3.6 – Screenshot denoting some of the diagnostic information displayed by Convert_NeXus.

This particular example shows information read from EMU MCS file r35000.ral and commentary from Temperature Log processing

3.4 Advanced

This section explains how to utilise the following functions available using Convert_NeXus:

- Selecting Modes.
- Selecting beam-lines and files from the command line.
- Specifying an Instrument Specific File.
- Including a User Information File.
- Entering Collimator Information.
- Entering Time-Zero Information.
- Specifying a Dead-Time File.
- Specifying a Grouping File.

All of these features are available as command line parameters when executing Convert_NeXus.

3.4.1 Selecting Modes

Convert_NeXus has the ability to run under three different modes, which effect the display of information during the conversion process. This gives the user the flexibility to run a Convert_NeXus program tailored to their needs.

Mode 0 – This is a stealth or silent mode. When running Convert_NeXus in Stealth mode, only fatal error messages, should they occur, will be displayed on screen. If the user wishes to run the application in Stealth mode, this MUST be indicated as a command line parameter along with the beam-line and selected file(s) (see the following sections for more information). This mode has the advantage that Convert_NeXus can be run as a script.

Mode 1 – This mode allows Menu options and Warnings to be displayed. The only difference is that no diagnostic information is displayed as Convert_NeXus executes.

Mode 2 – This is the default mode, allowing the user to view selected information read from the binary muon data files and diagnostics from critical processes.

<i>Option</i>	<i>Qualifiers</i>	<i>Arguments</i>	<i>Example</i>
MODE	2	MODE=2	
	1	MODE=1	
	0	MODE=0	(NB. Additional parameters must be included when using Stealth mode, 4.2 / 4.3)

Figure 3.7 – Mode Selection Table




Figure 3.8 – Mode Selection Screenshot

3.4.2 Selecting files from beam-lines 1-4

If the user wishes to run Convert_NeXus in stealth mode, using beam-lines 1-4, then the beam-line along with the first and last file to be converted must be supplied. If this format is not followed, Convert_NeXus will display an error message indicating which parameters were missing or incorrect. The format of these arguments is described below:

<i>Option</i>	<i>Qualifiers</i>	<i>Arguments</i>	<i>Example</i>	<i>Note</i>
BEAMLINE	1	BEAMLINE=1	EMU	
	2	BEAMLINE=2	MUSR (old data sequence)	
	3	BEAMLINE=3	MUSR (new data sequence)	
	4	BEAMLINE=4	MUT	
FIRSTFILE	Integer 10000 - 99999	FIRSTFILE=30000		
LASTFILE	Integer 10000 - 99999	LASTFILE=30005		

Figure 3.9 – Beam-line and File Range Selection Table




Figure 3.10 – Beam-line and File Range Selection Screenshot

3.4.3 Selecting files from Beam-line #5

As mentioned earlier, Beam-line 5 is handled slightly different to beam-lines 1-4. Instead of supplying the run numbers of the file to be processed, the user must instead supply the actual filename of the data file. The following table Figure 3.11 and screenshot Figure 3.12 show how to choose a file from beam-line 5 as a command-line parameter.

Option	Qualifiers	Arguments	Example	Note
BEAMLINE		5	BEAMLINE=5	TEST_DATA
TST		<filename>	TST=mut00477.raw	

Figure 3.11 – Beam-line and File Selection Table




Figure 3.12 – Beam-line and File Selection Screenshot

3.4.4 Specifying an Instrument Specific File

At present, the detector geometry of the instrument used to perform an experiment is not stored in the binary Muon Data Files. It is desirable to include this information within the NeXus files so that utilities that read these NeXus files are able to take into account the properties of the experiment apparatus. We have therefore developed an Instrument Specific File (ISF) for each of the Muon Instruments at ISIS.

When Convert_NeXus executes, it automatically selects an appropriate Instrument Specific File according to which spectrometer was used to perform the experiment. It is possible, however, to over-ride this option and manually include an Instrument Specific File of your choice. However, this is ill-advised!

The format for this command-line switch is detailed below:

<i>Option</i>	<i>Qualifiers</i>	<i>Arguments</i>	<i>Example</i>	<i>Note</i>
ISF		<filename>	ISF=EMU_specific.isf ISF=MUSR_specific.isf ISF=DEVA_specific.isf	

Figure 3.13 – ISF Selection Table




Figure 3.14 – ISF Selection Screenshot

3.4.5 Including a User Information File

Although the current ISIS muon data files contain certain environment conditions such as the temperature and the magnetic field, some experimenters wish to incorporate other information which is specific to their experiment (one example could be light frequencies if the measurement involves any stimulation by light). A mechanism has been devised to enable the user to store additional information in the NeXus data file. This is achieved by producing a User Information File (UIF) which is read during conversion, the file contains a description of the information to be included within the Nexus File.

To produce a User Information File see section 3.5.

The format for the command-line switch to include a User Information File is detailed below:

Option	Qualifiers	Arguments	Example	Note
UIF		<filename>	UIF=test_uif.uif	Example using absolute path and of UIF stored in a different directory.
			UIF=user\$disk:[df01.uif]test_uif.uif	

Figure 3.15 – UIF Selection Table



Figure 3.16 – UIF Selection Screenshot

3.4.6 Including Collimator Information

All of the Muon instruments have a slit that defines the horizontal size of the beam. At present there is no way of recording this information directly from the experiment. However it is possible to include this information in the NeXus file during the conversion. If the collimator type and aperture values are not supplied via the command line, Convert_NeXus defaults the Collimator type to "Slits" and the collimator aperture to "0". See below for command line inclusion format:

Option	Qualifiers	Arguments	Example	Notes
COL	TYPE	<type of collimator>	COL.TYPE=slits	
	APT	<aperture>	COL.APT=1.5	

Figure 3.17 – Collimator Specification Table




Figure 3.18 – Collimator Specification Screenshot

3.4.7 Entering T0 Information

Time_zero specifies the exact time of the centre of the muon pulse relative to the time the data acquisition electronics was triggered. This information is not currently stored in the MCS files, but is available from a file called *basetime.uda*. On execution, Convert_NeXus searches for this file in the current directory, and if successful, includes it within the NeXus file.

It is possible for the user to supersede this value and specify a particular Time Zero floating point value at the command line. If the supplied floating point number is invalid a warning message will be displayed and the number will not be included within the NeXus file. The format of this argument is illustrated below:

<i>Option</i>	<i>Qualifiers</i>	<i>Arguments</i>	<i>Example</i>
T0		<floating point number>	T0=2.45

Figure 3.19 – T0 Specification Table



Figure 3.20 – T0 Specification Screenshot

3.4.8 Specifying a Dead-Time File

Each detector can have a dead-time associated with it. When a detector on one of the muon spectrometers detects a positron there is a very short period following when further positrons are missed. This dead-time is calibrated and is usually incorporated into the experimental results to compensate for any missed counts. While processing a file, Convert_NeXus selects an appropriate dead-time file from the current directory depending on the instrument used to perform the experiment. If the user wishes they can override this and specify their own choice of dead-time file. Although this is not advisable, details of how to do this are shown below:

<i>Option</i>	<i>Qualifiers</i>	<i>Arguments</i>	<i>Example</i>	<i>Note</i>
DEADTIME		<filename>	DEADTIME=dtepar.dat	

Figure 3.21 – Dead-time Specification Table




Figure 3.22 – Dead-time Specification Screenshot

3.4.9 Specifying a grouping file from the command line

Histograms are usually grouped for analysis. Grouping information is extracted from an appropriate grouping file that is dependent on the orientation of the spectrometer's detector array. One example would be detectors that are in forward and backward positions in relation to the sample. Convert_NeXus includes grouping information for traditional longitudinal and transverse measurements. The user is able to specify their own grouping file from the command line by using the command line option detailed below:

<i>Option</i>	<i>Qualifiers</i>	<i>Arguments</i>	<i>Example</i>	<i>Note</i>
GROUPING		<filename>	GROUPING=trans.uda	

Figure 3.23 – Grouping Specification Table





Figure 3.24 – Grouping Specification Screenshot


3.5 Creating a User Information File

We have produced a very simplistic application called UIF_creator, which guides the user through a step-by-step process in creating a User Information File. This application is written in Fortran 90 and is available when you run the NeXus set-up file as discussed in Section 3.



1


To Invoke the application, type “uif_creator” at the command line



2

Once the application has commenced, it asks the user to enter a name for the element they wish to store e.g. “Light_frequency”.


This name must have no spaces. Once a correct entry has been accepted, a prompt appears asking the user to enter a value e.g. “1.5MHz”. This may be up to 30 characters long and will be stored as a string in the NeXus file.



3

once the name and the value of the element has been specified, the user should enter the actual run numbers that they wish this


information to be included in. When a User Information File is included in Convert_NeXus it cross references these run numbers with the run numbers of the files it is processing. If they match then the information is included. The run numbers must be 5 digits.



4

It is more than likely that the user will want a particular element stored in more than one run. For this reason UIF_creator


asks the user if they wish to enter another run number. If so the previous prompt appears again (as in step 3) asking the user to enter a run number.



5


Once the user has entered all run numbers and chosen not to enter any more, another prompt appears, enquiring if

the user would like to add another element to the UIF. If so the process is repeated from step 2.



6

When the user has finished adding elements to the UIF, a summary of the UIF elements is printed to the screen (as left). The user is then asked if they wish to save the UIF file.



7

Assuming the user is happy with the information, they can enter a filename and save the information. (Note: we recommend

that you adopt the file extension “.uif” to avoid any confusion).

3.5.1 Example User Information File

Below is a screen shot of a typical User Information File with a description of each entry. This may be useful for the experienced user.




Figure 3.25 – User Information File Screenshot

3.6 Trouble Shooting

This short chapter is a trouble shooting guide for the user looking for hints, tips or an understanding of some of the display messages produced by Convert_NeXus.

3.6.1 System Messages

Convert_NeXus will have three varieties of system messages. Error Messages, Warning Messages and Diagnostic Messages. These system messages will all conform to a specific outline format of: Error Code, followed by, Error Message. e.g. "ERROR : Beam line must be between 1 and 5".

3.6.1.1 Error Messages

Irrespective of which mode Convert_NeXus is running under, if the application encounters any errors it will display an error number and the corresponding message alongside. The possible error messages that can occur are listed below:

<i>Code</i>	<i>Message</i>	<i>Notes</i>
Error	Must enter 5 digit number	When entering Run number
Error	Last file must be $\geq$ first file	
Error	Beam line must be between 1 and 5	
Error	TST filename must be supplied for TST beam-line	When beam line = 5
Error	First file must be supplied in stealth mode	When verbose mode = 0
Error	Last file must be supplied in stealth mode	When verbose mode = 0

Figure 3.26 – Convert_NeXus Error Messages Table

3.6.1.2 Warning Messages

If Convert_NeXus is running under Mode 1 or mode 2 then it will display warning information as the application executes.

<i>Code</i>	<i>Message</i>	<i>Notes</i>
Warning	Specified UIF does not exist	UIF file supplied via the command line does not exist
Warning	Problem reading UIF	UIF is corrupt
Warning	User Information Not included in NeXus file	
Warning	Specified MCS file does not exist	Specified MCS file is missing
Warning	Problem reading MCS file	MCS file is corrupt
Warning	Default time zero file not found in current directory	

Warning	Supplied T0 value is invalid	Specified T0 value is invalid floating point number
Warning	Problem reading time zero file	T0 file is corrupt
Warning	Default ISF file not found	
Warning	Non EMU ISF being used for EMU experiment	Wrong ISF being used
Warning	Non MUSR ISF being used for MUSR experiment	Wrong ISF being used
Warning	Non MUT/DEVA ISF being used for MUT experiment	Wrong ISF being used
Warning	Cannot read ISF file	Corrupt ISF
Warning	Specified DEADTIME file does not exist	User specified dead-time file does not exist
Warning	Suitable dead-time file not found in current directory	Default dead-time file is missing from current directory
Warning	Non EMU dead-time file being used for EMU experiment	Wrong dead-time file being used
Warning	Non MUSR dead-time file being used for MUSR experiment	Wrong dead-time file being used
Warning	Non DEVA/MUT dead-time file being used for MUT experiment	Wrong dead-time file being used
Warning	Can not read dead-time file	Dead-time file is corrupt
Warning	Dead-time file unsuccessfully read	Summary of above
Warning	Dead-time file not included in NeXus file	Summary
Warning	Specified GROUPING file does not exist	User specified GROUPING file Not found
Warning	Default GROUPING file does not exist	Default grouping file does not exist
Warning	Transverse grouping file being used for longitudinal experiment	Wrong grouping file being used
Warning	Longitudinal grouping file being used for transverse experiment	Wrong grouping file being used
Warning	Cannot read grouping file	Corrupt grouping file
Warning	Grouping file unsuccessfully read	Summary of above
Warning	Default grouping file does not match histograms	If red/green mode or if file has differing histogram properties.
Warning	MACQlog file does not exist	Cannot find corresponding macq log file
Warning	MACQ file only has one value, therefore not included in NeXus file	Not worth including in NeXus file
Warning	Can not read macq file	Macq file is corrupt
Warning	TLOG file does not exist	Cannot find any corresponding tlog files
Warning	Problem reading TLOG file version ??	Version of TLOG is corrupt

Figure 3.27 – Convert_NeXus Warning Messages Table

3.6.1.3 Diagnostic Messages

If Convert_NeXus is running in Mode 2, diagnostic information will be displayed commenting on each of the major processes. All diagnostic information will take the form below:

<i>Code</i>	<i>Message</i>	<i>Notes</i>
Diagnostic	?? File is ??	e.g. UIF file is test_uif.uif
Diagnostic	?? File open	e.g. UIF file open
Diagnostic	?? File read successfully	e.g. UIF file read successfully
Diagnostic	Closing ?? file	e.g. Closing UIF file

Figure 3.28 – Convert_NeXus Diagnostic Messages Table

3.6.2 Hints and Tips

If you are running Convert_NeXus from your own user directory at ISIS it may be beneficial to increase your working set size (see local computing staff).

If converting large and/or frequent data sets it may be advisable to increase your disk quota (See local computing staff). Running out of disk space while converting files causes Convert_NeXus to halt.

4. CONVERT_NEXUS DESIGN DOCUMENT

4.1 General Structure

Convert_NeXus is written in Fortran 90 and uses the Fortran 90 and utility NeXus API. The required functionality has been abstracted into separate modules. These modules are not completely self-contained, but give the application a flexible design and allow these modules to be removed and/or replaced with relative ease. The various modules are listed below:

- MACQ Handler
- Temperature Handler
- Command Line Handler
- Histogram Handler
- UIF Handler
- Source Handler
- ISF Handler
- ISotime

These modules will be described in subsequent sections.




Figure 4.1 – Modular Program Structure Diagram

Below is some primitive Pseudo code, which depicts the flow of events through the Convert_NeXus program.

4.1.1 Convert_NeXus Pseudo Code

```
Get command line switches ()
Get Beamline ()
Get File Range ()

Read UIF ()
Read T0 file ()

Loop through file range
  If (file doesn't exist) then
    Display warning message
  Else
    Create NeXus file ()
    Open Nexus file ()
    Read MCS File ()
    Read MACQ log ()
    Read Temperature log ()
    Read Instrument Specific File ()
    Read Dead-time file ()
    Read Grouping file ()
    Read T0 information ()
    Convert Variables to NeXus Format () //e.g. ISO time
    If (differing histogram properties) then
      Split data groups ()
    End if
    Close NeXus file ()
  End if
End loop
```


4.2 Source Handler

The source handler currently has only one main method, which reads MCS files. The module has been designed so that in principle more than one data format can be read. The global variables below are the elements read from the MCS file. These can be accessed from the main Convert_NeXus application and directly included within the NeXus file.




Figure 4.2 – Source Handler Class Diagram

4.2.1 Source Handler Pseudo Code

```

Declare Equivalence block
Open File
Read header information into memory
Read data into memory
Index memory using equivalence block
Close File
  
```



```

DECterm 1
File Edit Commands Options Print Help

Dump of file MUSR$DISK0:[DATA.MUSR]R30000.RAL;1 on 21-AUG-2000 13:56:11.81
File ID (34479,1,0) End of file block 253 / Allocated 261

Virtual block number 1 (00000001), 512 (0200) bytes

00040005 00000000 00004452 00030402 ....RD..... 000000
03E80004 00007530 00000000 00003E80 .>.....Ou....è. 000010
03E803E8 03E803E8 03E803E8 002003E8 è. .è.è.è.è.è.è.è. 000020
03E803E8 03E803E8 03E803E8 03E803E8 è.è.è.è.è.è.è.è. 000030
03E803E8 03E803E8 03E803E8 03E803E8 è.è.è.è.è.è.è.è. 000040
03E803E8 03E803E8 03E803E8 03E803E8 è.è.è.è.è.è.è.è. 000050
00000000 00000000 00000000 03E803E8 è.è..... 000060
00000000 00000000 00000000 00000000 ..... 000070
47481388 00180000 03E80020 00000000 .... .è.....RG 000080
30302E30 32202000 20202020 20202031 1 . 20.00 000090
454E4F4E 00743030 2E303035 20200030 0. 500.00t.NONE 0000A0
00000000 00000000 00002020 20202020 ..... 0000B0
20205253 554D0000 00000000 00000000 .....MUSR 0000C0
38320000 00000000 00000000 20202020 .....28 0000D0
392D5941 4D2D3832 0039392D 59414D2D -MAY-99.28-MAY-9 0000E0
3A38333A 39313034 3A35353A 38310039 9.18:55:4019:38: 0000F0
00000000 00000000 00000000 00003833 38..... 000100
00000000 00000000 00000000 00000000 ..... 000110
RETURN/SPACE=More, PREV/NEXT=Scroll, INS/DEL=Pan, SELECT=00/132, CTRL/Z=Quit

```

Figure 4.3 – MCS File Screenshot

4.2.2 Notes

The entire data file is read into memory. This ensures that the data is placed in memory with the least amount of I/O as possible. By using a Fortran Equivalence statement it is possible to cross reference main memory to extract the individual and contiguous elements of the data file. The binary muon data file is split into two records. The first record is read into memory and the header information is indexed in the appropriate manner. The second record contains the actual histogram data. Its rank and dimensions can be determined from the header information. A dynamically allocatable two-dimensional array has therefore been created to store the histogram data.

4.3 MACQ Handler

The MACQ handler module consists of one method which takes a string value representing the filename of a particular MACQLOG file, and returns an integer representing the status of the method i.e. 1 on success or 0 on failure. Meanwhile, assuming the filename is valid, the method reads the desired event log information, processes the data and stores the values in the appropriate global variables, shown below.




Figure 4.4 – MACQ Handler Class Diagram

4.3.1 MACQ Handler Pseudo Code

```

Find existing macq files with matching run number
Open MACQ log file
Read header information into memory
While not found
{
    If (correct run and run times are within nexus run times)
    {
        Allocate space for data
        Read time info
        Read counts
        Close file
        Found = true
    }
}

```

MACQLOG - from: 09:41:14 To: 10:01:17

```

57 0.0000
1200.0000
0.0000
0.000000
(7(1X,F10.5)) (7(1X,F10.0))
0.00000 0.05408 0.05331 0.05371 0.05366 0.05402 0.05395
0.05373 0.05394 0.05336 0.05376 0.05370 0.05385 0.05417
0.05424 0.05374 0.05414 0.05386 0.05355 0.05422 0.05407
0.05411 0.05383 0.05320 0.05366 0.05391 0.05383 0.05371
0.05392 0.05344 0.05389 0.05373 0.05337 0.05376 0.05424
0.05413 0.05378 0.05357 0.05358 0.05320 0.05326 0.05390
0.05371 0.05406 0.05365 0.05360 0.05369 0.05380 0.05367
0.05354 0.05345 0.05383 0.05417 0.05366 0.05413 0.05450
0.05375
0. 22. 43. 65. 86. 107. 129.
150. 172. 193. 215. 236. 257. 279.
300. 321. 343. 364. 386. 407. 429.
450. 471. 493. 514. 536. 557. 578.
600. 621. 642. 664. 687. 709. 730.
751. 773. 794. 816. 837. 858. 880.

```

Buffer: R30001.MACQLOG | Write | Insert | Forward

EDT keypad defined (for more information, see help on EDT DIFFERENCES).

Figure 4.5 – MACQ File Screenshot

4.3.2 Notes

It is possible that more than one version of the MACQLOG file exists of a particular run. This may occur if a run is restarted due to unforeseen circumstances. It is therefore necessary to inspect all versions of a particular MACQLOG file and ensure that the run times lie within that of those stored in the NeXus file. There will, however, only be one correct version of the MACQLOG file. Once this is found no other MACQLOG files need be processed.

4.4 TLOG_Handler

The TLOG Handler module also consists of a single method which takes a string value representing the filename of a particular TLOG file, and returns an integer representing the status of the method i.e. 1 on success or 0 on failure. Meanwhile, assuming the filename is valid, the method reads the desired temperature log information, processes the data and stores the values in the appropriate global variables, shown below.




Figure 4.6 – TLOG Handle Class Diagram

4.4.1 TLOG Handler Pseudo Code

```

J=0
Get_versions
Write filenames to tlog_name array
TIME_A = nexus_run_start

For (k = 1, no_of_versions, k++)
{
    open tlog_names [k]
    read header info
    extract dates and times
    convert dates and times
    work out interval

    if (tlog times within run times)
    {
        Work out difference between TIME_A and tlog start

        For l = 1, reads, l++
        {
            j++
            if (diff > 6) then
            {
                if (j == 1) then
                {
                    tlog_placeholder [j] = 0
                    tlog_time [j] = (tlog_time - run_start_time) - 1 second
                    j++
                    tlog_placholder [j] = data [l]
                    tlog_time [j] = tlog_time [j-1] + 1 second
                }
            }
            else

```

```

{
    tlog_placeholder [j] = 0
    tlog_time [j] = (tlog_time [j - 1] + 1 second)
    j++
    tlog_placeholder [j] = 0
    tlog_time [j] = (tlog_start_time - run_start_time) - 1 second
    j++
    tlog_placeholder [j] = data [I]
    tlog_time [j] = tlog_time [j - 1] + 1 second

}end if

diff = 0

}
else
{
    if ( I == 1 ) then
    {
        tlog_placeholder [j] = data [I]
        tlog_time [j] = tlog_start_time - run_start_time
    }
    else
    {
        tlog_placeholder [j] = data [I]
        tlog_time [j] = tlog_time [j - 1] + interval_time
    }end if
}end if


}end for

TIME_A = tlog_end time

}end if

if ( at least one file has been converted ) then
{
    j++
    tlog_placeholder (j) = 0
    tlog_time (j) = tlog_time (j-1) + 1 second
    j++
    tlog_placeholder (j) = 0
    tlog_time (j) = run_stop_time - run_start_time
}end if
}end for

```



TLOG - from: 11-JUL-2000 19:17:52 to: 11-JUL-2000 19:32:35

159 5.0000

882.0000

5.0000

0.000000

(7(1X,F10.3)) (7(1X,F10.3))

4.999	4.999	4.999	4.999	4.999	4.999	5.000
5.001	5.001	5.001	5.001	5.001	5.001	5.001
5.000	5.000	4.999	4.999	4.999	5.000	5.000
5.000	5.000	4.998	4.998	4.998	4.999	5.000
5.000	5.000	5.000	5.000	5.001	5.001	5.001
5.001	5.001	5.000	5.000	4.999	4.998	4.998
4.998	4.998	4.998	4.998	4.998	4.998	4.999
4.999	4.999	4.999	5.000	5.000	5.000	5.001
5.001	5.001	5.001	5.001	5.001	5.001	5.001
5.001	5.000	5.000	5.000	5.000	5.000	4.998
4.998	5.000	4.999	4.999	4.999	4.999	4.999
4.999	4.999	4.999	4.999	4.999	4.999	4.998
4.998	4.999	5.001	5.001	5.001	5.001	5.001
5.001	5.001	5.001	5.001	5.001	5.001	4.999
4.999	4.999	4.999	4.999	4.999	4.999	4.999
4.999	4.999	4.999	4.999	4.999	4.999	4.999
4.999	5.000	5.000	5.000	5.000	5.000	5.000
5.001	5.001	5.001	5.001	5.000	5.000	5.000
5.000	5.000	5.000	4.999	4.998	4.998	4.998
4.999	4.999	5.000	5.000	5.000	5.000	5.000
5.000	5.000	5.000	5.000	5.000	5.000	5.001
5.001	5.002	5.002	5.001	5.001	5.001	5.001
5.000	4.999	4.999	4.999	4.999		

Buffer: R30000.TLOG | Write | Insert | Forward

Nothing to insert. Use Select to select a range of text.

Figure 4.7 – TLOG File Screenshot

4.4.2 Notes

This is a very complicated process. For legacy reasons, the temperature logs for a run frequently extend over several versions. This process involves identifying those versions that lie within the period of the run and then assembling the separate versions into a single temperature log representing the entire run period. Temperature time stamps are also recalculated to enhance their accuracy and readings outside the run period are discarded.

4.5 UIF Handler

The MACQ handler module consists of a single method which takes a string value representing the filename of a particular User Information File, and returns an integer representing the status of the method i.e. 1 on success or 0 on failure. Users can create User Information files using the UIF_Creator application. The method reads the UIF, allocating appropriate space and stores the information in the global variables that appear below.

Run_array	int [] []	Uif_item_value	String
Item_array	String [] []	Uif_data_name	String
Int read_uif (String filename)			

Figure 4.8 – UIF Handler

4.5.1 UIF Handler Pseudo Code

```

Open file
Read number of elements in user information file
Read number of run numbers in UIF

Allocate run array (sizeof number of run numbers,2)
Allocate item array (sizeof number of elements,2)

Counter=0

For l=1; number of elements ; l++
  Read item run count
  For j=1; item run count; j++
    Counter ++
    Read run number into 'run array' at pos (counter, 1)
    Run array (counter,2) = 1
  End for

  Read 'element name' into 'item array' at pos (l,1)
  Read 'element value' into 'item array' at pos (l,2)
End for

Close file

```

4.5.2 Format of UIF

```

2(number of elements) 6(number of run numbers)
2(numbers of run numbers within particular element) (start of first element)
30000(run numbers)
30004
Light_frequency(Nexus element name)
2.4MHz(element value)
4(number of run numbers within this element) (start of second element)
30002(run numbers)
30010

```


30102
 47890
 $RF^{(Nexus\ element\ name)}$
 $5MHz^{(element\ value)}$

NB. Text in brackets does not appear in file

4.5.3 Implementation

Array of items:

AttName(30char)	AttrValue(30char)
[1] Light_frequency	2.4 MHz
[2] RF	5 MHz

Array of Run Numbers:

[1] Run	Item
30000	1
30004	1
30002	2
30010	2
30102	2
47890	2




Figure 4.9 – UIF File Screenshot

4.5.4 Notes

The process of reading the User Information File is somewhat complex. This is because individual elements and their values have to be linked with a set of run numbers. This has been done by creating a 2-dimensional element array which contains the name of the element to be stored in the nexus file and its value. Also created is a 2-dimensional array to store the run numbers read from the UIF. Each run number is linked a particular element in the 'element array'. This is done by storing the element number in the second dimension of the run array.

4.6 ISF Handler

The ISF Handler encompasses a wide range of functionality, including:

- Instrument Specific Information (i.e. detector angles)
- Time Zero Information
- Dead-time Information
- Grouping Information

Each of these values are obtained from a different source, and separate methods have been produced to extract and process each value.




Figure 4.10 - ISF Handler Class Diagram

4.6.1 Instrument Specific File

If the Instrument Specific File is supplied via the command line

“*convert_nexus/isf=musr_specific.isf*” then the “*get_command_line*” function checks to see whether the file exists. If not, a warning message is displayed.

If an ISF file is not supplied then an appropriate ISF is chosen according to the variable “*NXM_instr_name*”. These instrument Specific files are named appropriately:

- Musr_specific.isf
- Emu_specific.isf
- Mut_specific.isf

If a suitable ISF file cannot be found, or an error occurs while reading the file, then the “detector_angles” are set to 0, and the available attribute of this element will be set to 0 (not

available). If the wrong ISF is supplied for an experiment, then a warning is displayed e.g. if an “emu_specific” ISF file is used for a MuSR experiment.

The coordinate system is set to “spherical”

4.6.1.1 Instrument Specific File Pseudo Code

Open ISF
Check whether ISF matches Instrument used to perform experiment
Read header information
Allocate space for detector angles
Read detector angles
Close ISF

muSr_specific					
32					
detector_no	x	y	z	p	
1	1.0	2.0	3.0	4.0	
2	5.0	6.0	7.0	8.0	
3	9.0	10.0	11.0	12.0	
4	13.0	14.0	15.0	16.0	
5	17.0	18.0	19.0	20.0	
6	21.0	22.0	23.0	24.0	
7	25.0	26.0	27.0	28.0	
8	29.0	30.0	31.0	32.0	
9	33.0	34.0	35.0	36.0	
10	37.0	38.0	39.0	40.0	
11	41.0	42.0	43.0	44.0	
12	45.0	46.0	47.0	48.0	
13	49.0	50.0	51.0	52.0	
14	53.0	54.0	55.0	56.0	
15	57.0	58.0	59.0	60.0	
16	61.0	62.0	63.0	64.0	
17	65.0	66.0	67.0	68.0	
18	69.0	70.0	71.0	72.0	
19	73.0	74.0	75.0	76.0	
20	77.0	78.0	79.0	80.0	
21	81.0	82.0	83.0	84.0	
22	85.0	86.0	87.0	88.0	
23	89.0	90.0	91.0	92.0	
24	93.0	94.0	95.0	96.0	
25	97.0	98.0	99.0	100.0	
26	101.0	102.0	103.0	104.0	
27	105.0	106.0	107.0	108.0	
28	109.0	110.0	111.0	112.0	
29	113.0	114.0	115.0	116.0	
30	117.0	118.0	119.0	120.0	
31	121.0	122.0	123.0	124.0	
32	125.0	126.0	127.0	128.0	
[End of file]					

Figure 4.11 – Instrument Specific File Screenshot

4.6.2 Deadtimes

If the deadtime file is supplied on the command line “convert_nexus/deadtime=dtpar.dat” then Convert_NeXus immediately checks to see whether the included file exists. If not, an error message is displayed and Convert_NeXus continues execution filling deadtimes with 0’s and setting the available attribute to 0.

If the dead-time file is not supplied, the program searches for a specific file in the current directory according to chosen beamline (i.e. “*dtepar.dat*” for EMU and “*dtpar.dat*” for others). If this file is not located then the dead-time array is filled with 0’s and the available attribute is set to 0.

Once a file has been located, the ‘*read_deadtimes*’ function is called from the ‘*isf_handler*’ module. This routine attempts to read the dead-time file. The routine returns 0 if unsuccessful or 1 if successful. If successful the dead-time values are written to the nexus file and the available attribute is set to 1.

The units are set to “*microseconds*”.

4.6.2.1 Dead-Times Pseudo Code

```

Open dead-time file
Read header information
Allocate space for detector dead-times i.e. no of detectors
Read dead-times
Close dead-time file

```

```

DECTerm 1
File Edit Commands Options Print Help
32
1 1.2233107E-02
2 1.3292350E-02
3 3.7105435E-03
4 4.5303139E-03
5 1.0547889E-02
6 1.1666382E-02
7 -1.4865699E-03
8 1.0790898E-02
9 1.1551796E-02
10 1.0585661E-02
11 6.3674781E-03
12 1.0772804E-02
13 4.2096488E-03
14 6.5120654E-03
15 1.2043515E-02
16 8.758552E-03
17 1.0470682E-02
18 8.0206711E-03
19 1.1466230E-02
20 3.9674528E-03
21 1.0470150E-02
22 1.5971461E-02
23 1.1935869E-02
24 6.8573016E-03
25 1.7256520E-03
26 7.8046639E-03
27 1.1020150E-02
28 1.3035779E-02
29 1.3754920E-02
30 8.1274696E-03
31 1.5183973E-02
32 1.2798513E-02
'Run 32705 11-JAN-00'
[End of file]

```

Figure 4.12 – Dead-Time File Screenshot

4.6.3 Groupings

If the grouping file is supplied via the command line “*convert_nexus/grouping=long.uda*” then the “*get_command_line*” function checks to see whether this file exists. If not an error message is displayed.

If the grouping file is not supplied then a suitable grouping file is chosen depending on the magnetic field orientation of the experiment (“*NXM_det_orientation*”); i.e. “*long.uda*” for longitudinal experiments and “*trans.uda*” for transverse experiments. These are searched for in the current directory.

If a grouping file cannot be found then the grouping values are set to 0’s and the available attribute is initialised to 0. This is also true if there is an error encountered while reading a grouping file.

If a longitudinal grouping file is used for a transverse experiment file, or visa versa, then a warning is displayed.

4.6.3.1 Grouping Pseudo Code

```
Open grouping file
Read header information
Allocate grouping array (size of (no of detectors))
If transverse field
{
    read grouping information
}

if longitudinal field
{
    read grouping information
    read alpha value
}

Close grouping file
```




Figure 4.13 – Grouping File Screenshot

4.6.4 Time Zero

When Convert_NeXus executes it searches for a file called “*basetime.uda*” in the current directory; this file contains a T0 value. If this file is available Convert_NeXus includes this value within the nexus file.

If this file is not available then the “*time_zero available*” attribute and the “*time_zero*” value are set to 0.

It is possible to over-ride this process, and include a value through the command line. If a valid floating point number is supplied then Convert_NeXus stores this value and ignores the file “*basetime.uda*”. If the supplied T0 is invalid a warning is supplied to the user and the file is converted without including any T0 value.

The units are set to “*microseconds*”.

4.6.4.1 Time Zero Pseudo Code

```

Open time zero information
Read header information
Read time zero information
Close time zero file

```





Figure 4.14 – Time Zero Screenshot

4.7 Command Line Handler

Virtually all the functionality of Convert_NeXus can be controlled by the command line interface. The Command Line Handler module takes care of the command line values supplied and sets the appropriate flags variables and error messages. This module is quite complex in that it involves many parameters, but consists mainly of standard Fortran string handling.




Figure 4.15 – Command Line Handler Class Diagram

4.7.1 Command Line Handler Pseudo Code

Extract mode
Extract user information file filename
Extract time zero value
Extract grouping file filename
Extract collimator aperture value
Extract collimator type value
Extract beamline
Extract filename/file range

If mode supplied → Ensure mode lies between 0 and 2
If beam-line supplied → Ensure beam-line value lies between 1 and 5
If first-file supplied → ensure value is an integer
If last-file supplied → ensure value is an integer
If time zero value supplied → ensure floating point number

If uif filename supplied → ensure file exists
If dead-time filename supplied → ensure file exists
If isf filename supplied → ensure file exists
If grouping file supplied → ensure file exists

If verbose mode == 0 → ensure beam-line & first-file and last-file || filename supplied

4.7.2 Note

Command line options and parameters supplied by the user are parsed by the `get_command_line()` method. This method sets the appropriate flags and variable values. If it encounters any errors along the way, the user is informed.

When running `Convert_NeXus` the user must choose a beam-line and the appropriate files to process. The methods to obtain these values reside in this module. A few other methods to validate run numbers and floating point numbers also exist in this module.

4.7.3 Collimator

Collimator information is not currently stored at ISIS. It is possible for the users to specify the “Collimator Type” and the “Collimator_aperture” when converting files. This is done by specifying command line parameters. These values are stored as strings and have a maximum length of 30 characters.

If, however, the user does not supply a value, the “collimator type” is set to “slits” and the “collimator aperture” is set to “0”. Because the values are stored as strings no validation takes place.

4.8 Histogram Handler

A module for handling histogram information is necessary because Convert_NeXus undertakes a lot of processing in splitting the histograms into groups with common characteristics.




Figure 4.16 – Histogram Handler Class Diagram

4.8.1 Get Data Groups Pseudo Code

```
Create hist_reference_list
For (l=1; l<= no of histograms; l++)
{
    Find 0 entry in hist_reference_list
    For (j=1; j<=no of histograms; j++)
    {
        If histogram properties (j) == entry
            Add 1 in hist_reference_list at entry j
    }
}
```

4.8.2 Split Data Groups Pseudo Code

```
Work out how many histograms in current group
Allocate space for data "temp_data"
For (l=1; l<=no of histograms; l++)
{
    if (hist_reference_list (i) == current group)
    {
        read histogram data into temp_data array
    }
}
```

4.8.3 Notes

The Get Data Groups method essentially creates a list with an entry for each histogram. Each entry contains a number which represents the group that the histogram belongs to. If all the histograms are the same, then they will all have the value 1.

The Split Data Groups method actually creates a temporary histogram data array (the size of the number of histograms in the current group x histogram length). It then reads the appropriate information into this array. Convert_NeXus is then free to directly include this data when processing the histogram_data group.

4.9 ISOtime

This module was written specifically for use by Convert_NeXus and is not meant to offer a definitive set ISO date and time functions.

Structures

Name : *ISOtime_struct*
entities : *hours* (int)
 : *minutes* (int)
 : *seconds* (int)

Name : *ISOdate_struct*
entities : *years* (int)
 : *months* (int)
 : *days* (int)

Operators +, -, .eq., .gt., .ge., .lt., .le.

These operators have been overloaded so that date structures can be compared lexically with each other and time structures can be compared lexically with each other.

Also dates can be subtracted from each other. This function is used to determine the difference between two times and will therefore give unexpected results if the first argument is less than the second. This is the same for the time structures.

Times can be added together. Although dates cannot.

Date functions

logical date_greater_than (ISOdate_struct, ISOdate_struct)
//logical function to determine if first date is greater than second date

logical date_greater_than_or_equal (ISOdate_struct, ISOdate_struct)
//logical function to determine if first date is greater than or equal to the second date

logical date_less_than (ISOdate_struct, ISOdate_struct)
//logical function to determine if first date is less than second date

logical date_less_than_or_equal (ISOdate_struct, ISOdate_struct)
//logical function to determine if first date is less than or equal to the second date

logical date_equal_to (ISOdate_struct, ISOdate_struct)
//logical function to determine if first date is equal to the second date

ISOdate_struct subtract_date (ISOdate_struct, ISOdate_struct)
//function to subtract second ISO date from the first, returning an ISOdate structure with the result

character(len=10) date_to_string (ISOdate_struct)
//function to return value of ISOdate in string format

ISOdate convertDate (character(len=9))
//function to convert a string date in this format "22-FEB-99" to ISOdate format

Time Functions

```

logical time_greater_than (ISotime_struct, ISotime_struct)
//logical function to determin if first time is greater than second time

logical time_greater_than_or_equal (ISotime_struct, ISotime_struct)
//logical function to determin if first time is greater than or equal to the second time

logical time_less_than (ISotime_struct, ISotime_struct)
//logical function to determin if first time is less than second time

logical time_less_than_or_equal (ISotime_struct, ISotime_struct)
//logical function to determin if first time is less than or equal to the second time

logical time_equal_to (ISotime_struct, ISotime_struct)
//logical function to determin if first time is equal to the second time

ISotime_struct convertTime (character(len=8))
//function which takes a string time in the format "22:45:12" and converts it into an ISotime structure

character(len=10) time_to_string (ISotime_struct)
//function to output time structure in string format

ISotime_struct add_ISotime (ISotime_struct, ISotime_struct)
//function which adds two time structures together and outputs a time structure with the result

ISotime_struct subtract_ISotime (ISotime_struct, ISotime_struct)
//function which subtracts second ISotime from first and outputs result in a time structure

ISotime_struct add_hours (ISotime_struct, int)
//function adds hours to time structure and outputs result in form of ISotime_struct

ISotime_struct add_minutes (ISotime_struct, int)
//function adds minutes to time structure and outputs result in form of ISotime_struct

ISotime_struct add_seconds (ISotime_struct, int)
//function adds seconds to time structure and outputs result in form of ISotime_struct

ISotime_struct subtract_hours (ISotime_struct, int)
//function subtracts hours from time structure and outputs result in form of ISotime_struct

ISotime_struct subtract_minutes (ISotime_struct, int)
//function subtracts minutes from time structure and outputs result in form of ISotime_struct

ISotime_struct subtract_seconds (ISotime_struct, int)
//function subtracts seconds from time structure and outputs result in form of ISotime_struct


```

```

miscellaneous functions

ISotime_struct build_time (int)
//functions which builds time structure when passed time in seconds

character(len=20) toString (int)
//function to convert in to string. Adds '0' to front if int is from 0..9

```


4.10 TESTING PROCEDURE FOR CONVERT_NEXUS

The functionality of various aspects of Convert_NeXus has been tested by evaluating the response of the program against input data files with known properties. The results are summarised in the following table:

4.10.1 MCS File

Test no	Start run	Stop run	Beam line	Comment	Command
1	-	-	5	File r31111.ral is corrupt. Convert_NeXus should discard this file from processing.	Convert_NeXus / /beamline=5 /tst=r31111.ral
2	-	-	5	File emu30000.ral has Red/green data. Convert_NeXus should convert this file easily	Convert_NeXus /beamline=5 /tst=emu30000.raw
3	-	-	5	File mut00477.raw has differing histogram properties. Convert_NeXus should separate these differing histograms into different NXdata groups	Convert_NeXus /beamline=5 /tst=mut00477.raw

4.10.2 Temperature_logs

Test No	Start run	Stop run	Beam line	Comment	Command
4	32482	32482	2	Run lasts over 3 days, tlog 10 versions long	Convert_NeXus /beamline=2 /firstfile=32482 /lastfile=32482
5	37500	37500	2	Normal Tlog	Convert_NeXus /beamline=2 /firstfile=37500 /lastfile=37500

4.10.3 MACQ files

Test No	Start run	Stop run	Beam line	Comment	Command
6	35074	35074	1	MACQ log file has only one value so not included, tlog doesn't exist, warning given	Convert_NeXus /beamline=1 /firstfile=35074 /lastfile=35074
7	37500	37500	2	Normal MACQ log	Convert_NeXus /beamline=2 /firstfile=37500 /lastfile=37500

4.10.4 Grouping file

Test No	Start run	Stop run	Beam line	Comment	Command
8	35000	35000	1	Specify a grouping file that doesn't exist. Should display warning message and continue conversion	Convert_NeXus /beamline=1 /firstfile=35000 /lastfile=35000 /grouping=xyz.xyz
9	37000	37000		Specify in effect a corrupt grouping file	Convert_NeXus /firstfile=37000 /lastfile=37000 /beamline=1 /grouping=emu30000.raw
10	35600	35600	2	Specify transverse grouping file for longitudinal experiment. Should convert but display a warning	Convert_NeXus /firstfile=35600 /lastfile=35600 /beamline=2 /grouping=trans.uda
11	35600	35600	2	Specify longitudinal grouping for same file. Should convert no problem	Convert_NeXus /firstfile=35600 /lastfile=35600 /beamline=2 /grouping=long.uda

4.10.5 Deadtime file

Test No	Start run	Stop run	Beam line	Comment	Command
12	35600	35600	2	Specify deadtime file that does not exist, Should display warning and continue conversion	Convert_NeXus /beamline=1 /firstfile=35600 /lastfile=35600 /deadtime=xyz.xyz
13	35600	35600	2	Specify in effect corrupt deadtime file. Should display warning and continue conversion	Convert_NeXus /Beamline=1 /firstfile=35600 /lastfile=35600 /deadtime=test_uif.uif
14	35600	35600	2	Specify valid deadtime file. Should process accordingly	Convert_NeXus /Beamline=1 /firstfile=35600 /lastfile=35600 /deadtime=dtpar.dat

4.10.6 User Information File

Test No	Start run	Stop run	Beam line	Comment	Command
15	30000	30003	2	UIF has elements referencing run 30000, 30001, 30002. Should include appropriate elements	Convert_NeXus /beamline=2 /firstfile=30000 /lastfile=30003 /uif=test_uif.uif
16	30000	30003	2	Specify non existent UIF, should display a warning and continue	Convert_NeXus /beamline=2 /firstfile=30000 /lastfile=30003 /uif=xyz.xyz
17	30000	30003	2	Specify corrupt uif. Should display a warning and continue conversion	Convert_NeXus /beamline=2 /firstfile=30000 /lastfile=30003 /uif=corrupt_uif.uif

4.10.7 Instrument Specific File

Test No	Start run	Stop run	Beam line	Comment	Command
18	35600	35600	2	Specify valid ISF for experiment	Convert_NeXus /firstfile=35600 /lastfile=35600 /beamline=2 /isf=musr_specific.isf
19	37000	37000	1	Specify MUSR ISF for EMU experiment. Should process but display warning	Convert_NeXus /firstfile=37000 /lastfile=37000 /beamline=1 /isf=musr_specific.isf
20	35600	35600	2	Specify non existent ISF. Should display warning and continue conversion	Convert_NeXus /firstfile=35600 /lastfile=35600 /beamline=2 /isf=xyz.xyz
21	35600	35600	2	Specify corrupt ISF. Should display warning and continue conversion	Convert_NeXus /firstfile=35600 /lastfile=35600 /beamline=2 /isf=test_uif.uif

4.10.8 Collimator Information

Test No	Start run	Stop run	Beam line	Comment	Command
22	30000	30000	2	Test to enter collimator information	Convert_NeXus /beamline=2 /firstfile=30000 /lastfile=30000 /col.type=slitties /col.apert=1.7
23	30000	30000	2	Test blank collimator information, should set to default and continue processing	Convert_NeXus /firstfile=30000 /lastfile=30000 /beamline=2 /col.type= /col.apert=

4.10.9 Time Zero Information

Test No	Start run	Stop run	Beam line	Comment	Command
24	30000	30000	2	Test to enter time_zero information	Convert_NeXus /beamline=2 /firstfile=30000 /lastfile=30000 /T0=1.45
25	30000	30000	2	Test with invalid floating point number	Convert_NeXus /firstfile=30000 /lastfile=30000 /beamline=2 /T0=hello

5. NeXus_READER

5.1 Introduction

Complementary to the Convert_NeXus application, NeXus read routines are provided to enable information stored in the NeXus muon data files to be read by user applications. Versions written in Fortran 77 and C are available. These read routines conform to the ISIS Muon Instrument Definition Version 1 and have been tested under a variety of operating systems including VMS and Windows NT.

5.1.1 Common features

Despite the idiosyncrasies of the C and Fortran programming languages, we have attempted to maintain a consistent style for the structure and semantics of the two versions of the ISIS Muon NeXus read routines. Variable names, constants, display messages and the order in which the read routines traverse the NeXus file have been synchronised. The main differences are the way each reader handles memory allocation and the method by which the user can access the desired information.

5.1.2 ISIS Muon NeXus Instrument Definition

The NeXus read routines have to follow a standard. We have produced an Instrument Definition which describes the structure and content of our NeXus file. We are currently at version 1, and the read routines traverse the NeXus files according to this specification. The ISIS Muon Instrument Definition is listed in Chapter 2. This definition denotes the optional and essential elements in accordance with the μ SR community Instrument Definition discussed during the NeXus workshop held at PSI, Switzerland in Spring 2001. Also specified are the units, coordinate systems and data types (including the rank of data types). It is important to note that although we have released version 1 of our Instrument Definition and accompanying software, this is by no means definitive and future modifications are envisaged.

5.2 C Version

The C NeXus read routine has been developed in Microsoft Visual C++ on the Windows platform; the code is ANSI compatible. The routine has been tested under both Windows NT and the Open VMS operating systems.

The C NeXus read routine comprise of the following files:

- Nexus_reader.h – C header file which defines a structure which contains all the Muon data, defines a number of constants and methods which are used in the NeXus reader.
- Nexus_reader.c – contains the code that actually reads the Nexus file.
- NAPI.h – ANSI C NeXus API Header File.

5.2.1 Nexus_reader.h

```

struct NXM_MUONUUF
{
    NXname      uif_element_name;
    NXname      uif_element_value;
};

struct NXM_MUONIDF
{
    /* NXrun */
    NXname      run_program_name;
    NXname      run_program_version;
    int         run_IDF_version;
    int         run_number;
    .
    .
    /* NXuser */
    NXname      user_name;
    NXname      user_experiment_number;
    int         user_uif_array_length;
    struct NXM_MUONUUF * user_uif_array;
    .
    .
    /* NXdata histogram_data_1 */
    int *       histogram_counts;
    int         histogram_number;
    int         histogram_length;
    int         histogram_t0_bin;
    int         histogram_first_good_bin;
    int         histogram_last_good_bin;
    int         histogram_resolution;
    NXname      histogram_resolution_units;
    float *     histogram_raw_time;
    NXname      histogram_raw_time_units;
    float *     histogram_corrected_time;
    NXname      histogram_corrected_time_units;
    int         histogram_grouping_available;
    int *       histogram_grouping;
    int         histogram_alpha_available;
    float *     histogram_alpha;
    .
    .
    /* NXlog events_log_1 */
    int         events_log_available;
    NXname      events_log_name;
    float *     events_log_values;
    NXname      events_log_values_units;
    float *     events_log_time;
    .
    .
    NXname      events_log_time_units;
};

int NXMread(char *nxfname, struct NXM_MUONIDF **nxpdata);
void NXMfreememory(struct NXM_MUONIDF *nxpdata);

```

Here a structure is defined for the UIF elements.

Next a structure is defined to store the contents of the nexus file. This structure may contain many UIF structures.

Note, for brevity, some of the element declarations are missing from this figure.

This code denotes the methods that are available in the NeXus read routine, along with their input and return parameters.

This header file consists of the 'muon_def' structure which will act as a container for the data read from the nexus file. This method has the advantage that the nexus read routine can read the data from the nexus file into this structure and pass the self contained structure neatly back to the user's program. The header file also defines the functions contained within the 'c_nexus_reader.c' file.

5.2.2 Nexus_reader.c

```
#include <stdio.h>
#include "napi.h"
#include "nexus_reader.h"
```

*These are the files
that NeXus_Reader.c
are dependant*

```
/* switch for debugging information */
#define _DEBUG_ 1
```

```
/* information locating read routine */
#define _READVERSION_ 1 /* Version number of read routine */
#define _IDFVERSION_ 1 /* Compatible version of Instrument Definition File */
#define _ANALYSIS_ "muonTD" /* Compatible analysis */
#define _LAB_ "ISIS" /* Compatible lab */
```

*These are constants
used by the read
routine*

```
/* structure defining file information */
struct NXM_MUONIDF muon;
```

*Here is where the
structure that defines
the contents of the
NeXus file is defined*

```
/* macro for adding debug information */
#define NXMdebug(format, value) \
    if (_DEBUG_) printf(format, value)
```

```
/* macro handling function calling, returns with NX_ERROR if function call failed */
#define NXMerrorhandler(function) \
    if (function != NX_OK) return NX_ERROR
```

```
static int NXMdatasize(int dataType)
{
    int size;

    switch (dataType)
    {
        case NX_CHAR:
            size = SIZE_CHAR8;
            break;
        case NX_FLOAT32:
            size = SIZE_FLOAT32;
            break;
        case NX_FLOAT64:
            size = SIZE_FLOAT64;
            break;
        case NX_INT8:
            size = SIZE_INT8;
            break;
        case NX_UINT8:
            size = SIZE_UINT8;
            break;
        case NX_INT16:
            size = SIZE_INT16;
            break;
        case NX_UINT16:
            size = SIZE_UINT16;
            break;
        case NX_INT32:
            size = SIZE_INT32;
            break;
        case NX_UINT32:
            size = SIZE_UINT32;
            break;
        default:
            size = 0;
            break;
    }

    return size;
}
```

*Function returns size
of supplied data type
in bytes. This make
the read routine
scalable across
different platforms*


```

static int NXMgetstringdata(NXhandle file_id, char *str)
{
    int i, rank, type, dims[32];
    NXname data_value;

    NXErrorHandler((NXgetinfo(file_id, &rank, dims, &type)));
    if ((type != NX_CHAR) || (rank > 1) || (dims[0] >= VGNAMELENMAX))
    {
        printf("\nFatal error in routine 'getstringdata'\n\n");
        exit(0);
    }

    NXErrorHandler((NXgetdata(file_id, data_value)));

    for (i = 0; i < dims[0]; i++)
        *(str + i) = *(data_value + i);
    *(str + i) = '\0';

    return 1;
}

```

This method handles the reading of string data from the NeXus file, adding a null terminator for compatibility with C string handling

```

static int NXMgetstringattr(NXhandle file_id, char *name, char *str)
{
    int i, attlen, atttype;
    NXname data_value;

    attlen = VGNAMELENMAX - 1;
    atttype = NX_CHAR;
    NXErrorHandler((NXgetattr(file_id, name, data_value, &attlen, &atttype)));

    for (i = 0; i < attlen; i++)
        *(str + i) = *(data_value + i);
    *(str + i) = '\0';

    return 1;
}

```

This method handles the reading of string data attributes from the NeXus file, adding a null terminator for compatibility with C String handling

```

static void *NXMmemoryhandler(NXhandle file_id)
{
    int i, rank, type, dims[32], size;
    char *data_ptr;

    NXErrorHandler((NXgetinfo(file_id, &rank, dims, &type)));
    size = dims[0];
    for (i = 1; i < rank; i++)
        size = size * dims[i];
    size = size * NXMdatasize(type);
    if ((data_ptr = malloc (size)) == NULL)
    {
        printf("\nMemory allocation error: can't allocate %d bytes\n\n", size);
        exit(0);
    }
}

```

This method allocates memory for a NeXus data item, returning a pointer to the allocated memory block

```

static void NXMfree(void *data_ptr)
{
    if (data_ptr != NULL) free(data_ptr);
}

```

Method frees memory from NeXus data items provided the pointer argument is non NULL

```

int NXMread(char *nxfname, struct NXM_MUONIDF **nxfdata)
{
    NXhandle file_id;
    NXname group_name, class_name, data_name;
    int i, attlen, atttype, numitems, data_type, uif_count;

    *nxfdata = &muon;

    /* open NeXus file for reading */
    NXMerrorhandler((NXopen(nxfname, NXACC_READ, &file_id)));
    NXMdebug("%s\n", "Opening NeXus file...");

    /* open root group 'NXentry' */
    NXMerrorhandler((NXopengroup(file_id, "run", "NXentry")));
    NXMdebug("%s\n", "Opening 'NXentry'...");

    /* display (in debug mode) NeXus reader version and IDF compatible version number */
    NXMdebug("\nNeXus reader version %d\n", _READVERSION_);
    NXMdebug("Compatible with Instrument Definition File version %d\n", _IDFVERSION_);

    /* read IDF version */
    NXMerrorhandler((NXopendata(file_id, "IDF_version")));
    NXMerrorhandler((NXgetdata(file_id, &muon.run_IDF_version)));
    NXMerrorhandler((NXclosedata(file_id)));
    NXMdebug("IDF version: %d\n", muon.run_IDF_version);

    /* check reader compatibility */
    if ((muon.run_IDF_version != _IDFVERSION_) ||
        (strcmp(muon.run_analysis, _ANALYSIS_) || (strcmp(muon.run_lab, _LAB_)))
    {
        printf("\nFatal error: Read routine incompatible with Instrument Definition used
                to write data\n");
        printf("Found, IDF Version: %d, Analysis: %s, Lab: %s\n", muon.run_IDF_version,
                muon.run_analysis, muon.run_lab);
        printf("Expecting: IDF Version: %d, Analysis: %s, Lab: %s\n\n", _IDFVERSION_,
                _ANALYSIS_, _LAB_);
        exit(0);
    }

    /* close root group 'NXentry' */
    NXMerrorhandler((NXclosegroup(file_id)));
    NXMdebug("%s\n", "Close 'NXentry'");

    /* close NeXus file */
    NXMerrorhandler((NXclose (&file_id)));
    NXMdebug("%s\n", "Close NeXus file");

    return NX_OK;
}

```

Returns data from the NeXus file 'nxfname' in structure 'nxfdata'

Here the read routine checks that the NeXus file is compatible with the reader.

The read routine continues to traverse the NeXus file and closes the file on completion.

```

/* NXMfreememory - free memory claimed in NeXus file */
void NXMfreememory(struct NXM_MUONIDF *nxfdata)
{
    /* NXlog (temperature) */
    NXMfree(nxfdata->temperature_log_values);
    NXMfree(nxfdata->temperature_log_time);
    /* NXlog (events) */
    NXMfree(nxfdata->events_log_values);
    NXMfree(nxfdata->events_log_time);
    /* NXdata */
    NXMfree(nxfdata->histogram_counts);
    NXMfree(nxfdata->histogram_raw_time);
    NXMfree(nxfdata->histogram_corrected_time);
    NXMfree(nxfdata->histogram_grouping);
    NXMfree(nxfdata->histogram_alpha);
    /* NXuser */
    NXMfree(nxfdata->user_uif_array);
    /* NXsample */
    NXMfree(nxfdata->sample_magnetic_field_vector);
    /* NXdetector */
    NXMfree(nxfdata->detector_angles);
    NXMfree(nxfdata->detector_deadtimes);
}

```

This method frees the memory of each of the allocatable data items by calling NXMfree() on each of them in turn.

5.2.3 Test_nexus_reader.c

```
#include "napi.h"
#include "nexus_reader.h"
```

These are the files that the user must include in their own program.

```
void main()
{
    struct NXM_MUONIDF *muon_ptr;
    int *histogram_ptr;
    int i, j, histogram;
    char nxfname[80];

    /* read filename from keyboard */
    printf("Please enter nexus filename : ");
    scanf("%s", nxfname);
```

Here the user must declare a pointer that will reference the NXM_MUONIDF structure and a pointer that will reference the histogram data within that structure.

```
/* read NeXus file */
(void) NXMread(nxfname, &muon_ptr);
```

To read the NeXus file, the user must pass the NXMread routine a filename and a pointer to a NXM_MUONIDF structure.

```
/* print contents of NeXus file */
printf("\n\nNeXus File Contents\n\n");
printf("Program name: %s\n", muon_ptr->run_program_name);
printf("Program version: %s\n", muon_ptr->run_program_version);
.
.
.
/* NXdata */
printf("Number of histograms: %d\n", muon_ptr->histogram_number);
printf("Time zero bin: %d\n", muon_ptr->histogram_t0_bin);
printf("Histogram length: %d\n", muon_ptr->histogram_length);
printf("First good bin: %d\n", muon_ptr->histogram_first_good_bin);
printf("Last good bin: %d\n", muon_ptr->histogram_last_good_bin);
printf("Histogram resolution: %d\n", muon_ptr->histogram_resolution);
.
.
.
printf("\nHistogram information - (Raw Time, Corrected Time, Counts)
(time in %s):\n", muon_ptr->histogram_corrected_time_units);
histogram = 0; /* choose which histogram you want (starting at zero) */
histogram_ptr = &(muon_ptr->histogram_counts[histogram * muon_ptr->
```

The supplied pointer now references the NeXus data and the user can access the information as left.

```
histogram_length]);

* get the start of desired histogram */
for (i = 0; i < (muon_ptr->histogram_length); i+=50) /* extract data */
{
    printf("(%f, %f, %d)\n", *(muon_ptr->histogram_raw_time + i),
        *(muon_ptr->histogram_corrected_time + i), *(histogram_ptr+=50));
}
printf("\n");
.
.
.
/* access logged information ... events */
if (muon_ptr->events_log_available > 0)
```

This code demonstrates pointer manipulation of individual histogram by printing every 50th value in the 1st histogram with time stamp.

```
{
    printf("\nEvent Log - Time (%s), Recorded events (%s)\n", muon_ptr->
        events_log_time_units, muon_ptr->events_log_values_units);
    for (i = 0; i < (muon_ptr->events_log_available); i++)
    {
        printf("(%f, %f)\n", muon_ptr->events_log_time[i], muon_ptr->
            events_log_values[i]);
    }
}
else
    printf("\nEvent log not available\n");
```

This sample checks to see if any event log values are available, if so it prints them to the screen.

```
/* free memory claimed by read routine */
NXMfreememory(muon_ptr);
}
```

To prevent memory leaks the user should call the NXMfreememory routine to reclaim the memory occupied by the NeXus data.

5.2.4 Compilation Instructions


To recompile this application, the following files will be required along with a standard ANSI C compiler:

- NAPI.h
- Nexus_reader.c
- Nexus_reader.h
- Test_nexus_reader.c (or suitably modified for a user application)

In these instructions it is assumed that the HDF libraries and NeXus API's are currently installed. If these are not yet installed, please refer to the NeXus Data Format Homepage (http://www.neutron.anl.gov/nexus/NeXus_API.html#install) for detailed instructions.

Step 1 - Compile the above ".c" files in sequence

Step 2 - Link Test_nexus_reader, Nexus_reader, NeXus API (where NeXus API exists)

A screenshot of a DECterm window titled "DECterm 1". The window has a menu bar with "File", "Edit", "Commands", "Options", "Print", and "Help". The main text area contains the following commands:

```
$set verify
$define hdf_root axplib$disk:[4_1r3_alphavms7_1.] ! where you installed HDF, trailing "." is needed
$define nexus_root user$disk:[df01.] ! where the nexus library exists
$ CC/INCLUDE=axplib$disk:[4_1r3_alphavms7_1.include] NEXUS-READER.C
$ CC/INCLUDE=HDF_ROOT:[INCLUDE] test_nexus_reader.c
$ LINK test_nexus_reader,nexus_reader,nexus_root:[nexus]NEXUS.OLB/LIB,HDF_ROOT:[LIB]MFHDF/LIB,DF/LIB,LIBZ/LIB,LIBJPEG/LIB
```

Figure 5.1 - NeXus_Reader Compilation Instructions

5.2.5 Limitations

The C NeXus read routine is restricted in that it currently only reads one NXdata group. At present, however, this is acceptable for the ISIS muon data files since it is very rare that an experiment has differing histogram properties. Also, the reader only reads one set of temperature log data and one set of event log data.

5.3 F77 Reader

The F77 NeXus read-routine has been developed on the VMS operating system using a standard DEC Fortran 77 compiler. We have decided to produce a Fortran 77 version because this will be compatible both with Fortran 77 and Fortran 90 applications. The read routine comprises of the following files:

- NAPIF.INC
- NEXUS_READER.FOR
- MUON_DEF_F77.INC or MUON_DEF_F90.INC (for either Fortran 77 or Fortran 90 respectively)

5.3.1 MUON_DEF

There are two versions of the MUON_DEF: MUON_DEF_F77.INC and MUON_DEF_F90.INC. This is purely because of the syntactic peculiarities of the F77 and F90 programming languages. The user must choose a version according to the language their application is written in. The MUON_DEF essentially declares a common block for character data, integer data and float data, defining the variables that will be read from the NeXus file. These variables will then be accessible from the users application.

```
C      constant to define length of strings
      INTEGER NXMname
      PARAMETER (NXMname = 60)

C      switch for debugging information
      LOGICAL DEBUG
      PARAMETER (DEBUG = .TRUE.)

C      version number of the read routine
      INTEGER READVERSION
      PARAMETER (READVERSION = 1)

C      compatible version of instrument definition file
      INTEGER IDFVERSION
      PARAMETER (IDFVERSION = 1)

C      compatible analysis
      CHARACTER*(NXMname) ANALYSIS
      PARAMETER (ANALYSIS = 'muonTD')

C      compatible lab
      CHARACTER*(NXMname) LAB
      PARAMETER (LAB = 'ISIS')

C      uif array length
      INTEGER UIFLENGTH
      PARAMETER (UIFLENGTH = 30)

C      number of detectors
      INTEGER NUMDETECTORS
      PARAMETER (NUMDETECTORS = 32)
```

It is necessary to define certain constants for both versions of the NeXus Reader, but the Fortran 77's inability to support dynamic memory allocation means that we must declare the upper boundaries for each array.

```

C ##### NXrun #####
CHARACTER*(NXMname) NXM_run_program_name
CHARACTER*(NXMname) NXM_run_program_version
INTEGER NXM_run_idf_version
INTEGER NXM_run_number
CHARACTER*(NXMname) NXM_run_title

```

Here each of the variables that will be utilised by the users program is defined.

Each variable has been named carefully to give a self explanatory description of what each element holds.

```

C ##### NXdata histogram_data_1 #####
INTEGER NXM_histogram_counts (MAXHISNUM, MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_counts_units
INTEGER NXM_histogram_number
INTEGER NXM_histogram_length
INTEGER NXM_histogram_t0_bin
INTEGER NXM_histogram_first_good_bin
INTEGER NXM_histogram_last_good_bin
INTEGER NXM_histogram_resolution
CHARACTER*(NXMname) NXM_histogram_resolution_units

```

Once all the variables have been defined, they must be included in a common block.

The common block is used because the reader extracts the data from the NeXus file and places it in memory, and by defining this common block in the MUON_DEF header file, the users application can access this same area of memory.

Note: To avoid confusion, a common block has been declared for each data type, i.e. integer, real and character. This reduces the likelihood of memory overwrite.

```

C ##### NXlog event_log_1 #####
INTEGER NXM_event_log_available
CHARACTER*(NXMname) NXM_event_log_name
REAL NXM_event_log_values (MAXLOGLEN)
CHARACTER*(NXMname) NXM_event_log_values_units
REAL NXM_event_log_time (MAXLOGLEN)
CHARACTER*(NXMname) NXM_event_log_time_units

```

```

C ----- define common blocks here -----

```

```

C block for integer data

```

```

COMMON /NXM_integer_data_block/ NXM_run_idf_version, NXM_run_number,
NXM_run_switching_states, NXM_user_uif_array_length,
&NXM_sample_mfield_vec_available, NXM_detector_number,
NXM_detector_angles_available, NXM_deadtims_available,
&NXM_beam_daereads, NXM_beam_frames, NXM_histogram_counts, NXM_histogram_number,
NXM_histogram_length,
&NXM_histogram_t0_bin, NXM_histogram_first_good_bin, NXM_histogram_last_good_bin,
NXM_histogram_resolution,
&NXM_histogram_grp_available, NXM_histogram_grouping,
NXM_histogram_alpha_available, NXM_histogram_t0_available,
&NXM_temp_log_available, NXM_event_log_available

```

```

C block for real data

```

```

COMMON /NXM_real_data_block/ NXM_sample_temperature, NXM_sample_magnetic_field,
NXM_sample_mfield_vector,
&NXM_detector_angles, NXM_detector_deadtims, NXM_beam_total_counts,
NXM_histogram_raw_time,
&NXM_histogram_corrected_time, NXM_histogram_alpha, NXM_histogram_time_zero,
NXM_temp_log_values, NXM_temp_log_time,
&NXM_event_log_values, NXM_event_log_time

```

```

C block for character data

```

```

COMMON /NXM_character_data_block/ NXM_run_program_name, &NXM_run_program_version,
&NXM_run_title, NXM_run_notes, NXM_run_analysis, NXM_run_lab, NXM_run_beamline,
&NXM_run_start_time, NXM_run_stop_time, NXM_run_duration, NXM_user_name,
&NXM_user_experiment_number, NXM_user_uif_element_name, NXM_user_uif_element_value,
&NXM_sample_name, NXM_sample_temperature_units, NXM_sample_magnetic_field_units,
&NXM_sample_shape, NXM_sample_magnetic_field_state, NXM_sample_mfield_vector_coord,
&NXM_sample_environment, NXM_instrument_name, NXM_detector_orientation,
&NXM_detector_angles_coord, NXM_detector_deadtims_units, NXM_collimator_type,
&NXM_collimator_aperture, NXM_beam_total_counts_units, NXM_histogram_counts_units,
&NXM_histogram_resolution_units, NXM_histogram_raw_time_units,
&NXM_histogram_ctime_units, NXM_histogram_time_zero_units, NXM_temp_log_name,
&NXM_temp_log_values_units, NXM_temp_log_time_units, NXM_event_log_name,
&NXM_event_log_values_units, NXM_event_log_time_units

```

5.3.2 NEXUS_READER.FOR

The F77 NeXus consists of one subroutine “NXMread” which traverses the NeXus file and places its contents in memory (via common blocks). All array sizes are declared statically using predefined constants.

```

INTEGER FUNCTION NXMread (nxfname)
    include 'NAPIF.INC'

C    constant to define length of strings
    INTEGER NXMname
    PARAMETER (NXMname = 60)

C    switch for debugging information
    LOGICAL DEBUG
    PARAMETER (DEBUG = .FALSE.)

C    version number of the read routine
    INTEGER READVERSION
    PARAMETER (READVERSION = 1)

C    compatible version of instrument definition file
    INTEGER IDVERSION
    PARAMETER (IDVERSION = 1)

C    maximum number of histograms
    INTEGER MAXHISNUM
    PARAMETER (MAXHISNUM = 80)

.
.
C    ----- define NeXus elements -----
C    ##### NXrun #####
    CHARACTER*(NXMname) NXM_run_program_name
    CHARACTER*(NXMname) NXM_run_program_version
    INTEGER NXM_run_idf_version
    INTEGER NXM_run_number
    CHARACTER*(NXMname) NXM_run_title
    CHARACTER*(NXMname) NXM_run_notes
    CHARACTER*(NXMname) NXM_run_analysis
    CHARACTER*(NXMname) NXM_run_lab
    CHARACTER*(NXMname) NXM_run_beamline

.
.
C    ##### NXdata histogram_data_1 #####
    INTEGER NXM_histogram_counts (MAXHISNUM, MAXHISLEN)
    CHARACTER*(NXMname) NXM_histogram_counts_units
    INTEGER NXM_histogram_number
    INTEGER NXM_histogram_length
    INTEGER NXM_histogram_t0_bin
    INTEGER NXM_histogram_first_good_bin
    INTEGER NXM_histogram_last_good_bin
    INTEGER NXM_histogram_resolution
    CHARACTER*(NXMname) NXM_histogram_resolution_units
    REAL NXM_histogram_raw_time (MAXHISLEN)

.
.
C    ##### NXlog temperature_log_1 #####
    INTEGER NXM_temp_log_available
    CHARACTER*(NXMname) NXM_temp_log_name
    REAL NXM_temp_log_values (MAXLOGLEN)
    CHARACTER*(NXMname) NXM_temp_log_values_units
    REAL NXM_temp_log_time (MAXLOGLEN)
    CHARACTER*(NXMname) NXM_temp_log_time_units

```

The reader takes as an input parameter a string representing the NeXus filename.

The Reader also defines a number of constants. As mentioned before, the F77 header must declare upper boundary limits on arrays because we cannot dynamically allocate space in Fortran 77.

We must declare variables to store the data read from the NeXus file. These names have been kept identical to those used in the MUON_DEF.

```

C ----- define common blocks here -----
C block for integer data
COMMON /NXM_integer_data_block/ NXM_run_idf_version, NXM_run_number,
&NXM_run_switching_states, NXM_user_uif_array_length,
&NXM_sample_mfield_vec_available, NXM_detector_number,
&NXM_detector_angles_available, NXM_deadtrees_available,
&NXM_beam_daereads, NXM_beam_frames, NXM_histogram_counts, NXM_histogram_number,
&NXM_histogram_length, NXM_histogram_t0_bin, NXM_histogram_first_good_bin,
&NXM_histogram_last_good_bin, NXM_histogram_resolution,
&NXM_histogram_grp_available, NXM_histogram_grouping,
&NXM_histogram_alpha_available, NXM_histogram_t0_available,
&NXM_temp_log_available, NXM_event_log_available
.
.
C @@@@@@@@@@@@@@@@@@ Open Nexus File @@@@@@@@@@@@@@@@@@
if (NXopen (nxfname, NXACC_READ, file_id) .NE. NX_OK) stop
if (debug) print *, 'Opening NeXus file...'

C @@@@@@@@@@@@@@@@@@ Open 'Run' Group @@@@@@@@@@@@@@@@@@
if (NXopengroup (file_id, 'run', 'NXentry') .NE. NX_OK) stop
if (debug) print *, 'Opening "NXentry"...'

C display (in debug mode) Nexus reader version and IDF compatible
version number
if (debug) then
  print *, ' '
  print *, 'Nexus reader version ', READVERSION
  print *, 'Compatible with ISIS Muon Instrument Definition
& version ', IDFVERSION
  print *, ' '
end if

C read idf version
if (NXopendata (file_id, 'idf_version') .NE. NX_OK) stop
if (NXgetdata (file_id, NXM_run_idf_version) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'IDF version : ', NXM_run_idf_version

C read analysis
if (NXopendata (file_id, 'analysis') .NE. NX_OK) stop
if (NXgetchardata (file_id, NXM_run_analysis) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Analysis : ', NXM_run_analysis

C read lab
if (NXopendata (file_id, 'lab') .NE. NX_OK) stop
if (NXgetchardata (file_id, NXM_run_lab) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Lab : ', NXM_run_lab

C read beamline
if (NXopendata (file_id, 'beamline') .NE. NX_OK) stop
if (NXgetchardata (file_id, NXM_run_beamline) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Beamline : ', NXM_run_beamline
.
.
C @@@@@@@@@@@@@@@@@@ close NXentry group @@@@@@@@@@@@@@@@@@
if (NXclosegroup (file_id) .NE. NX_OK) stop
if (debug) print *, ' '
if (debug) print *, 'Closing NXentry...'

C @@@@@@@@@@@@@@@@@@ close NeXus file @@@@@@@@@@@@@@@@@@
if (NXclose (file_id) .NE. NX_OK) stop
if (debug) print *, ' '
if (debug) print *, 'Closing NeXus file...'

NXMread = 1
END

```

After declaring the variables to store the NeXus data, we must create a set of common blocks identical to those in the MUON_DEF. This will share the memory between the reader and the user' application.

The Reader then navigates the each group in the NeXus file, and places the appropriate data in the variables defined above.

5.3.3 TEST_NEXUS_READER.F77

```
program test_nexus_reader
include 'muon_def_F77.inc'
character*60 fname
integer status, i, j, k, histogram
```

The users application must include the appropriate MUON_DEF.

```
WRITE (*,*) 'Please Enter Nexus File Name'
READ (*,*) fname
print *, ' '
print *, ' '
status = NXMread (fname)
```

Then they must call the NXMread routine supplying a valid NeXus filename.

```
print *, 'NeXus file Contents...'
print *, ' '
```

```
C  NXrun -----
print *, 'Program name : ',NXM_run_program_name
print *, 'Program version : ',NXM_run_program_version
print *, 'Run number : ', NXM_run_number
print *, 'Title : ', NXM_run_title
print *, 'Notes : ', NXM_run_notes
print *, 'Start time : ', NXM_run_start_time
print *, 'Stop time : ', NXM_run_stop_time
print *, 'Duration : ', NXM_run_duration
print *, 'Switching states : ', NXM_run_switching_states
```

The data is then read into memory and is easily accessible, as shown in the example code.

```
C  NXuser-----
print *, 'User name : ', NXM_user_name
print *, 'Experiment number : ', NXM_user_experiment_number
```

```
C  NXdata-----
print *, 'Number of histograms : ', NXM_histogram_number
print *, 'Time histogram length : ', NXM_histogram_length
print *, 'Time zero bin : ', NXM_histogram_t0_bin
print *, 'First good bin : ', NXM_histogram_first_good_bin
print *, 'Last good bin : ', NXM_histogram_last_good_bin
print *, 'Histogram resolution : ', NXM_histogram_resolution, ' in
',NXM_histogram_resolution_units
print *, 'Histogram units : ', NXM_histogram_counts_units
```

```
.
```

```
C  NXlog-----
if (NXM_temp_log_available.gt.1) then
  print *, ' '
  print *, 'Temperature Log - time (seconds), temperature (kelvin)'
  do i=1, NXM_temp_log_available
    print *, '(', NXM_temp_log_time(i), ',', NXM_temp_log_values(i), ')'
  end do
end if
```

```
C  NXlog-----
if (NXM_event_log_available.gt.0) then
  print *, ' '
  print *, 'Event Log - time (seconds), events (counts)'
  do i=1, NXM_event_log_available
    print *, '(', NXM_event_log_time(i), ',', NXM_event_log_values(i), ')'
  end do
end if
end
```

5.3.4 Compilation Instructions


To recompile this application, the following files will be required along with a standard DEC Fortran compiler:

- NAPIF.INC
- NEXUS_READER.FOR
- MUON_DEF_F77.INC or MUON_DEF_F90.INC (for either Fortran 77 or Fortran 90 respectively)
- TEST_NEXUS_READER.FOR or TEST_NEXUS_READER.F90 (or suitably modified for a user application)

In these instructions it is assumed that the HDF libraries and NeXus API's are currently installed. If these are not yet installed, please refer to the NeXus Data Format Homepage (http://www.neutron.anl.gov/nexus/NeXus_API.html#install) for detailed instructions.

Step 1 - Compile the above Fortran files in sequence

Step 2 - Link TEST_NEXUS_READER, NEXUS_READER, NeXus API (where NeXus API exists)

A screenshot of a DECterm window titled "DECterm 1". The window has a menu bar with "File", "Edit", "Commands", "Options", "Print", and "Help". The main text area contains the following Fortran code:

```
$set verify
$define hdf_root axplib$disk:[4_1r3_alphavms7_1.] ! where you installed HDF, trailing "." is needed
$define nexus_root user$disk:[df01.] ! where the nexus library exists
$FOR/ extend nexus_reader.for
$FOR/ extend test_nexus_reader.for
$link test_nexus_reader, nexus_reader, axplib$disk:[nexus]nexus_f77/opt
```

Figure 5.2 – Fortran NeXus_Reader Compilation Instructions

5.3.5 Limitations

The Fortran NeXus reader has a number of limitations. Of particular importance is the inability of the language to support dynamic memory allocation, which means that we have had to statically declare the sizes of the various arrays. The Reader is also restricted to reading only one NXdata group. At present, however, this is acceptable for the ISIS muon data files because it is very rare that an experiment has differing histogram properties. Also, the reader only reads one set of temperature log data and one set of event log data.

5.4 Error Messages

The error messages produced from each reader have been made as consistent as possible.

Below is a table listing each error message and describing why each message has appeared.

	<i>Message</i>	<i>Meaning</i>
Error	"Fatal Error : Read routine incompatible with instrument definition used to write data"	The Instrument definition is expected to change over time. The situation may occur where an old version of a NeXus file is read by an up to date NeXus read routine.
Warning	"Magnetic field vector not available"	This message appears if the reader cannot find any magnetic field vector values
Warning	"Detector angles not available"	This message appears if the reader cannot find any detector angle values
Warning	"Detector dead-times not available"	This message appears if the reader cannot find any dead-times
Warning	"Histogram time zero not available"	This message appears if the reader cannot find a time zero value
Warning	"No groups available"	This message appears if the reader cannot find any grouping information
Warning	"No alpha pairs available"	This message appears if no alpha pairs are present
Diagnostic	<element name> " : " <element value>	As the reader progresses through the nexus file it displays the information that it reads e.g. "lab : ISIS"
Diagnostic	<element name> " units : " <element units>	Some elements in the nexus file have units associated with them. E.g. "Histogram_resolution units : picoseconds"
Diagnostic	<element name> " coordinate system : " <element coordinate system>	Some elements in the NeXus file have a coordinate system associated. E.g. "Detector Angles coordinate system : spherical"
Diagnostic	"reading " <element name>	When reading some of the more complex elements, the read routine displays a message to let the user know that it is reading a value or set of values
Diagnostic	"opening " <group name>	As the read routine traverses the NeXus file it notifies the user when it opens a group.
Diagnostic	"closing " <group name>	Likewise the reader notifies the user when the reader has closed a group

Figure 5.3 – NeXus_Reader Error Messages Table

6. TMOGGER

6.1 Tmogger User Manual

Tmogger is a hybrid of the existing applications “Tlogger” and “MACQlogger”. Tmogger enables the user to plot temperature and event log information individually or together on the same page. Tmogger is written in Fortran 90 and uses the Fortran 77 NeXus read routine to interpret the user supplied NeXus files.

Advantages of Tmogger:

- Can display Temperature Log data.
- Can display Event log data.
- Temperature log data has been collated and time stamped from all the different versions of the temperature log file.
- Event log data has been time-stamped.
- Can display all the data available in the NeXus file.

6.1.1 How To Compile Tmogger


To recompile this application, the source code for the following files will be required along with a standard Fortran 90 compiler:

- Nexus_reader.for
- Tmogger_module.f90
- Tmogger.f90
- Napif.inc
- Muon_def.inc

In these instructions it is assumed that the HDF libraries and NeXus API's are currently installed. If these are not yet installed please refer to the NeXus Data Format Homepage (http://www.neutron.anl.gov/nexus/NeXus_API.html#Install) for detailed instructions.

Step 1 – Compile above Fortran files in sequence

Step 2 – Link Tmogger, Tmogger_module, NeXus F77 API (where Nexus Library exists)

A screenshot of a DECterm window titled "DECterm 1". The window has a menu bar with "File", "Edit", "Commands", "Options", "Print", and "Help". The main text area contains the following commands:

```
$set verify
$define hdf_root axplib$disk:[4_lr3_alphavms7_1.] ! where you installed HDF, trailing "." is needed
$define nexus_root user$disk:[df01.] ! where the nexus library exists
$For/extend nexus_reader.for
$f90/noop/debug tmogger_module.f90
$f90/noop/debug tmogger.f90
$link tmogger,tmogger_module,nexus_reader, axplib$disk:[nexus]nexus_f77/opt, pgplot/opt
```

Figure 6.1 – Tmogger Compilation Instructions Screenshot

6.1.2 How To Set-Up Tmogger

This manual guides the user through a step by step process of utilising the features of Tmogger.

To run the application the user must first initialise the application. This is achieved by typing: “@ muon\$disk:[muonmgr.software.utilities.NeXus]NeXus_setup” at the ISIS VMS command prompt (See Figure 3.2)

The above command will execute a command file, which will install a suite of muon NeXus applications in the users current directory. The applications installed are:


- *Convert_NeXus* – Application which converts ISIS Muon MCS data files into NeXus format.
- *NXbrowse* – Open source NeXus file browser.
- *NeXus_UDA* - μ SR Data Analysis program which can interpret ISIS Muon NeXus files.
- *F77_NeXus_Reader* – ISIS Muon NeXus read routine written in F77 (allows users to read NeXus data into their own programs).
- *C_NeXus_Reader* - ISIS Muon NeXus read routine written in C (allows users to read NeXus data into their own programs).
- *Tmogger* – Application which plots temperature and event logs stored in ISIS Muon NeXus files.
- *UIF_creator* – Application which produces UIF files for Convert_NeXus.

It should also be noted that the command file copies across any additional information required by these applications including a selection of test data, and will also initialise any logical variables needed.

We will only concern ourselves with the ‘Tmogger’ application here. For more information on the other available software please refer to the Muon NeXus homepage (www.isis.rl.ac.uk/muons/nexus).

6.1.3 Getting Started

Now that the suite of NeXus applications have been set-up, Tmogger can be run by typing “*TMOGGER*” at the command line.



```
DECterm 1
File Edit Commands Options Print Help
$ TMOGGER
Welcome, you are running TMOGGER version 1
This program plots temperature and/or event log information
from the supplied NeXus file
If you have any problems or queries, please contact
Damian.Flannery@rl.ac.uk


1. T logger          - Temperature Values
2. MACQ logger       - Event Values
3. T logger and MACQ logger - Temperature and Event Values

Please choose option :
```

1

When the user executes Tmogger a welcome message appears on the screen along with a short menu. This menu asks the user which values they

would like to plot on the screen. Option 1 plots temperature log values on the screen. Option 2 plots event log values on the screen. Option 3 plots both temperature and event logs together on the screen.



```
DECterm 1
File Edit Commands Options Print Help
This program plots temperature and/or event log information
from the supplied NeXus file
If you have any problems or queries, please contact
Damian.Flannery@rl.ac.uk


1. T logger          - Temperature Values
2. MACQ logger       - Event Values
3. T logger and MACQ logger - Temperature and Event Values

Please choose option :
3

Please enter NeXus filename :
35000.nxs
```

2

Once the user has chosen which values to plot, Tmogger requests a NeXus filename. (Note. This is relative from the current directory)



```

DECTerm 1
File Edit Commands Options Print Help
Please enter NeXus filename :
35000.nxs
Opening NeXus file...
Opening "NXentry"...

Nexus reader version      1
Compatible with ISIS Muon Instrument Definition version      1

IDf version :              1
Analysis : muonTD
Lab : ISIS
Beamline : EMU
Program Name : MCS
Program Version : RD
Run Number : 35000
Title : Temperature = 10 K, Magnetic Field = 0 Gauss, Sample = LACUO
Notes : slits=10mm, 2f
Start Time : 2000-07-14 11:14:50
Stop Time : 2000-07-14 11:38:29
Duration : 00:23:39
Switching States : 1

Opening NXuser...
Data name: name
Data value : PAVIA
Data name: experiment_number
Data value : 11293
UIF array length : 0
Closing NXuser..

Opening NXsample...
Sample name : LACUOZ30
Sample temperature : 10.00000
Temperature units :
kelvin
Sample magnetic field : 0.0000000E+00
Field units : gauss
Sample shape : n/a
Field state : 0.0000000E+00
Field units : [ 0.0000000E+00 , 0.0000000E+00 , 1.000000 ]
Field Vector coordinate system :
cartesian
Sample environment :
n/a
Closing NXsample


Opening NXinstrument...

```

3

Assuming that the user has supplied a valid NeXus filename and the NeXus file conforms to the ISIS Muon Instrument Definition version 1, the contents of the NeXus file are dumped to the screen.

This information is read by the Fortran 77 version of the NeXus reader and prints the contents of the file as it traverses each element.



```

DECTerm 1
File Edit Commands Options Print Help
Histogram number : 32
Histogram length : 1000
Histogram to bin : 19
Histogram first good bin : 26
Histogram last good bin : 1000
Histogram units : counts
Histogram resolution : 16000
Histogram resolution units :
picoseconds
Reading histogram time zero...
Histogram time zero : 0.2780000 value found
Histogram time zero units :
picoseconds
Reading raw time...
Histogram raw time units :
microseconds
Reading corrected time...
Histogram corrected time units :
microseconds
Reading histogram groupings...
2 groups found
Reading alpha for grouped pairs...
1 alpha pairs found
Closing NXdata...

Opening NXlog..
Logged name : sample temperature
260 values available
Reading temperature values...
First temperature value logged : 10.02500
Reading time values...
Logging started at time : 4.000000
Closing NXlog..

Opening NXlog..
Logged name : ISIS beam
14 values available
Reading event values...
First event value logged : 0.5627000
Reading time values...
Logging started at time : 103.0000
Closing NXlog..

Closing NXentry..
Closing NeXus file...
Graphics device/type (? to see list, default /NULL): /xwindow

```


4

Once the diagnostics have been printed to the screen, a prompt appears asking the user which graphics device they would like to choose. Type "/xwindow" here. This will plot the desired data in a separate window. (For a list of other devices available, type "?". For further information view pgplot documentation)

5


After the user has selected “/xwindow” as the display device, a pgplot window appears with the plotted data.


This is the temperature and event log data for run 35000.nxs



6


If the user had chosen to plot just temperature log data in step 2, this is what the user would see.





7

If the user had chosen to just plot event log data in step 2, this is what the user would see.



8

To discard the plotting window and continue, the user must select the VMS prompt window and press <RETURN>

(as above). Tmogger will then ask the user if they require a laser hardcopy. Please select "y" or "n". If the user selects no, Tmogger exits. If the user selects yes, Tmogger produces a postscript file called pgplot.ps. The user must then print this file to the printer of their choice.

6.2 Tmogger Design

The Tmogger application has been written in Fortran 90 and uses the F77 NeXus read routine. Tmogger has one supporting module, “Tmogger_module” which is designed to handle the plotting of treated data. This is described in more detail in the following section. The application firstly asks the user which data they wish to plot and secondly the filename of the NeXus file where they wish to extract this information. Tmogger then performs the appropriate calculations for the selected data – for example RMS deviation, Average Events/Temperature, minimum and maximum values etc. It then supplies this data to the appropriate plotting subroutine.




Figure 6.3 – A diagram of external module association

6.2.1 Tmogger Module

The Tmogger consists of three main methods, which all involve plotting data. Within these routines are calls to PGPLOT which is a Fortran written graphics plotting package that Tmogger links to.

Plot_macq_data (..) → *Plots Event Log Data*

Plot_dat (..) → *Plots Temperature Log Data*

Plot_tlog_and_macq_data (..) → *Plots Temperature and event Log Data*

and a function method which obtains which information the user wishes to plot.

Get_application_choice () → *Asks the user which log they wish to plot*

These methods are described in more detail below:

6.2.2 Get_application_choice ()

Input arguments : *none*

Output arguments : application_choice int

Displays a menu on the screen asking the user which log information they wish to plot. The function then reads the input from the user and returns this value.

6.2.3 Plot_macq_data ()

Input arguments:

MIN	Real	MAX	Real
YMIN	Real	YMAX	Real
READS	Int	RUN	Int
TIME	Real []	EVENTS	Real []
AVERAGE	Real	RMS	Real
TITLE	String	STIME	String
FTIME	String		

Output Arguments: *none*

This method receives a number of parameters, listed above, which contain the data to be plotted, along with any other labels, statements or calculations. The temperature readings are then plotted as a function of time.

6.2.4 Plot_dat ()

Input arguments:

XMIN	real	XMAX	Real
YMIN	Real	YMAX	Real
READS	Int	RUN	Int
TIME	Real []	TEMP	Real []
AVERAGE	Real	RMS	Real
TITLE	String	STIME	String
FTIME	String	DUR	String

Output arguments: none

This method also receives a number of parameters, as listed above, which contain the data to be plotted, along with any other labels, statements or calculations. The event readings are then plotted as a function of time.

6.2.5 Plot_tlog_and_macq_data ()

Input arguments:

T_XMIN	Real	M_XMIN	Real
T_XMAX	Real	M_XMAX	Real
T_YMIN	Real	M_YMIN	Real
T_YMAX	Real	M_YMAX	Real
T_READS	Int	M_READS	Int
RUN	Int	T_TIME	Real [] []
M_TIME	Real [] []	EVENTS	Real []
TEMP	Real []	T_AVERAGE	Real
M_AVERAGE	Real	T_RMS	Real
M_RMS	Real	TITLE	String
SET_TEMP	Real	SAMPLE_NAME	String
SET_FIELD	Real	S_TIME	String
F_TIME	String		

Output arguments: *none*

This method also receives the parameters listed above, which contain the data to be plotted, along with any other labels, statements or calculations for both temperature and event logs. The event readings are then plotted as a function of time, beneath the temperature log values which are also plotted as a function of time.

7. NEXUS_UDA

7.1 Introduction

UDA is an easy-to-use μ SR data analysis program written in Fortran which allows scientists to plot histograms and perform analysis of their data. A new version of UDA has been developed so that it can read NeXus files as well as the current binary data files produced by the present data acquisition system (MCS). What follows is a description of how to run this new version of UDA, and also a short tutorial on how to read a NeXus file into UDA. For a more detailed article on UDA see <http://www.isis.rl.ac.uk/muons/emu/eguide/chapter6/index.htm>.

7.2 How To Set Up UDA

To run the application the user must first initialise the application. This is can be achieved by typing: “@ muon\$disk:[muonmgr.software.utilities.NeXus]NeXus_setup” at the ISIS VMS command prompt (See Figure 3.2).

The above command will execute a command file, which will install a suite of muon NeXus applications in the users current directory. The applications installed are:

- *Convert_NeXus* – Application which converts ISIS Muon MCS data files into NeXus format.
- *NXbrowse* – Open source NeXus file browser.
- *NeXus_UDA* - μ SR Data Analysis program which can interpret ISIS Muon NeXus files.
- *F77_NeXus_Reader* – ISIS Muon NeXus read routine written in F77 (allows users to read NeXus data into their own programs).
- *C_NeXus_Reader* - ISIS Muon NeXus read routine written in C (allows users to read NeXus data into their own programs).
- *Tmogger* – Application which plots temperature and event logs stored in ISIS Muon NeXus files.
- *UIF_creator* – Application which produces UIF files for Convert_NeXus.

It should also be noted that the command file copies across any additional information needed by these applications including a selection of test data, and will also initialise any logical variables needed.


We will only concern ourselves with the 'NeXus_UDA' application here. For more information on other available software please refer to the Muon NeXus homepage (www.isis.rl.ac.uk/muons/nexus).

7.3 How to use UDA

Below is a short step by step tutorial describing how to read a NeXus file into UDA. This is not a tutorial detailing the capabilities of UDA. For a more detailed article on UDA see <http://www.isis.rl.ac.uk/muons/emu/eguide/chapter6/index.htm>.

To execute NeXus_UDA just type "NEXUS_UDA"

Step 1 – Choose NeXus option on the main menu.




Step (1)

After Successfully running UDA v14.0. You will be presented with the opening menu (left).

Using the cursor keys or by pressing 'N', select the Nexus option and press enter

Step 2 – Enter the run number for the file you wish to display. Note. This version of UDA points to EMU MCS binary data files, and will look in the current directory for NeXus files.




Step (2)

At the bottom of the screen, a prompt will appear asking the user to enter a run number.

Here type a run number of a nexus file e.g. 30000. And then press Enter.

Step 3 – View the header information. Now the data is in memory, you can use UDA's capabilities to analyse the information.



Step (3)

At this point, UDA reads the nexus file and prints the header information to the screen.

Press Enter to continue back to the main menu where you can perform the standard analytical and plotting functions.

Note: If the user types the run number of a file that doesn't exist; an error message will appear and control will return back to the UDA main menu.

If you have any trouble plotting data using UDA, you may need to set your terminal display. You can do this by typing "set display/create/transport=tcpip/node=<your ip address>".

8. CONCLUSIONS AND PLANS FOR THE FUTURE

Although adequate for the present data storage requirements, the ISIS muon Instrument Definition will certainly evolve with time as new instrumentation is developed. The key advantage of adopting a structured and self describing format, such as NeXus, is that information can be added without affecting the ability of existing subroutines to read (a subset of) the new data files. Hopefully, other laboratories will work to create their own Instrument Definitions and adopt NeXus in some way, perhaps initially as an intermediate data format with appropriate file converters.

The intention is to continue the development of the NeXus utilities as presented above, in the belief that NeXus offers a viable approach for evolving a common μ SR file format that will enable user to share data and analysis programs between facilities. New developments will be posted on the NeXus pages of our web site (www.isis.rl.ac.uk/muons/nexus).

There are many possibilities as to future work with NeXus, however the killer applications are yet to appear. In the short term a project is underway that will allow NeXus files to be served over the web, possibly enabling simple histograms and logged data sets to be plotted using a Java applet installed in the clients web browser. This mechanism will ensure that users will not have to download, compile or install any software, and potentially offers a truly cross-platform solution for data manipulation.

FURTHER READING

SFJ Cox, Introduction to μ SR : What, How, Where? RAL-94-033, 1994

SFJ Cox, AM Stoneham, Muon Beams, Used for Studying the Solid State, RAL-92-002, 1992

GH Eaton, MA Clarke-Gayther, CA Scott, SFJ Cox, SH Kilcoyne, SP Cottrell, SR Brown and WG Williams, The Muon Beamline at ISIS, RAL-94-077, 1994

T.J. Pearson, PGPLOT Graphics Subroutine Library User Manual, 1989

M Metcalf, J Reid, Fortran 90/95 explained 2nd ed, 1999

ISIS Muon Homepage, <http://www.isis.rl.ac.uk/muons>

NeXus Data Format Homepage, <http://www.neutron.anl.gov/nexus/>

BD Hahn, Fortran 90 for Scientists and Engineers, 1997

NCSA Hierarchical Data Format Homepage, <http://hdf.ncsa.uiuc.edu/>

VMS Library information, <http://www.openvms.digital.com:8000/index.html>

ACKNOWLEDGEMENTS

The author would like to thank Dr. SP Cottrell, Dr. PJC King, Prof. SFJ Cox, Dr. GH Eaton, Dr. JS Lord, Dr. WG Williams, Dr. F Pratt and Dr AD Hillier (ISIS Muon Facility, RAL) for their services in educating a Software Engineer in the field of μ SR!

Dr FA Akeroyd (ISIS Facility, RAL) for his help with Fortran, NeXus and VMS.

Also thanks to T.M Riseman (Birmingham University) for all her useful feedback during the progression of the project.

Special thanks to Dr. SP Cottrell who has dedicated much of his time in working closely on this project.

APPENDIX A: Nexus_reader.h

```
/*
 * Structure definition for NeXus muon file information
 */

struct NXM_MUONUIF
{
    NXname          uif_element_name;
    NXname          uif_element_value;
};

struct NXM_MUONIDF
{
    /* NXrun */
    NXname          run_program_name;
    NXname          run_program_version;
    int             run_idf_version;
    int             run_number;
    NXname          run_title;
    NXname          run_notes;
    NXname          run_analysis;
    NXname          run_lab;
    NXname          run_beamline;
    NXname          run_start_time;
    NXname          run_stop_time;
    NXname          run_duration;
    int             run_switching_states;
    /* NXuser */
    NXname          user_name;
    NXname          user_experiment_number;
    int             user_uif_array_length;
    struct NXM_MUONUIF * user_uif_array;
    /* NXsample */
    NXname          sample_name;
    float           sample_temperature;
    NXname          sample_temperature_units;
    float           sample_magnetic_field;
    NXname          sample_magnetic_field_units;
    NXname          sample_shape;
    NXname          sample_magnetic_field_state;
    int             sample_magnetic_field_vector_available;
    float *         sample_magnetic_field_vector;
    NXname          sample_magnetic_field_vector_coord;
    NXname          sample_environment;
    /* NXinstrument */
    NXname          instrument_name;
    /* NXdetector */
    int             detector_number;
    NXname          detector_orientation;
    int             detector_angles_available;
    float *         detector_angles;
    NXname          detector_angles_coord;
    int             detector_deadtimes_available;
    float *         detector_deadtimes;
    NXname          detector_deadtimes_units;
    /* NXcollimator */
    NXname          collimator_type;
    NXname          collimator_aperture;
    /* NXbeam */
    float           beam_total_counts;
    NXname          beam_total_counts_units;
    int             beam_daereads;
    int             beam_frames;
    /* NXdata histogram_data_1 */
    NXname          histogram_units;
    int *           histogram_counts;
    int             histogram_number;
    int             histogram_length;
    int             histogram_t0_bin;
    int             histogram_first_good_bin;
    int             histogram_last_good_bin;
};
```

```

    int          histogram_resolution;
    NXname        histogram_resolution_units;
    float         histogram_time_zero;
    int           histogram_time_zero_available;
    NXname        histogram_time_zero_units;
    float *       histogram_raw_time;
    NXname        histogram_raw_time_units;
    float *       histogram_corrected_time;
    NXname        histogram_corrected_time_units;
    int           histogram_grouping_available;
    int *         histogram_grouping;
    int           histogram_alpha_available;
    float *       histogram_alpha;
    /* NXlog temperature_log_1 */
    int           temperature_log_available;
    NXname        temperature_log_name;
    float *       temperature_log_values;
    NXname        temperature_log_values_units;
    float *       temperature_log_time;
    NXname        temperature_log_time_units;
    /* NXlog events_log_1 */
    int           events_log_available;
    NXname        events_log_name;
    float *       events_log_values;
    NXname        events_log_values_units;
    float *       events_log_time;
    NXname        events_log_time_units;
};

/*
 * Function definitions for reading NeXus muon file information
 */

int NXMread(char *nxfname, struct NXM_MUONIDF **nxfddata);
void NXMfreememory(struct NXM_MUONIDF *nxfddata);

```

APPENDIX B: Nexus_reader.c

```
/*
 * Read routine for NeXus files based on an Instrument Definition Function (IDF)
 * for ISIS Muon Time Differential instrument
 *
 *
 * Version 1
 *
 *
 * To use the read routine call:
 *
 *   int NXMread(char *nxfname, struct NXM_MUONIDF **nxpdata)
 *
 * where:
 *
 *   'nxfname' is a pointer to the NeXus filename
 *   'nxpdata' returns a pointer to a pointer to a structure containing the file information
 *
 * the value NX_OK is returned on successful read else NX_ERROR is returned
 *
 * The structure NXM_MUONIDF and function are defined in the file 'nexus_reader.h'
 *
 *
 * After use, the workspace claimed by the routine should be freed using the call:
 *
 *   void NXMfreememory(struct NXM_MUONIDF *nxpdata)
 *
 * where:
 *
 *   'nxpdata' is a pointer to the structure containing the file information
 *
 *
 * Known limitations:
 * Assumes the NXrun group contains a single NXdata group ('histogram_data_1') and
 * two NXlog groups ('temperature_log_1' and 'events_log_1')
 *
 *
 * Tested: 1) Visual C++ 5.0, Windows NT 4.0
 *         2) OpenVMS 7.2
 *
 *
 * D.W.Flannery and S.P.Cottrell
 */

#include <stdio.h>
#include "napi.h"
#include "nexus_reader.h"

#define _DEBUG_ 1 /* switch for debugging information */

/* information locating read routine */
#define _READVERSION_ 1 /* Version number of read routine */
#define _IDFVERSION_ 1 /* Compatible version of Instrument Definition File */
#define _ANALYSIS_ "muonTD" /* Compatible analysis */
#define _LAB_ "ISIS" /* Compatible lab */

/* structure defining file information */
struct NXM_MUONIDF muon;
```

```

/* macro for adding debug information */
#define NXMdebug(format, value) \
    if (_DEBUG_) printf(format, value)

/* macro handling function calling, returns with NX_ERROR if function call failed */
#define NXMerrorhandler(function) \
    if (function != NX_OK) return NX_ERROR

static int NXMdatasize(int dataType)

/*
 * function returning size (in bytes) of NeXus data type
 */
{
    int size;

    switch (dataType)
    {
        case NX_CHAR:
            size = SIZE_CHAR8;
            break;
        case NX_FLOAT32:
            size = SIZE_FLOAT32;
            break;
        case NX_FLOAT64:
            size = SIZE_FLOAT64;
            break;
        case NX_INT8:
            size = SIZE_INT8;
            break;
        case NX_UINT8:
            size = SIZE_UINT8;
            break;
        case NX_INT16:
            size = SIZE_INT16;
            break;
        case NX_UINT16:
            size = SIZE_UINT16;
            break;
        case NX_INT32:
            size = SIZE_INT32;
            break;
        case NX_UINT32:
            size = SIZE_UINT32;
            break;
        default:
            size = 0;
            break;
    }

    return size;
}

static int NXMgetstringdata(NXhandle file_id, char *str)

/*
 * read a rank 1 data item of type NX_CHAR from a NeXus file
 * add a null terminator character for compatibility with C string handling
 *
 * assumes space reserved for string of length 'VGNAMELENMAX', as defined by 'NXname'
 */
{
    int i, rank, type, dims[32];
    NXname data_value;

    NXMerrorhandler((NXgetinfo(file_id, &rank, dims, &type)));

```

```

    if ((type != NX_CHAR) || (rank > 1) || (dims[0] >= VGNAMELENMAX))
    {
        printf("\nFatal error in routine 'getstringdata'\n\n");
        exit(0);
    }

    NXErrorHandler((NXgetdata(file_id, data_value)));

    for (i = 0; i < dims[0]; i++)
        *(str + i) = *(data_value + i);
    *(str + i) = '\0';

    return 1;
}

static int NXMgetstringattr(NXhandle file_id, char *name, char *str)
/*
 * read an attribute of type NX_CHAR from a NeXus file
 * add a Null terminator character for compatibility with C string handling
 *
 * assumes space reserved for string of length 'VGNAMELENMAX', as defined by 'NXname'
 */
{
    int i, attlen, atttype;
    NXname data_value;

    attlen = VGNAMELENMAX - 1;
    atttype = NX_CHAR;
    NXErrorHandler((NXgetattr(file_id, name, data_value, &attlen, &atttype)));

    for (i = 0; i < attlen; i++)
        *(str + i) = *(data_value + i);
    *(str + i) = '\0';

    return 1;
}

static void *NXMmemoryhandler(NXhandle file_id)
/*
 * allocate memory for a NeXus data item, calls NXgetinfo to determine size
 *
 * function returns a pointer to memory block
 */
{
    int i, rank, type, dims[32], size;
    char *data_ptr;

    NXErrorHandler((NXgetinfo(file_id, &rank, dims, &type)));
    size = dims[0];
    for (i = 1; i < rank; i++)
        size = size * dims[i];
    size = size * NXMdatasize(type);
    if ((data_ptr = malloc (size)) == NULL)
    {
        printf("\nMemory allocation error: can't allocate %d bytes\n\n", size);
        exit(0);
    }
}

static void NXMfree(void *data_ptr)
/*
 *
 * free memory from NeXus data items provided pointer argument is non NULL
 */

```

```

{
    if (data_ptr != NULL) free(data_ptr);
}

int NXMread(char *nxfname, struct NXM_MUONIDF **nxfdata)

/*
 * return data from NeXus file 'nxfname' in structure 'nxfdata'
 *
 * function returns NX_OK on success, otherwise NX_ERROR
 */

{
    NXhandle file_id;
    NXname group_name, class_name, data_name;
    int i, attlen, atttype, numitems, data_type, uif_count;

    *nxfdata = &muon;

    /* open NeXus file for reading */
    NXMerrorhandler((NXopen(nxfname, NXACC_READ, &file_id)));
    NXMdebug("%s\n", "Opening NeXus file...");

    /* open root group 'NXentry' */
    NXMerrorhandler((NXopengroup(file_id, "run", "NXentry")));
    NXMdebug("%s\n", "Opening 'NXentry'...");

    /* display (in debug mode) NeXus reader version and IDF compatible version
    number */
    NXMdebug("\nNeXus reader version %d\n", _READVERSION_);
    NXMdebug("Compatible with Instrument Definition File version %d\n\n",
        _IDFVERSION_);

    /* read IDF version */
    NXMerrorhandler((NXopendata(file_id, "idf_version")));
    NXMerrorhandler((NXgetdata(file_id, &muon.run_idf_version)));
    NXMerrorhandler((NXclosedata(file_id)));
    NXMdebug("IDF version: %d\n", muon.run_idf_version);

    /* read 'analysis' */
    NXMerrorhandler((NXopendata(file_id, "analysis")));
    NXMerrorhandler((NXMgetstringdata(file_id, muon.run_analysis)));
    NXMerrorhandler((NXclosedata(file_id)));
    NXMdebug("Analysis: %s\n", muon.run_analysis);

    /* read 'lab' */
    NXMerrorhandler((NXopendata(file_id, "lab")));
    NXMerrorhandler((NXMgetstringdata(file_id, muon.run_lab)));
    NXMerrorhandler((NXclosedata(file_id)));
    NXMdebug("Lab: %s\n", muon.run_lab);

    /* check reader compatibility */
    if ((muon.run_idf_version != _IDFVERSION_) || (strcmp(muon.run_analysis,
        _ANALYSIS_) || (strcmp(muon.run_lab, _LAB_)))
    {
        printf("\nFatal error: Read routine incompatible with Instrument
            Definition used to write data\n");
        printf("Found, IDF Version: %d, Analysis: %s, Lab: %s\n",
            muon.run_idf_version, muon.run_analysis, muon.run_lab);
        printf("Expecting: IDF Version: %d, Analysis: %s, Lab: %s\n\n",
            _IDFVERSION_, _ANALYSIS_, _LAB_);
        exit(0);
    }

    /* read 'beamline' */
    NXMerrorhandler((NXopendata(file_id, "beamline")));
    NXMerrorhandler((NXMgetstringdata(file_id, muon.run_beamline)));
    NXMerrorhandler((NXclosedata(file_id)));
    NXMdebug("Beamline: %s\n", muon.run_beamline);

    /* read 'program_name' and 'program_version' */
    NXMerrorhandler((NXopendata(file_id, "program_name")));
    NXMerrorhandler((NXMgetstringdata(file_id, muon.run_program_name)));
    NXMerrorhandler((NXMgetstringattr(file_id, "version",
        muon.run_program_version)));

```

```

NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Program name: %s\n", muon.run_program_name);
NXMdebug("Program version: %s\n", muon.run_program_version);

/* read 'number' */
NXMErrorHandler((NXopendata(file_id, "number")));
NXMErrorHandler((NXgetdata(file_id, &muon.run_number)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Run number: %d\n", muon.run_number);

/* read 'title' */
NXMErrorHandler((NXopendata(file_id, "title")));
NXMErrorHandler((NXMgetstringdata(file_id, muon.run_title)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Title: %s\n", muon.run_title);

/* read 'notes' */
NXMErrorHandler((NXopendata(file_id, "notes")));
NXMErrorHandler((NXMgetstringdata(file_id, muon.run_notes)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Notes: %s\n", muon.run_notes);

/* read 'start_time' */
NXMErrorHandler((NXopendata(file_id, "start_time")));
NXMErrorHandler((NXMgetstringdata(file_id, muon.run_start_time)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Start time: %s\n", muon.run_start_time);

/* read 'stop_time' */
NXMErrorHandler((NXopendata(file_id, "stop_time")));
NXMErrorHandler((NXMgetstringdata(file_id, muon.run_stop_time)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Stop time: %s\n", muon.run_stop_time);

/* read 'duration' */
NXMErrorHandler((NXopendata(file_id, "duration")));
NXMErrorHandler((NXMgetstringdata(file_id, muon.run_duration)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Duration: %s\n", muon.run_duration);

/* read 'switching_states' (rgmode) */
NXMErrorHandler((NXopendata(file_id, "switching_states")));
NXMErrorHandler((NXgetdata(file_id, &muon.run_switching_states)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Switching States: %d\n", muon.run_switching_states);

/* open subgroup 'NXuser' */
NXMErrorHandler((NXopengroup(file_id, "user", "NXuser")));
NXMdebug("%s\n", "Opening 'NXuser'...");

/* read user information - two fixed entries and the information from the UIF
file */
NXMErrorHandler((NXinitgroupdir(file_id)));
NXMErrorHandler((NXgetgroupinfo(file_id, &numitems, group_name, class_name)));
NXMdebug("Group Info: %s", class_name);
NXMdebug(" ", %s, group_name);
NXMdebug(" ", %d\n", numitems);
if ((numitems - 2) > 0)
{
    if ((muon.user_uif_array = (struct NXM_MUONUIF *) malloc((numitems - 2) *
        sizeof(struct NXM_MUONUIF))) == NULL)
    {
        printf("\nUnable to allocate memory for User Information
            array\n\n");
        exit(0);
    }
}
else
{
    muon.user_uif_array = (struct NXM_MUONUIF *) NULL;
}
uif_count = 0;
for (i = 0; i < numitems; i++)
{
    NXMErrorHandler((NXgetnextentry(file_id, data_name, class_name,
        &data_type)));
    NXMErrorHandler((NXopendata(file_id, data_name)));

```

```

NXMdebug("Data name: %s\n", data_name);
if (strcmp(data_name, "name") == 0)
{
    NXMerrorhandler((NXMgetstringdata(file_id, muon.user_name)));
    NXMdebug("User name: %s\n", muon.user_name);
}
else if (strncmp(data_name, "experiment_number", 2) == 0)
{
    NXMerrorhandler((NXMgetstringdata(file_id,
    muon.user_experiment_number)));
    NXMdebug("Experiment number: %s\n", muon.user_experiment_number);
}
else
{
    NXMerrorhandler((NXMgetstringdata(file_id, (char *)
    muon.user_uif_array[uif_count].uif_element_value)));
    strcpy((char *) muon.user_uif_array[uif_count].uif_element_name,
    (char *) data_name);
    NXMdebug("Data value %d: ", uif_count);
    NXMdebug("%s\n",
    muon.user_uif_array[uif_count].uif_element_value);
    uif_count++;
}
NXMerrorhandler((NXclosedata(file_id)));
}
muon.user_uif_array_length = uif_count;
NXMdebug("UIF array length: %d\n", muon.user_uif_array_length);

/* close subgroup 'NXuser' */
NXMerrorhandler((NXclosegroup(file_id)));
NXMdebug("%s\n", "Close 'NXuser'");

/* open subgroup 'NXsample' */
NXMerrorhandler((NXopengroup(file_id, "sample", "NXSample")));
NXMdebug("%s\n", "Opening 'NXsample'...");

/* read 'name' */
NXMerrorhandler((NXopendata(file_id, "name")));
NXMerrorhandler((NXMgetstringdata(file_id, muon.sample_name)));
NXMerrorhandler((NXclosedata(file_id)));
NXMdebug("Sample name: %s\n", muon.sample_name);

/* read 'temperature' */
NXMerrorhandler((NXopendata(file_id, "temperature")));
NXMerrorhandler((NXgetdata(file_id, &muon.sample_temperature)));
NXMerrorhandler((NXMgetstringattr(file_id, "units",
muon.sample_temperature_units)));
NXMerrorhandler((NXclosedata(file_id)));
NXMdebug("Temperature label: %f\n", muon.sample_temperature);
NXMdebug("Temperature units: %s\n", muon.sample_temperature_units);

/* read 'magnetic_field' */
NXMerrorhandler((NXopendata(file_id, "magnetic_field")));
NXMerrorhandler((NXgetdata(file_id, &muon.sample_magnetic_field)));
NXMerrorhandler((NXMgetstringattr(file_id, "units",
muon.sample_magnetic_field_units)));
NXMerrorhandler((NXclosedata(file_id)));
NXMdebug("Magnetic field: %f\n", muon.sample_magnetic_field);
NXMdebug("Field units: %s\n", muon.sample_magnetic_field_units);

/* read 'shape' */
NXMerrorhandler((NXopendata(file_id, "shape")));
NXMerrorhandler((NXMgetstringdata(file_id, muon.sample_shape)));
NXMerrorhandler((NXclosedata(file_id)));
NXMdebug("Shape: %s\n", muon.sample_shape);

/* read 'magnetic_field_state' */
NXMerrorhandler((NXopendata(file_id, "magnetic_field_state")));
NXMerrorhandler((NXMgetstringdata(file_id, muon.sample_magnetic_field_state)));
NXMerrorhandler((NXclosedata(file_id)));
NXMdebug("Magnetic field state: %s\n", muon.sample_magnetic_field_state);

/* read 'magnetic_field_vector' and 'coordinate_system' attribute */
NXMerrorhandler((NXopendata(file_id, "magnetic_field_vector")));
attlen = SIZE_INT32;
atttype = NX_INT32;
NXMerrorhandler((NXgetattr(file_id, "available",

```



```

        &muon.sample_magnetic_field_vector_available, &attlen, &atttype));
if (muon.sample_magnetic_field_vector_available != 0)
{
    muon.sample_magnetic_field_vector = (float *) NXMemoryhandler(file_id);
    NXErrorHandler((NXgetdata(file_id, muon.sample_magnetic_field_vector)));
    NXMdebug("Magnetic field vector: [%f, ",
        muon.sample_magnetic_field_vector[0]);
    NXMdebug("%f, ", muon.sample_magnetic_field_vector[1]);
    NXMdebug("%f]\n", muon.sample_magnetic_field_vector[2]);
    NXErrorHandler((NXMgetstringattr(file_id, "coordinate_system",
        muon.sample_magnetic_field_vector_coord)));
    NXMdebug("Coordinate system: %s\n",
        muon.sample_magnetic_field_vector_coord);
}
else
{
    muon.sample_magnetic_field_vector = (float *) NULL;
    NXMdebug("%s\n", "Magnetic field vector not available\n");
}
NXErrorHandler((NXclosedata(file_id)));

/* read 'environment' */
NXErrorHandler((NXopendata(file_id, "environment")));
NXErrorHandler((NXMgetstringdata(file_id, muon.sample_environment)));
NXErrorHandler((NXclosedata(file_id)));
NXMdebug("Sample environment: %s\n", muon.sample_environment);

/* close subgroup 'NXsample' */
NXErrorHandler((NXclosegroup(file_id)));
NXMdebug("%s\n", "Close 'NXsample'");

/* open subgroup 'NXinstrument' */
NXErrorHandler((NXopengroup(file_id, "instrument", "NXinstrument")));
NXMdebug("%s\n", "Open 'NXinstrument'...");

/* read 'name' */
NXErrorHandler((NXopendata(file_id, "name")));
NXErrorHandler((NXMgetstringdata(file_id, muon.instrument_name)));
NXErrorHandler((NXclosedata(file_id)));
NXMdebug("Instrument name: %s\n", muon.instrument_name);

/* open subgroup 'NXdetector' */
NXErrorHandler((NXopengroup(file_id, "detector", "NXdetector")));
NXMdebug("%s\n", "Open 'NXdetector'...");

/* read 'number' */
NXErrorHandler((NXopendata(file_id, "number")));
NXErrorHandler((NXgetdata(file_id, &muon.detector_number)));
NXErrorHandler((NXclosedata(file_id)));
NXMdebug("Number of detectors: %d\n", muon.detector_number);

/* read 'orientation' */
NXErrorHandler((NXopendata(file_id, "orientation")));
NXErrorHandler((NXMgetstringdata(file_id, muon.detector_orientation)));
NXErrorHandler((NXclosedata(file_id)));
NXMdebug("Orientation of detectors: %s\n", muon.detector_orientation);

/* read 'angles' and 'coordinate_system' attribute */
NXErrorHandler((NXopendata(file_id, "angles")));
attlen = SIZE_INT32;
atttype = NX_INT32;
NXErrorHandler((NXgetattr(file_id, "available",
    &muon.detector_angles_available, &attlen, &atttype)));
if (muon.detector_angles_available != 0)
{
    muon.detector_angles = (float *) NXMemoryhandler(file_id);
    NXErrorHandler((NXgetdata(file_id, muon.detector_angles)));
    NXMdebug("%s\n", "Reading detector angles");
    NXErrorHandler((NXMgetstringattr(file_id, "coordinate_system",
        muon.detector_angles_coord)));
    NXMdebug("Coordinate system: %s\n", muon.detector_angles_coord);
}
else
{
    muon.detector_angles = (float *) NULL;
    NXMdebug("%s\n", "not available");
}

```

```

NXMErrorHandler((NXclosedata(file_id)));

/* read 'available' attribute and, if necessary, the 'deadtimes' and 'units' */
NXMErrorHandler((NXopendata(file_id, "deadtimes")));
NXMdebug("%s", "Reading detector deadtimes ... ");
attlen = SIZE_INT32;
atttype = NX_INT32;
NXMErrorHandler((NXgetattr(file_id, "available",
    &muon.detector_deadtimes_available, &attlen, &atttype)));
if (muon.detector_deadtimes_available > 0)
{
    muon.detector_deadtimes = (float *) NXMemoryhandler(file_id);
    NXMErrorHandler((NXgetdata(file_id, muon.detector_deadtimes)));
    NXMErrorHandler((NXMgetstringattr(file_id, "units",
        muon.detector_deadtimes_units)));
    NXMdebug("%s\n", "available");
}
else
{
    muon.detector_deadtimes = (float *) NULL;
    NXMdebug("%s\n", "not available");
}

/* close subgroup 'NXdetector' */
NXMErrorHandler((NXclosegroup(file_id)));
NXMdebug("%s\n", "Close 'NXinstrument'");

/* open subgroup 'NXcollimator' */
NXMErrorHandler((NXopengroup(file_id, "collimator", "NXcollimator")));
NXMdebug("%s\n", "Open 'NXcollimator'...");

/* read 'type' */
NXMErrorHandler((NXopendata(file_id, "type")));
NXMErrorHandler((NXMgetstringdata(file_id, muon.collimator_type)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Orientation of detectors: %s\n", muon.collimator_type);

/* read 'aperture' */
NXMErrorHandler((NXopendata(file_id, "aperture")));
NXMErrorHandler((NXMgetstringdata(file_id, muon.collimator_aperture)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Orientation of detectors: %s\n", muon.collimator_aperture);

/* close subgroup 'NXcollimator' */
NXMErrorHandler((NXclosegroup(file_id)));
NXMdebug("%s\n", "Close 'NXcollimator'");

/* open subgroup 'NXbeam' */
NXMErrorHandler((NXopengroup(file_id, "beam", "NXbeam")));
NXMdebug("%s\n", "Open 'NXbeam'...");

/* read 'total_counts' */
NXMErrorHandler((NXopendata(file_id, "total_counts")));
NXMErrorHandler((NXgetdata(file_id, &muon.beam_total_counts)));
NXMErrorHandler((NXMgetstringattr(file_id, "units",
    muon.beam_total_counts_units)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Total counts: %f\n", muon.beam_total_counts);
NXMdebug("Units: %s\n", muon.beam_total_counts_units);

/* read 'daereads' */
NXMErrorHandler((NXopendata(file_id, "daereads")));
NXMErrorHandler((NXgetdata(file_id, &muon.beam_daereads)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("DAE Reads: %d\n", muon.beam_daereads);

/* read 'frames' */
NXMErrorHandler((NXopendata(file_id, "frames")));
NXMErrorHandler((NXgetdata(file_id, &muon.beam_frames)));
NXMErrorHandler((NXclosedata(file_id)));
NXMdebug("Frames: %d\n", muon.beam_frames);

/* close subgroup 'NXbeam' */
NXMErrorHandler((NXclosegroup(file_id)));
NXMdebug("%s\n", "Close 'NXbeam'");

/* close subgroup 'NXdetector' */

```

```

NXMerrorhandler((NXclosegroup(file_id)));
NXMdebug("%s\n", "Close 'NXdetector'");

/* open subgroup 'NXdata' */
NXMerrorhandler((NXopengroup(file_id, "histogram_data_1", "NXdata")));
NXMdebug("%s\n", "Open 'NXdata'...");

/* read 'counts' */
NXMerrorhandler((NXopendata(file_id, "counts")));
muon.histogram_counts = (int *) NXMmemoryhandler(file_id);
NXMerrorhandler((NXgetdata(file_id, muon.histogram_counts)));
/* read 'counts' attributes */
NXMerrorhandler((NXMgetstringattr(file_id, "units", muon.histogram_units)));
attlen = SIZE_INT32;
atttype = NX_INT32;
NXMerrorhandler((NXgetattr(file_id, "number", &muon.histogram_number, &attlen,
    &atttype)));
NXMerrorhandler((NXgetattr(file_id, "length", &muon.histogram_length, &attlen,
    &atttype)));
NXMerrorhandler((NXgetattr(file_id, "t0_bin", &muon.histogram_t0_bin, &attlen,
    &atttype)));
NXMerrorhandler((NXgetattr(file_id, "first_good_bin",
    &muon.histogram_first_good_bin, &attlen, &atttype)));
NXMerrorhandler((NXgetattr(file_id, "last_good_bin",
    &muon.histogram_last_good_bin, &attlen, &atttype)));
NXMerrorhandler((NXclosedata(file_id)));
NXMdebug("%s\n", "Read histogram data");

/* read 'histogram_resolution' */
NXMerrorhandler((NXopendata(file_id, "resolution")));
NXMerrorhandler((NXgetdata(file_id, &muon.histogram_resolution)));
NXMerrorhandler((NXMgetstringattr(file_id, "units",
    muon.histogram_resolution_units)));
NXMerrorhandler((NXclosedata(file_id)));
NXMdebug("Histogram resolution: %d\n", muon.histogram_resolution);
NXMdebug("Units: %s\n", muon.histogram_resolution_units);

/* read 'histogram_time_zero' */
NXMerrorhandler((NXopendata(file_id, "time_zero")));
attlen = SIZE_INT32;
atttype = NX_INT32;
NXMerrorhandler((NXgetattr(file_id, "available",
    &muon.histogram_time_zero_available, &attlen, &atttype)));
if (muon.histogram_time_zero_available > 0)
{
    NXMerrorhandler((NXgetdata(file_id, &muon.histogram_time_zero)));
    NXMerrorhandler((NXMgetstringattr(file_id, "units",
        muon.histogram_time_zero_units)));
    NXMdebug("Histogram time zero: %f\n", muon.histogram_time_zero);
    NXMdebug("Units: %s\n", muon.histogram_time_zero_units);
}
else
{
    muon.histogram_time_zero = 0;
    NXMdebug("%s\n", "no time zero data found");
}
NXMerrorhandler((NXclosedata(file_id)));

/* read 'raw_time' */
NXMerrorhandler((NXopendata(file_id, "raw_time")));
muon.histogram_raw_time = (float *) NXMmemoryhandler(file_id);
NXMerrorhandler((NXgetdata(file_id, muon.histogram_raw_time)));
/* read 'raw_time' attributes */
NXMerrorhandler((NXMgetstringattr(file_id, "units",
    muon.histogram_raw_time_units)));
NXMerrorhandler((NXclosedata(file_id)));
NXMdebug("%s\n", "Read Raw Time");
NXMdebug("Units: %s\n", muon.histogram_raw_time_units);

/* read 'corrected_time' */
NXMerrorhandler((NXopendata(file_id, "corrected_time")));
muon.histogram_corrected_time = (float *) NXMmemoryhandler(file_id);
NXMerrorhandler((NXgetdata(file_id, muon.histogram_corrected_time)));
/* read 'corrected_time' attributes */
NXMerrorhandler((NXMgetstringattr(file_id, "units",
    muon.histogram_corrected_time_units)));
NXMerrorhandler((NXclosedata(file_id)));

```

```

NXMdebug("%s\n", "Read Corrected Time");
NXMdebug("Units: %s\n", muon.histogram_corrected_time_units);

/* read 'grouping' if attribute 'available' is non-zero */
NXMerrorhandler((NXopendata(file_id, "grouping")));
NXMdebug("%s", "Reading Histogram grouping ... ");
attlen = SIZE_INT32;
atttype = NX_INT32;
NXMerrorhandler((NXgetattr(file_id, "available",
    &muon.histogram_grouping_available, &attlen, &atttype)));
if (muon.histogram_grouping_available > 0)
{
    NXMdebug("%d groups found\n", muon.histogram_grouping_available);
    muon.histogram_grouping = (int *) NXMemoryhandler(file_id);
    NXMerrorhandler((NXgetdata(file_id, muon.histogram_grouping)));
}
else
{
    muon.histogram_grouping = (int *) NULL;
    NXMdebug("%s\n", "no grouping data found");
}

/* read 'alpha' if attribute 'available' is non-zero */
NXMerrorhandler((NXopendata(file_id, "alpha")));
NXMdebug("%s", "Reading alpha for grouped pairs ... ");
attlen = SIZE_INT32;
atttype = NX_INT32;
NXMerrorhandler((NXgetattr(file_id, "available",
    &muon.histogram_alpha_available, &attlen, &atttype)));
if (muon.histogram_alpha_available > 0)
{
    NXMdebug("%d alpha pairs found\n", muon.histogram_alpha_available);
    muon.histogram_alpha = (float *) NXMemoryhandler(file_id);
    NXMerrorhandler((NXgetdata(file_id, muon.histogram_alpha)));
}
else
{
    muon.histogram_alpha = (float *) NULL;
    NXMdebug("%s\n", "no alpha pairs found");
}

/* close subgroup 'NXdata' */
NXMerrorhandler((NXclosegroup(file_id)));
NXMdebug("%s\n", "Close 'NXdata'");

/* open subgroup 'NXlog' (temperature) */
NXMerrorhandler((NXopengroup(file_id, "temperature_log_1", "NXlog")));
NXMdebug("%s\n", "Open 'NXlog' (temperature)...");

/* read 'temperature name' and the 'available' attribute */
NXMerrorhandler((NXopendata(file_id, "name")));
NXMerrorhandler((NXMgetstringdata(file_id, muon.temperature_log_name)));
attlen = SIZE_INT32;
atttype = NX_INT32;
NXMerrorhandler((NXgetattr(file_id, "available",
    &muon.temperature_log_available, &attlen, &atttype)));
NXMerrorhandler((NXclosedata(file_id)));
NXMdebug("Logged name: %s, ", muon.temperature_log_name);
NXMdebug("%d values available\n", muon.temperature_log_available);

/* if 'available' then read ... */
if (muon.temperature_log_available > 0)
{
    /* 'temperature values' and 'length' and 'units' attributes */
    NXMerrorhandler((NXopendata(file_id, "values")));
    muon.temperature_log_values = (float *) NXMemoryhandler(file_id);
    NXMerrorhandler((NXgetdata(file_id, muon.temperature_log_values)));
    NXMerrorhandler((NXMgetstringattr(file_id, "units",
        muon.temperature_log_values_units)));
    NXMerrorhandler((NXclosedata(file_id)));
    NXMdebug("First temperature value logged: %f\n",
        muon.temperature_log_values[0]);

    /* 'temperature time' and 'units' attribute */
    NXMerrorhandler((NXopendata(file_id, "time")));
    muon.temperature_log_time = (float *) NXMemoryhandler(file_id);
    NXMerrorhandler((NXgetdata(file_id, muon.temperature_log_time)));
}

```

```

        NXMerrorhandler((NXMgetstringattr(file_id, "units",
            muon.temperature_log_time_units)));
        NXMerrorhandler((NXclosedata(file_id)));
        NXMdebug("%s\n", "Read time data");
        NXMdebug("Logging started at time %f\n", muon.temperature_log_time[0]);
    }
    else
    {
        muon.temperature_log_values = (float *) NULL;
        muon.temperature_log_time = (float *) NULL;
    }

    /* close subgroup 'NXlog' (temperature) */
    NXMerrorhandler((NXclosegroup(file_id)));
    NXMdebug("%s\n", "Close 'NXlog' (temperature)");

    /* open subgroup 'NXlog' (events) */
    NXMerrorhandler((NXopengroup(file_id, "event_log_1", "NXlog")));
    NXMdebug("%s\n", "Open 'NXlog' (events)...");

    /* read 'events name' */
    NXMerrorhandler((NXopendata(file_id, "name")));
    NXMerrorhandler((NXMgetstringdata(file_id, muon.events_log_name)));
    attlen = SIZE_INT32;
    atttype = NX_INT32;
    NXMerrorhandler((NXgetattr(file_id, "available", &muon.events_log_available,
        &attlen, &atttype)));
    NXMerrorhandler((NXclosedata(file_id)));
    NXMdebug("Logged name: %s, ", muon.events_log_name);
    NXMdebug("%d values available\n", muon.events_log_available);

    /* if 'available' then read ... */
    if (muon.events_log_available > 0)
    {
        /* 'events values' and 'units' attributes */
        NXMerrorhandler((NXopendata(file_id, "values")));
        muon.events_log_values = (float *) NXMemoryhandler(file_id);
        NXMerrorhandler((NXgetdata(file_id, muon.events_log_values)));
        NXMerrorhandler((NXMgetstringattr(file_id, "units",
            muon.events_log_values_units)));
        NXMerrorhandler((NXclosedata(file_id)));
        NXMdebug("First events value logged: %f\n", muon.events_log_values[0]);

        /* read 'events time' and 'units' attribute */
        NXMerrorhandler((NXopendata(file_id, "time")));
        muon.events_log_time = (float *) NXMemoryhandler(file_id);
        NXMerrorhandler((NXgetdata(file_id, muon.events_log_time)));
        NXMerrorhandler((NXMgetstringattr(file_id, "units",
            muon.events_log_time_units)));
        NXMerrorhandler((NXclosedata(file_id)));
        NXMdebug("%s\n", "Read time data");
        NXMdebug("Logging started at time %f\n", muon.events_log_time[0]);
    }
    else
    {
        muon.events_log_values = (float *) NULL;
        muon.events_log_time = (float *) NULL;
    }

    /* close subgroup 'NXlog' (events) */
    NXMerrorhandler((NXclosegroup(file_id)));
    NXMdebug("%s\n", "Close 'NXlog' (events)");

    /* close root group 'NXentry' */
    NXMerrorhandler((NXclosegroup(file_id)));
    NXMdebug("%s\n", "Close 'NXentry'");

    /* close NeXus file */
    NXMerrorhandler((NXclose (&file_id)));
    NXMdebug("%s\n", "Close NeXus file");

    return NX_OK;
}

```

```

/* NXMfreememory - free memory claimed in NeXus file */

```

```

void NXMfreememory(struct NXM_MUONIDF *nxfdata)
{
    /* NXlog (temperature) */
    NXMfree(nxfdata->temperature_log_values);
    NXMfree(nxfdata->temperature_log_time);
    /* NXlog (events) */
    NXMfree(nxfdata->events_log_values);
    NXMfree(nxfdata->events_log_time);
    /* NXdata */
    NXMfree(nxfdata->histogram_counts);
    NXMfree(nxfdata->histogram_raw_time);
    NXMfree(nxfdata->histogram_corrected_time);
    NXMfree(nxfdata->histogram_grouping);
    NXMfree(nxfdata->histogram_alpha);
    /* NXuser */
    NXMfree(nxfdata->user_uif_array);
    /* NXsample */
    NXMfree(nxfdata->sample_magnetic_field_vector);
    /* NXdetector */
    NXMfree(nxfdata->detector_angles);
    NXMfree(nxfdata->detector_deadtimes);
}

```

APPENDIX C: Test_nexus_reader.c

```
/*
 * Demonstration program using the NeXus ISIS Read Routine
 *
 *
 * Tested: 1) Visual C++ 5.0, Windows NT 4.0
 *         2) OpenVMS 7.2
 *
 *
 * D.W.Flannery and S.P.Cottrell
 */

#include "napi.h"
#include "nexus_reader.h"

void main()
{
    struct NXM_MUONIDF *muon_ptr;
    int *histogram_ptr;
    int i, j, histogram;
    char nxfname[80];

    /* input filename from keyboard */
    printf("Enter NeXus filename: ");
    scanf("%s", nxfname);
    printf("\nReading NeXus file: %s\n", nxfname);

    /* read NeXus file */
    if(NXMread(nxfname, &muon_ptr) != NX_OK)
    {
        printf("\nError reading NeXus file %s\n", nxfname);
        exit(0);
    }

    /* print contents of NeXus file */
    printf("\n\nNeXus File Contents\n\n");
    printf("Program name: %s\n", muon_ptr->run_program_name);
    printf("Program version: %s\n", muon_ptr->run_program_version);
    printf("Run number: %d\n", muon_ptr->run_number);
    printf("Title: %s\n", muon_ptr->run_title);
    printf("Notes: %s\n", muon_ptr->run_notes);
    printf("Analysis: %s\n", muon_ptr->run_analysis);
    printf("Start_time: %s\n", muon_ptr->run_start_time);
    printf("Stop_time: %s\n", muon_ptr->run_stop_time);
    printf("Duration: %s\n", muon_ptr->run_duration);
    printf("Switching States: %d\n", muon_ptr->run_switching_states);
    /* NXuser */
    printf("Username: %s\n", muon_ptr->user_name);
    printf("Experiment number: %s\n", muon_ptr->user_experiment_number);
    /* NXsample */
    printf("Sample: %s\n", muon_ptr->sample_name);
    printf("Temperature: %f %s\n", muon_ptr->sample_temperature, muon_ptr->
    >sample_temperature_units);
    printf("Magnetic field: %f %s\n", muon_ptr->sample_magnetic_field, muon_ptr->
    >sample_magnetic_field_units);
    printf("Shape: %s\n", muon_ptr->sample_shape);
    printf("Magnetic field state: %s\n", muon_ptr->sample_magnetic_field_state);
    if (muon_ptr->sample_magnetic_field_vector_available != 0)
    {
        printf("Magnetic field vector: [%f, %f, %f]\n", muon_ptr->sample_magnetic_field_vector[0],
            muon_ptr->sample_magnetic_field_vector[1],
            muon_ptr->sample_magnetic_field_vector[2]);
        printf("Coordinate system: %s\n", muon_ptr->sample_magnetic_field_vector_coord);
    }
    else
        printf("Magnetic field vector not available\n");
    printf("Sample environment: %s\n", muon_ptr->sample_environment);
    /* NXinstrument */
    printf("Instrument name: %s\n", muon_ptr->instrument_name);
}
```

```

/* NXdetector */
printf("Number of detectors: %d\n", muon_ptr->detector_number);
printf("Detector orientation: %s\n", muon_ptr->detector_orientation);
/* print detector angle information */
if (muon_ptr->detector_angles_available != 0)
{
    printf("Position of first detector: [%f, %f, %f]\n", muon_ptr->detector_angles[0],
        muon_ptr->detector_angles[1],
        muon_ptr->detector_angles[2]);
    printf("Coordinate system: %s\n", muon_ptr->detector_angles_coord);
    printf("Solid angle of 1st detector: %f\n", muon_ptr->detector_angles[3]);
}
else
    printf("Detector angles not available\n");
/* print detector deadtime information */
if (muon_ptr->detector_deadtimes_available != 0)
{
    printf("\nDeadtimes available (in %s) ...\n", muon_ptr->detector_deadtimes_units);
    for (i = 0; i < muon_ptr->detector_number; i++)
        printf("Deadtime of detector %d: %f\n", (i + 1), muon_ptr->detector_deadtimes[i]);
}
else
    printf("\nDeadtime information not available\n");

/* NXcollimator */
printf("Collimator type: %s\n", muon_ptr->collimator_type);
printf("Collimator aperture: %s\n", muon_ptr->collimator_aperture);
/* NXbeam */
printf("Events: %f %s\n", muon_ptr->beam_total_counts, muon_ptr->beam_total_counts_units);
printf("Frames: %d\n", muon_ptr->beam_frames);
printf("DAEreads: %d\n", muon_ptr->beam_daereads);
/* NXdata */
printf("Number of histograms: %d\n", muon_ptr->histogram_number);
printf("Time zero bin: %d\n", muon_ptr->histogram_t0_bin);
printf("Histogram length: %d\n", muon_ptr->histogram_length);
printf("First good bin: %d\n", muon_ptr->histogram_first_good_bin);
printf("Last good bin: %d\n", muon_ptr->histogram_last_good_bin);
printf("Histogram resolution: %d\n", muon_ptr->histogram_resolution);
printf("Histogram units: %s\n", muon_ptr->histogram_units);

/* print time zero information if available */
if (muon_ptr->histogram_time_zero_available > 0)
    printf("Histogram time zero: %f %s\n", muon_ptr->histogram_time_zero,
        muon_ptr->histogram_time_zero_units);

/* print grouping information if available */
if (muon_ptr->histogram_grouping_available > 0)
{
    printf("\nFound %d histogram groups\n", muon_ptr->histogram_grouping_available);
    for (i = 1; i <= muon_ptr->histogram_grouping_available; i++)
    {
        printf("Group %d histograms: ", i);
        for (j = 0; j < muon_ptr->histogram_number; j++)
            if ((*muon_ptr->histogram_grouping + j) == i)
                printf("%d, ", (j + 1));
        printf("\n");
    }
}
/* print alpha values for paired groups if available */
if (muon_ptr->histogram_alpha_available > 0)
{
    printf("alpha defined for %d pairs\n", muon_ptr->histogram_alpha_available);
    for (i = 0; i < muon_ptr->histogram_alpha_available; i++)
        printf("alpha for group pair %d (groups %f and %f): %f\n", (i + 1),
            *(muon_ptr->histogram_alpha + (i * 3)),
            *(muon_ptr->histogram_alpha + (i * 3) + 1),
            *(muon_ptr->histogram_alpha + (i * 3) + 2));
}

/* print User Information */
printf("\n");
for (i = 0; i < muon_ptr->user_uif_array_length; i++)
{
    printf("User Information %d:\n", i + 1);
    printf("    Data Name: %s\n", muon_ptr->user_uif_array[i].uif_element_name);
    printf("    Data Value: %s\n", muon_ptr->user_uif_array[i].uif_element_value);
}

```



```

/* demonstrate pointer manipulation of individual histograms ... print every 50th value
in the 1st histogram with time stamp */

printf("\nHistogram information - (Raw Time, Corrected Time, Counts) (time in %s):\n",
muon_ptr->histogram_corrected_time_units);
histogram = 0; /* choose which histogram you want (starting at zero) */
histogram_ptr = &(muon_ptr->histogram_counts[histogram * muon_ptr->histogram_length]);
for (i = 0; i < (muon_ptr->histogram_length); i+=50) /* extract data */
{
    printf("%f, %f, %d\n", *(muon_ptr->histogram_raw_time + i),
        *(muon_ptr->histogram_corrected_time + i),
        *(histogram_ptr+=50));
}
printf("\n");

/* access logged information ... temperature */
if (muon_ptr->temperature_log_available > 0)
{
    printf("\nTemperature Log - Time (%s), Temperature (%s)\n",
        muon_ptr->temperature_log_time_units,
        muon_ptr->temperature_log_values_units);
    for (i = 0; i < (muon_ptr->temperature_log_available); i++)
    {
        printf("(%f, %f)\n", muon_ptr->temperature_log_time[i],
            muon_ptr->temperature_log_values[i]);
    }
}
else
    printf("\nTemperature log not available\n");

/* access logged information ... events */
if (muon_ptr->events_log_available > 0)
{
    printf("\nEvent Log - Time (%s), Recorded events (%s)\n",
        muon_ptr->events_log_time_units, muon_ptr->events_log_values_units);
    for (i = 0; i < (muon_ptr->events_log_available); i++)
    {
        printf("(%f, %f)\n", muon_ptr->events_log_time[i],
            muon_ptr->events_log_values[i]);
    }
}
else
    printf("\nEvent log not available\n");

/* free memory claimed by read routine */
NXMfreememory(muon_ptr);
}

```


APPENDIX D: MUON_DEF_F77.INC

```
C      include file for NeXus_reader.For
C      written by Damian Flannery

C      ----- define constants -----

C      constant to define length of strings
      INTEGER NXMname
      PARAMETER (NXMname = 60)

C      switch for debugging information
      LOGICAL DEBUG
      PARAMETER (DEBUG = .TRUE.)

C      version number of the read routine
      INTEGER READVERSION
      PARAMETER (READVERSION = 1)

C      compatible version of instrument definition file
      INTEGER IDFVERSION
      PARAMETER (IDFVERSION = 1)

C      compatible analysis
      CHARACTER*(NXMname) ANALYSIS
      PARAMETER (ANALYSIS = 'muonTD')

C      compatible lab
      CHARACTER*(NXMname) LAB
      PARAMETER (LAB = 'ISIS')

C      uif array length
      INTEGER UIFLENGTH
      PARAMETER (UIFLENGTH = 30)

C      number of detectors
      INTEGER NUMDETECTORS
      PARAMETER (NUMDETECTORS = 32)

C      maximum number of histograms
      INTEGER MAXHISNUM
      PARAMETER (MAXHISNUM = 80)

C      maximum histogram length
      INTEGER MAXHISLEN
      PARAMETER (MAXHISLEN = 2048)

C      maximum log length
      INTEGER MAXLOGLEN
      PARAMETER (MAXLOGLEN = 80000)

C      ##### NXrun #####
      CHARACTER*(NXMname) NXM_run_program_name
      CHARACTER*(NXMname) NXM_run_program_version
      INTEGER NXM_run_idf_version
      INTEGER NXM_run_number
      CHARACTER*(NXMname) NXM_run_title
      CHARACTER*(NXMname) NXM_run_notes
      CHARACTER*(NXMname) NXM_run_analysis
      CHARACTER*(NXMname) NXM_run_lab
      CHARACTER*(NXMname) NXM_run_beamline
      CHARACTER*(NXMname) NXM_run_start_time
      CHARACTER*(NXMname) NXM_run_stop_time
      CHARACTER*(NXMname) NXM_run_duration
      INTEGER NXM_run_switching_states

C      ##### NXuser #####
      CHARACTER*(NXMname) NXM_user_name
      CHARACTER*(NXMname) NXM_user_experiment_number
      INTEGER NXM_user_uif_array_length
```

```

CHARACTER*(NXMname) NXM_user_uif_element_name (UIFLENGTH)
CHARACTER*(NXMname) NXM_user_uif_element_value (UIFLENGTH)

C      ##### NXsample #####
CHARACTER*(NXMname) NXM_sample_name
REAL NXM_sample_temperature
CHARACTER*(NXMname) NXM_sample_temperature_units
REAL NXM_sample_magnetic_field
CHARACTER*(NXMname) NXM_sample_magnetic_field_units
CHARACTER*(NXMname) NXM_sample_shape
CHARACTER*(NXMname) NXM_sample_magnetic_field_state
INTEGER NXM_sample_mfield_vec_available
REAL NXM_sample_mfield_vector(3)
CHARACTER*(NXMname) NXM_sample_mfield_vector_coord
CHARACTER*(NXMname) NXM_sample_environment

C      ##### NXinstrument #####
CHARACTER*(NXMname) NXM_instrument_name

C      ##### NXdetector #####
INTEGER NXM_detector_number
CHARACTER*(NXMname) NXM_detector_orientation
INTEGER NXM_detector_angles_available
REAL NXM_detector_angles (4, NUMDETECTORS)
CHARACTER*(NXMname) NXM_detector_angles_coord
INTEGER NXM_deadtimses_available
REAL NXM_detector_deadtimses (NUMDETECTORS)
CHARACTER*(NXMname) NXM_detector_deadtimses_units

C      ##### NXcollimator #####
CHARACTER*(NXMname) NXM_collimator_type
CHARACTER*(NXMname) NXM_collimator_aperture

C      ##### NXbeam #####
REAL NXM_beam_total_counts
CHARACTER*(NXMname) NXM_beam_total_counts_units
INTEGER NXM_beam_daereads
INTEGER NXM_beam_frames

C      ##### NXdata histogram_data_1 #####
INTEGER NXM_histogram_counts (MAXHISNUM, MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_counts_units
INTEGER NXM_histogram_number
INTEGER NXM_histogram_length
INTEGER NXM_histogram_t0_bin
INTEGER NXM_histogram_first_good_bin
INTEGER NXM_histogram_last_good_bin
INTEGER NXM_histogram_resolution
CHARACTER*(NXMname) NXM_histogram_resolution_units
REAL NXM_histogram_raw_time (MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_raw_time_units
REAL NXM_histogram_corrected_time (MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_ctime_units
INTEGER NXM_histogram_grp_available
INTEGER NXM_histogram_grouping (NUMDETECTORS)
INTEGER NXM_histogram_alpha_available
REAL NXM_histogram_alpha (3,4)
REAL NXM_histogram_time_zero
CHARACTER*(NXMname) NXM_histogram_time_zero_units
INTEGER NXM_histogram_t0_available

C      ##### NXlog temperature_log_1 #####
INTEGER NXM_temp_log_available
CHARACTER*(NXMname) NXM_temp_log_name
REAL NXM_temp_log_values (MAXLOGLEN)
CHARACTER*(NXMname) NXM_temp_log_values_units
REAL NXM_temp_log_time (MAXLOGLEN)
CHARACTER*(NXMname) NXM_temp_log_time_units

C      ##### NXlog event_log_1 #####
INTEGER NXM_event_log_available
CHARACTER*(NXMname) NXM_event_log_name
REAL NXM_event_log_values (MAXLOGLEN)
CHARACTER*(NXMname) NXM_event_log_values_units
REAL NXM_event_log_time (MAXLOGLEN)
CHARACTER*(NXMname) NXM_event_log_time_units

```

```

C      ----- define common blocks here -----

C      block for integer data

COMMON /NXM_integer_data_block/ NXM_run_idf_version, NXM_run_number,
&NXM_run_switching_states, NXM_user_uif_array_length,
&NXM_sample_mfield_vec_available, NXM_detector_number,
&NXM_detector_angles_available, NXM_deadtimes_available,
&NXM_beam_daereads, NXM_beam_frames, NXM_histogram_counts, NXM_histogram_number,
&NXM_histogram_length, NXM_histogram_t0_bin, NXM_histogram_first_good_bin,
&NXM_histogram_last_good_bin, NXM_histogram_resolution,
&NXM_histogram_grp_available, NXM_histogram_grouping,
&NXM_histogram_alpha_available, NXM_histogram_t0_available,
&NXM_temp_log_available, NXM_event_log_available

C      block for real data

COMMON /NXM_real_data_block/ NXM_sample_temperature, NXM_sample_magnetic_field,
&NXM_sample_mfield_vector, NXM_detector_angles, NXM_detector_deadtimes,
&NXM_beam_total_counts, NXM_histogram_raw_time, NXM_histogram_corrected_time,
&NXM_histogram_alpha, NXM_histogram_time_zero, NXM_temp_log_values,
&NXM_temp_log_time, NXM_event_log_values, NXM_event_log_time

C      block for character data

COMMON /NXM_character_data_block/ NXM_run_program_name, NXM_run_program_version,
&NXM_run_title, NXM_run_notes, NXM_run_analysis, NXM_run_lab, NXM_run_beamline,
&NXM_run_start_time, NXM_run_stop_time, NXM_run_duration, NXM_user_name,
&NXM_user_experiment_number, NXM_user_uif_element_name, NXM_user_uif_element_value,
&NXM_sample_name, NXM_sample_temperature_units, NXM_sample_magnetic_field_units,
&NXM_sample_shape, NXM_sample_magnetic_field_state,
&NXM_sample_mfield_vector_coord, NXM_sample_environment, NXM_instrument_name,
&NXM_detector_orientation, NXM_detector_angles_coord, NXM_detector_deadtimes_units,
&NXM_collimator_type, NXM_collimator_aperture, NXM_beam_total_counts_units,
&NXM_histogram_counts_units, NXM_histogram_resolution_units,
&NXM_histogram_raw_time_units, NXM_histogram_ctime_units,
&NXM_histogram_time_zero_units, NXM_temp_log_name, NXM_temp_log_values_units,
&NXM_temp_log_time_units, NXM_event_log_name, NXM_event_log_values_units,
&NXM_event_log_time_units

```


APPENDIX E: MUON_DEF_F90.INC

```
!      F90 include file for NeXus_reader.for
!      written by Damian Flannery

!      ----- define constants -----

!      constant to define length of strings
      INTEGER NXMname
      PARAMETER (NXMname = 60)

!      switch for debugging information
      LOGICAL DEBUG
      PARAMETER (DEBUG = .TRUE.)

!      version number of the read routine
      INTEGER READVERSION
      PARAMETER (READVERSION = 1)

!      compatible version of instrument definition file
      INTEGER IDFVERSION
      PARAMETER (IDFVERSION = 1)

!      compatible analysis
      CHARACTER*(NXMname) ANALYSIS
      PARAMETER (ANALYSIS = 'muonTD')

!      compatible lab
      CHARACTER*(NXMname) LAB
      PARAMETER (LAB = 'ISIS')

!      uif array length
      INTEGER UIFLENGTH
      PARAMETER (UIFLENGTH = 30)

!      number of detectors
      INTEGER NUMDETECTORS
      PARAMETER (NUMDETECTORS = 32)

!      maximum number of histograms
      INTEGER MAXHISNUM
      PARAMETER (MAXHISNUM = 80)

!      maximum histogram length
      INTEGER MAXHISLEN
      PARAMETER (MAXHISLEN = 2048)

!      maximum log length
      INTEGER MAXLOGLEN
      PARAMETER (MAXLOGLEN = 80000)

!      ##### NXrun #####
      CHARACTER*(NXMname) NXM_run_program_name
      CHARACTER*(NXMname) NXM_run_program_version
      INTEGER NXM_run_idf_version
      INTEGER NXM_run_number
      CHARACTER*(NXMname) NXM_run_title
      CHARACTER*(NXMname) NXM_run_notes
      CHARACTER*(NXMname) NXM_run_analysis
      CHARACTER*(NXMname) NXM_run_lab
      CHARACTER*(NXMname) NXM_run_beamline
      CHARACTER*(NXMname) NXM_run_start_time
      CHARACTER*(NXMname) NXM_run_stop_time
      CHARACTER*(NXMname) NXM_run_duration
      INTEGER NXM_run_switching_states

!      ##### NXuser #####
      CHARACTER*(NXMname) NXM_user_name
      CHARACTER*(NXMname) NXM_user_experiment_number
      INTEGER NXM_user_uif_array_length
      CHARACTER*(NXMname) NXM_user_uif_element_name (UIFLENGTH)
```

```

CHARACTER*(NXMname) NXM_user_uif_element_value (UIFLENGTH)

! ##### NXsample #####
CHARACTER*(NXMname) NXM_sample_name
REAL NXM_sample_temperature
CHARACTER*(NXMname) NXM_sample_temperature_units
REAL NXM_sample_magnetic_field
CHARACTER*(NXMname) NXM_sample_magnetic_field_units
CHARACTER*(NXMname) NXM_sample_shape
CHARACTER*(NXMname) NXM_sample_magnetic_field_state
INTEGER NXM_sample_mfield_vec_available
REAL NXM_sample_mfield_vector(3)
CHARACTER*(NXMname) NXM_sample_mfield_vector_coord
CHARACTER*(NXMname) NXM_sample_environment

! ##### NXinstrument #####
CHARACTER*(NXMname) NXM_instrument_name

! ##### NXdetector #####
INTEGER NXM_detector_number
CHARACTER*(NXMname) NXM_detector_orientation
INTEGER NXM_detector_angles_available
REAL NXM_detector_angles (4, NUMDETECTORS)
CHARACTER*(NXMname) NXM_detector_angles_coord
INTEGER NXM_deadtimes_available
REAL NXM_detector_deadtimes (NUMDETECTORS)
CHARACTER*(NXMname) NXM_detector_deadtimes_units

! ##### NXcollimator #####
CHARACTER*(NXMname) NXM_collimator_type
CHARACTER*(NXMname) NXM_collimator_aperture

! ##### NXbeam #####
REAL NXM_beam_total_counts
CHARACTER*(NXMname) NXM_beam_total_counts_units
INTEGER NXM_beam_daereads
INTEGER NXM_beam_frames

! ##### NXdata histogram_data_1 #####
INTEGER NXM_histogram_counts (MAXHISNUM, MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_counts_units
INTEGER NXM_histogram_number
INTEGER NXM_histogram_length
INTEGER NXM_histogram_t0_bin
INTEGER NXM_histogram_first_good_bin
INTEGER NXM_histogram_last_good_bin
INTEGER NXM_histogram_resolution
CHARACTER*(NXMname) NXM_histogram_resolution_units
REAL NXM_histogram_raw_time (MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_raw_time_units
REAL NXM_histogram_corrected_time (MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_ctime_units
INTEGER NXM_histogram_grp_available
INTEGER NXM_histogram_grouping (NUMDETECTORS)
INTEGER NXM_histogram_alpha_available
REAL NXM_histogram_alpha (3,4)
REAL NXM_histogram_time_zero
CHARACTER*(NXMname) NXM_histogram_time_zero_units
INTEGER NXM_histogram_t0_available

! ##### NXlog temperature_log_1 #####
INTEGER NXM_temp_log_available
CHARACTER*(NXMname) NXM_temp_log_name
REAL NXM_temp_log_values (MAXLOGLEN)
CHARACTER*(NXMname) NXM_temp_log_values_units
REAL NXM_temp_log_time (MAXLOGLEN)
CHARACTER*(NXMname) NXM_temp_log_time_units

! ##### NXlog event_log_1 #####
INTEGER NXM_event_log_available
CHARACTER*(NXMname) NXM_event_log_name
REAL NXM_event_log_values (MAXLOGLEN)
CHARACTER*(NXMname) NXM_event_log_values_units
REAL NXM_event_log_time (MAXLOGLEN)
CHARACTER*(NXMname) NXM_event_log_time_units

```



```

!      ----- define common blocks here -----
!
!      block for integer data

      COMMON /NXM_integer_data_block/
&NXM_run_idf_version,NXM_run_number,NXM_run_switching_states,&
&NXM_user_uif_array_length,NXM_sample_mfield_vec_available,&
&NXM_detector_number,&NXM_detector_angles_available,NXM_deadtimes_available,&
&NXM_beam_daereads, NXM_beam_frames, NXM_histogram_counts, NXM_histogram_number,&
&NXM_histogram_length,NXM_histogram_t0_bin,NXM_histogram_first_good_bin,&
&NXM_histogram_last_good_bin, NXM_histogram_resolution,&
&NXM_histogram_grp_available,NXM_histogram_grouping,&
&NXM_histogram_alpha_available,NXM_histogram_t0_available,&
&NXM_temp_log_available, NXM_event_log_available

!      block for real data

      COMMON /NXM_real_data_block/ NXM_sample_temperature,&
&NXM_sample_magnetic_field,NXM_sample_mfield_vector,&
&NXM_detector_angles,NXM_detector_deadtimes,NXM_beam_total_counts,&
&NXM_histogram_raw_time,NXM_histogram_corrected_time,NXM_histogram_alpha,&
&NXM_histogram_time_zero,NXM_temp_log_values,NXM_temp_log_time,&
&NXM_event_log_values, NXM_event_log_time

!      block for character data

      COMMON /NXM_character_data_block/ NXM_run_program_name,&
&NXM_run_program_version, NXM_run_title, NXM_run_notes,&
&NXM_run_analysis,NXM_run_lab,NXM_run_beamline,NXM_run_start_time,&
&NXM_run_stop_time, NXM_run_duration,NXM_user_name,&
&NXM_user_experiment_number, NXM_user_uif_element_name,&
&NXM_user_uif_element_value, NXM_sample_name,&
&NXM_sample_temperature_units, NXM_sample_magnetic_field_units,&
&NXM_sample_shape, NXM_sample_magnetic_field_state,&
&NXM_sample_mfield_vector_coord, NXM_sample_environment, NXM_instrument_name,&
&NXM_detector_orientation, NXM_detector_angles_coord,&
&NXM_detector_deadtimes_units, NXM_collimator_type, NXM_collimator_aperture,&
&NXM_beam_total_counts_units, NXM_histogram_resolution_units,&
&NXM_histogram_raw_time_units, NXM_histogram_ctime_units,&
&NXM_histogram_time_zero_units,NXM_temp_log_name,NXM_temp_log_values_units,&
&NXM_temp_log_time_units, NXM_event_log_name,&
&NXM_event_log_values_units, NXM_event_log_time_units

```


APPENDIX F: NEXUS_READER.FOR

```

        INTEGER FUNCTION NXMread (nxfname)

        include 'NAPIF.INC'

C      ----- define constants -----
C
C      constant to define length of strings
        INTEGER NXMname
        PARAMETER (NXMname = 60)
C
C      switch for debugging information
        LOGICAL DEBUG
        PARAMETER (DEBUG = .FALSE.)
C
C      version number of the read routine
        INTEGER READVERSION
        PARAMETER (READVERSION = 1)
C
C      compatible version of instrument definition file
        INTEGER IDVERSION
        PARAMETER (IDVERSION = 1)
C
C      compatible analysis
        CHARACTER*(NXMname) ANALYSIS
        PARAMETER (ANALYSIS = 'muonTD')
C
C      compatible lab
        CHARACTER*(NXMname) LAB
        PARAMETER (LAB = 'ISIS')
C
C      uif array length
        INTEGER UIFLENGTH
        PARAMETER (UIFLENGTH = 30)
C
C      number of detectors
        INTEGER NUMDETECTORS
        PARAMETER (NUMDETECTORS = 32)
C
C      maximum number of histograms
        INTEGER MAXHISNUM
        PARAMETER (MAXHISNUM = 80)
C
C      maximum histogram length
        INTEGER MAXHISLEN
        PARAMETER (MAXHISLEN = 2048)
C
C      maximum log length
        INTEGER MAXLOGLEN
        PARAMETER (MAXLOGLEN = 80000)

C      ----- define NeXus elements -----
C
C      ##### NXrun #####
        CHARACTER*(NXMname) NXM_run_program_name
        CHARACTER*(NXMname) NXM_run_program_version
        INTEGER NXM_run_idf_version
        INTEGER NXM_run_number
        CHARACTER*(NXMname) NXM_run_title
        CHARACTER*(NXMname) NXM_run_notes
        CHARACTER*(NXMname) NXM_run_analysis
        CHARACTER*(NXMname) NXM_run_lab
        CHARACTER*(NXMname) NXM_run_beamline
        CHARACTER*(NXMname) NXM_run_start_time
        CHARACTER*(NXMname) NXM_run_stop_time
        CHARACTER*(NXMname) NXM_run_duration
        INTEGER NXM_run_switching_states
C
C      ##### NXuser #####
        CHARACTER*(NXMname) NXM_user_name
        CHARACTER*(NXMname) NXM_user_experiment_number
        INTEGER NXM_user_uif_array_length
        CHARACTER*(NXMname) NXM_user_uif_element_name (UIFLENGTH)
        CHARACTER*(NXMname) NXM_user_uif_element_value (UIFLENGTH)
C
C      ##### NXsample #####
        CHARACTER*(NXMname) NXM_sample_name
        REAL NXM_sample_temperature
        CHARACTER*(NXMname) NXM_sample_temperature_units

```

```

REAL NXM_sample_magnetic_field
CHARACTER*(NXMname) NXM_sample_magnetic_field_units
CHARACTER*(NXMname) NXM_sample_shape
CHARACTER*(NXMname) NXM_sample_magnetic_field_state
INTEGER NXM_sample_mfield_vec_available
REAL NXM_sample_mfield_vector(3)
CHARACTER*(NXMname) NXM_sample_mfield_vector_coord
CHARACTER*(NXMname) NXM_sample_environment

C ##### NXinstrument #####
CHARACTER*(NXMname) NXM_instrument_name

C ##### NXdetector #####
INTEGER NXM_detector_number
CHARACTER*(NXMname) NXM_detector_orientation
INTEGER NXM_detector_angles_available
REAL NXM_detector_angles (4, NUMDETECTORS)
CHARACTER*(NXMname) NXM_detector_angles_coord
INTEGER NXM_deadtimses_available
REAL NXM_detector_deadtimses (NUMDETECTORS)
CHARACTER*(NXMname) NXM_detector_deadtimses_units

C ##### NXcollimator #####
CHARACTER*(NXMname) NXM_collimator_type
CHARACTER*(NXMname) NXM_collimator_aperture

C ##### NXbeam #####
REAL NXM_beam_total_counts
CHARACTER*(NXMname) NXM_beam_total_counts_units
INTEGER NXM_beam_daereads
INTEGER NXM_beam_frames

C ##### NXdata histogram_data_1 #####
INTEGER NXM_histogram_counts (MAXHISNUM, MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_counts_units
INTEGER NXM_histogram_number
INTEGER NXM_histogram_length
INTEGER NXM_histogram_t0_bin
INTEGER NXM_histogram_first_good_bin
INTEGER NXM_histogram_last_good_bin
INTEGER NXM_histogram_resolution
CHARACTER*(NXMname) NXM_histogram_resolution_units
REAL NXM_histogram_raw_time (MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_raw_time_units
REAL NXM_histogram_corrected_time (MAXHISLEN)
CHARACTER*(NXMname) NXM_histogram_ctype_units
INTEGER NXM_histogram_grp_available
INTEGER NXM_histogram_grouping (NUMDETECTORS)
INTEGER NXM_histogram_alpha_available
REAL NXM_histogram_alpha (3,4)
REAL NXM_histogram_time_zero
CHARACTER*(NXMname) NXM_histogram_time_zero_units
INTEGER NXM_histogram_t0_available

C ##### NXlog temperature_log_1 #####
INTEGER NXM_temp_log_available
CHARACTER*(NXMname) NXM_temp_log_name
REAL NXM_temp_log_values (MAXLOGLEN)
CHARACTER*(NXMname) NXM_temp_log_values_units
REAL NXM_temp_log_time (MAXLOGLEN)
CHARACTER*(NXMname) NXM_temp_log_time_units

C ##### NXlog event_log_1 #####
INTEGER NXM_event_log_available
CHARACTER*(NXMname) NXM_event_log_name
REAL NXM_event_log_values (MAXLOGLEN)
CHARACTER*(NXMname) NXM_event_log_values_units
REAL NXM_event_log_time (MAXLOGLEN)
CHARACTER*(NXMname) NXM_event_log_time_units

C ----- define local variables here -----

INTEGER file_id (NXHANDLESIZE)
CHARACTER*(NXMname) nxfname

INTEGER NLEN
INTEGER XRANK
INTEGER XDIM
INTEGER XTYPE
INTEGER ITEMN
INTEGER I, J, K, ATTLEN, RGNUM
CHARACTER*20 GNAME
CHARACTER*20 GCLASS
CHARACTER*8 TESTX

INTEGER NTINT, NTYPE
INTEGER NCOUNT

```

```

INTEGER numitems, datatype, uif_count
CHARACTER*(NXMname) groupname, classname, dataname

INTEGER temp_data (163840)

! ----- define common blocks here -----
!
! block for integer data

COMMON /NXM_integer_data_block/
&NXM_run_idf_version, NXM_run_number, NXM_run_switching_states, &
&NXM_user_uif_array_length, NXM_sample_mfield_vec_available, &
&NXM_detector_number, &NXM_detector_angles_available, NXM_deadtimses_available, &
&NXM_beam_daereads, NXM_beam_frames, NXM_histogram_counts, NXM_histogram_number, &
&NXM_histogram_length, NXM_histogram_t0_bin, NXM_histogram_first_good_bin, &
&NXM_histogram_last_good_bin, NXM_histogram_resolution, &
&NXM_histogram_grp_available, NXM_histogram_grouping, &
&NXM_histogram_alpha_available, NXM_histogram_t0_available, &
&NXM_temp_log_available, NXM_event_log_available

! block for real data

COMMON /NXM_real_data_block/ NXM_sample_temperature, &
&NXM_sample_magnetic_field, NXM_sample_mfield_vector, &
&NXM_detector_angles, NXM_detector_deadtimses, NXM_beam_total_counts, &
&NXM_histogram_raw_time, NXM_histogram_corrected_time, NXM_histogram_alpha, &
&NXM_histogram_time_zero, NXM_temp_log_values, NXM_temp_log_time, &
&NXM_event_log_values, NXM_event_log_time

! block for character data

COMMON /NXM_character_data_block/ NXM_run_program_name, &
&NXM_run_program_version, NXM_run_title, NXM_run_notes, &
&NXM_run_analysis, NXM_run_lab, NXM_run_beamline, NXM_run_start_time, &
&NXM_run_stop_time, NXM_run_duration, NXM_user_name, &
&NXM_user_experiment_number, NXM_user_uif_element_name, &
&NXM_user_uif_element_value, NXM_sample_name, &
&NXM_sample_temperature_units, NXM_sample_magnetic_field_units, &
&NXM_sample_shape, NXM_sample_magnetic_field_state, &
&NXM_sample_mfield_vector_coord, NXM_sample_environment, NXM_instrument_name, &
&NXM_detector_orientation, NXM_detector_angles_coord, &
&NXM_detector_deadtimses_units, NXM_collimator_type, NXM_collimator_aperture, &
&NXM_beam_total_counts_units, NXM_histogram_resolution_units, &
&NXM_histogram_raw_time_units, NXM_histogram_ctime_units, &
&NXM_histogram_time_zero_units, NXM_temp_log_name, NXM_temp_log_values_units, &
&NXM_temp_log_time_units, NXM_event_log_name, &
&NXM_event_log_values_units, NXM_event_log_time_units

C @@@@@@@@@@@@@@@@@@ Open Nexus File @@@@@@@@@@@@@@@@@@
if (NXopen (nxfname, NXACC_READ, file_id) .NE. NX_OK) stop
if (debug) print *, 'Opening NeXus file...'

C @@@@@@@@@@@@@@@@@@ Open 'Run' Group @@@@@@@@@@@@@@@@@@
if (NXopengroup (file_id, 'run', 'NXentry') .NE. NX_OK) stop
if (debug) print *, 'Opening "NXentry"...'

C display (in debug mode) Nexus reader version and IDF compatible version number
if (debug) then
  print *, ' '
  print *, 'Nexus reader version ', READVERSION
  print *, 'Compatible with ISIS Muon Instrument Definition version ', IDFVERSION
  print *, ' '
end if

C read idf version
if (NXopendata (file_id, 'idf_version') .NE. NX_OK) stop
if (NXgetdata (file_id, NXM_run_idf_version) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'IDF version : ', NXM_run_idf_version

C read analysis
if (NXopendata (file_id, 'analysis') .NE. NX_OK) stop
if (NXgetchardata (file_id, NXM_run_analysis) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Analysis : ', NXM_run_analysis

```

```

C      read lab
      if (NXopendata (file_id, 'lab') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_run_lab) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Lab : ', NXM_run_lab

C      check reader compatibility
      if ((NXM_run_idf_version .NE. IDFVERSION) .OR. ((NXM_run_analysis .NE. ANALYSIS) .OR. (NXM_run_lab .NE. LAB))) then
        print *, ' '
        print *, 'Fatal Error : Read routine incompatible with instrument definition used to write data'
        print *, 'Found IDF version: ',NXM_run_idf_version,', Analysis: ',NXM_run_analysis,', Lab: ', NXM_run_lab
        print *, 'Expecting IDF version: ',IDFVERSION,', Analysis: ',ANALYSIS,', Lab: ',LAB
        print *, ' '
        stop
      end if

C      read beamline
      if (NXopendata (file_id, 'beamline') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_run_beamline) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Beamline : ', NXM_run_beamline

      ATTLEN=3
      datatype = NX_CHAR
C      read program name and program version
      if (NXopendata (file_id, 'program_name') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_run_program_name) .NE. NX_OK) stop
      if (NXgetcharattr (file_id, 'version', NXM_run_program_version, ATTLEN, datatype) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Program Name : ', NXM_run_program_name
      if (debug) print *, 'Program Version : ', NXM_run_program_version

C      read run number
      if (NXopendata (file_id, 'number') .NE. NX_OK) stop
      if (NXgetdata (file_id, NXM_run_number) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Run Number : ', NXM_run_number

C      read title
      if (NXopendata (file_id, 'title') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_run_title) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Title : ', NXM_run_title

C      read notes
      if (NXopendata (file_id, 'notes') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_run_notes) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Notes : ', NXM_run_notes

C      read start_time
      if (NXopendata (file_id, 'start_time') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_run_start_time) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Start Time : ', NXM_run_start_time

C      read stop_time
      if (NXopendata (file_id, 'stop_time') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_run_stop_time) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Stop Time : ', NXM_run_stop_time

C      read duration
      if (NXopendata (file_id, 'duration') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_run_duration) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Duration : ', NXM_run_duration

C      read switching states (rgmode0
      if (NXopendata (file_id, 'switching_states') .NE. NX_OK) stop
      if (NXgetdata (file_id, NXM_run_switching_states) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Switching States : ', NXM_run_switching_states

C      @@@@@@@@@@@@@@@ open NXuser group @@@@@@@@@@@@@@@@@@
      if (NXopengroup (file_id, 'user', 'NXuser') .NE. NX_OK) stop
      if (debug) print *, ' '
      if (debug) print *, 'Opening NXuser...'

      if (NXinitgroupdir (file_id) .NE. NX_OK) stop
      if (NXgetgroupinfo (file_id, numitems, groupname, classname) .NE. NX_OK) stop
C      if (debug) print *, 'Group Info: ', trim(classname), ', ', trim(groupname), ', ', numitems

      uif_count = 0
      do i=1, numitems

```

```

if (NXgetnextentry (file_id, dataname, classname, datatype) .NE. NX_OK) stop
if (NXopendata(file_id, dataname) .NE. NX_OK) stop
if (debug) print *, 'Data name: ', dataname

if (dataname.EQ.'name') then
  if (NXgetchardata(file_id, NXM_user_name) .NE. NX_OK) stop
  if (debug) print *, 'Data value : ', NXM_user_name
else if (dataname.EQ.'experiment_number') then
  if (NXgetchardata (file_id, NXM_user_experiment_number) .NE. NX_OK) stop
  if (debug) print *, 'Data value : ', NXM_user_experiment_number
else
  uif_count = uif_count + 1
  if (NXgetchardata (file_id, NXM_user_uif_element_value (uif_count)) .NE. NX_OK) stop
  NXM_user_uif_element_name (uif_count) = dataname
  if (debug) print *, 'Data value : ', NXM_user_uif_element_value(uif_count)
end if
end do

NXM_user_uif_array_length = uif_count
if (debug) print *, 'UIF array length : ', NXM_user_uif_array_length

C      @@@@@@@@@@@@@@ close NXuser group @@@@@@@@@@@@@@@@@@
if (NXclosegroup (file_id) .NE. NX_OK) stop
if (debug) print *, 'Closing NXuser...'

C      @@@@@@@@@@@@@@ open NXsample group @@@@@@@@@@@@@@@@@@
if (NXopengroup (file_id, 'sample', 'NXSample') .NE. NX_OK) stop
if (debug) print *, ' '
if (debug) print *, 'Opening NXsample...'

C      read sample name
if (NXopendata (file_id, 'name') .NE. NX_OK) stop
if (NXgetchardata (file_id, NXM_sample_name) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Sample name : ', NXM_sample_name

ATTLEN=7
datatype = NX_CHAR
C      read temperature and units
if (NXopendata (file_id, 'temperature') .NE. NX_OK) stop
if (NXgetcharattr (file_id, 'units', NXM_sample_temperature_units, ATTLEN, datatype) .NE. NX_OK) stop
if (NXgetdata (file_id, NXM_sample_temperature) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Sample temperature : ', NXM_sample_temperature
if (debug) print *, 'Temperature units : ', NXM_sample_temperature_units

ATTLEN=6
datatype = NX_CHAR
C      read magnetic field and units
if (NXopendata (file_id, 'magnetic_field') .NE. NX_OK) stop
if (NXgetcharattr (file_id, 'units', NXM_sample_magnetic_field_units, ATTLEN, datatype) .NE. NX_OK) stop
if (NXgetdata (file_id, NXM_sample_magnetic_field) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Sample magnetic field : ', NXM_sample_magnetic_field
if (debug) print *, 'Field units : ', NXM_sample_magnetic_field_units

C      read shape
if (NXopendata (file_id, 'shape') .NE. NX_OK) stop
if (NXgetchardata (file_id, NXM_sample_shape) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Sample shape : ', NXM_sample_shape

C      read magnetic field state
if (NXopendata (file_id, 'magnetic_field_state') .NE. NX_OK) stop
if (NXgetchardata (file_id, NXM_sample_magnetic_field_state) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Field state : ', NXM_sample_magnetic_field

ATTLEN=1
datatype = NX_INT32
C      read magnetic field vector and coordinate system attribute
if (NXopendata (file_id, 'magnetic_field_vector') .NE. NX_OK) stop
if (NXgetattr (file_id, 'available', NXM_sample_mfield_vec_available, ATTLEN, datatype) .NE. NX_OK) stop
if (NXM_sample_mfield_vec_available.NE.0) then
  ATTLEN=20
  datatype = NX_CHAR
  if (NXgetcharattr (file_id, 'coordinate_system', NXM_sample_mfield_vector_coord, ATTLEN, datatype) .NE. NX_OK)
&      stop
  if (NXgetdata (file_id, NXM_sample_mfield_vector) .NE. NX_OK) stop

  if (debug) print *, 'Field units : [', NXM_sample_mfield_vector(1), ', ', NXM_sample_mfield_vector(2)
&      ', ', NXM_sample_mfield_vector(3), ']'
  if (debug) print *, 'Field Vector coordinate system : ', NXM_sample_mfield_vector_coord
else

```

```

        if (debug) print *, 'Magnetic field vector not available'
    end if
    if (NXclosedata (file_id) .NE. NX_OK) stop

C    read environment
    if (NXopendata (file_id, 'environment') .NE. NX_OK) stop
    if (NXgetchardata (file_id, NXM_sample_environment) .NE. NX_OK) stop
    if (NXclosedata (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Sample environment : ', NXM_sample_environment

C    ##### close NXsample group #####
    if (NXclosegroup (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Closing NXsample'

C    ##### open NXinstrument group #####
    if (NXopengroup (file_id, 'instrument', 'NXinstrument') .NE. NX_OK) stop
    if (debug) print *, ' '
    if (debug) print *, 'Opening NXinstrument...'

C    read instrument name
    if (NXopendata (file_id, 'name') .NE. NX_OK) stop
    if (NXgetchardata (file_id, NXM_instrument_name) .NE. NX_OK) stop
    if (NXclosedata (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Instrument name : ', NXM_instrument_name

C    ##### open NXdetector group #####
    if (NXopengroup (file_id, 'detector', 'NXdetector') .NE. NX_OK) stop
    if (debug) print *, ' '
    if (debug) print *, 'Opening NXdetector...'

C    read detector number
    if (NXopendata (file_id, 'number') .NE. NX_OK) stop
    if (NXgetdata (file_id, NXM_detector_number) .NE. NX_OK) stop
    if (NXclosedata (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Detector number : ', NXM_detector_number

C    read detector orientation
    if (NXopendata (file_id, 'orientation') .NE. NX_OK) stop
    if (NXgetchardata (file_id, NXM_detector_orientation) .NE. NX_OK) stop
    if (NXclosedata (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Detector orientation : ', NXM_detector_orientation

C    read angles and coordinate system attribute
    if (NXopendata (file_id, 'angles') .NE. NX_OK) stop
    ATTLEN=1
    datatype = NX_INT32
    if (NXgetattr (file_id, 'available', NXM_detector_angles_available, ATTLEN, datatype) .NE. NX_OK) stop
    if (NXM_detector_angles_available .NE. 0) then
        if (NXgetdata (file_id, NXM_detector_angles) .NE. NX_OK) stop
        if (debug) print *, 'Reading detector angles...'
        ATTLEN=20
        datatype = NX_CHAR
        if (NXgetcharattr (file_id, 'coordinate_system', NXM_detector_angles_coord, ATTLEN, datatype) .NE. NX_OK)
&        stop
        if (debug) print *, 'Detector angles coordinate system : ', NXM_detector_angles_coord
    else
        if (debug) print *, 'Detector angles not available'
    end if
    if (NXclosedata (file_id) .NE. NX_OK) stop

C    read deadtimes and coordinate system attribute
    if (NXopendata (file_id, 'deadtimes') .NE. NX_OK) stop
    ATTLEN=1
    datatype = NX_INT32
    if (NXgetattr (file_id, 'available', NXM_deadtmes_available, ATTLEN, datatype) .NE. NX_OK) stop
    if (NXM_deadtmes_available .NE. 0) then
        if (debug) print *, 'Reading detector deadtimes...'
        if (NXgetdata (file_id, NXM_detector_deadtmes) .NE. NX_OK) stop
        ATTLEN=20
        datatype = NX_CHAR
        if (NXgetcharattr (file_id, 'units', NXM_detector_deadtmes_units, ATTLEN, datatype) .NE. NX_OK) stop
        if (debug) print *, 'Detector deadtime units : ', NXM_detector_deadtmes_units
    else
        if (debug) print *, 'Detector deadtimes not available'
    end if
    if (NXclosedata (file_id) .NE. NX_OK) stop

C    ##### close NXdetector group #####
    if (NXclosegroup (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Closing NXdetector...'

```



```

C      @@@@@@@@@@@@@@ open NXcollimator group @@@@@@@@@@@@@@@@@@
      if (NXopengroup (file_id, 'collimator', 'NXcollimator') .NE. NX_OK) stop
      if (debug) print *, ' '
      if (debug) print *, 'Opening NXcollimator...'

C      read collimator type
      if (NXopendata (file_id, 'type') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_collimator_type) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Collimator type : ', NXM_collimator_type

C      read collimator aperture
      if (NXopendata (file_id, 'aperture') .NE. NX_OK) stop
      if (NXgetchardata (file_id, NXM_collimator_aperture) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Collimator aperture : ', NXM_collimator_aperture

C      @@@@@@@@@@@@@@ close NXcollimator group @@@@@@@@@@@@@@@@@@
      if (NXclosegroup (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Closing NXcollimator...'

C      @@@@@@@@@@@@@@ open NXbeam group @@@@@@@@@@@@@@@@@@
      if (NXopengroup (file_id, 'beam', 'NXbeam') .NE. NX_OK) stop
      if (debug) print *, ' '
      if (debug) print *, 'Opening NXbeam...'

C      read total counts and units
      if (NXopendata (file_id, 'total_counts') .NE. NX_OK) stop
      if (NXgetdata (file_id, NXM_beam_total_counts) .NE. NX_OK) stop
      ATTLEN=4
      datatype = NX_CHAR
      if (NXgetcharattr (file_id, 'units', NXM_beam_total_counts_units, ATTLEN, datatype) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Beam counts : ', NXM_beam_total_counts

C      read daereads
      if (NXopendata (file_id, 'daereads') .NE. NX_OK) stop
      if (NXgetdata (file_id, NXM_beam_daereads) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Beam daereads : ', NXM_beam_daereads

C      read frames
      if (NXopendata (file_id, 'frames') .NE. NX_OK) stop
      if (NXgetdata (file_id, NXM_beam_frames) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Beam frames : ', NXM_beam_frames

C      @@@@@@@@@@@@@@ close NXbeam group @@@@@@@@@@@@@@@@@@
      if (NXclosegroup (file_id) .NE. NX_OK) stop
      if (debug) print *, 'Closing NXbeam...'

C      @@@@@@@@@@@@@@ close NXinstrument group @@@@@@@@@@@@@@@@@@
      if (NXclosegroup (file_id) .NE. NX_OK) stop
      if (debug) print *, ' '
      if (debug) print *, 'Closing NXinstrument...'

C      @@@@@@@@@@@@@@ open NXdata group @@@@@@@@@@@@@@@@@@
      if (NXopengroup (file_id, 'histogram_data_1', 'NXdata') .NE. NX_OK) stop
      if (debug) print *, ' '
      if (debug) print *, 'Opening NXdata...'

C      read counts
      if (NXopendata (file_id, 'counts') .NE. NX_OK) stop
      if (debug) print *, 'Reading histogram_counts'
      if (NXgetdata (file_id, temp_data) .NE. NX_OK) stop

      ATTLEN=1
      datatype = NX_INT32
      if (NXgetattr (file_id, 'number', NXM_histogram_number, ATTLEN, datatype) .NE. NX_OK) stop
      if (NXgetattr (file_id, 'length', NXM_histogram_length, ATTLEN, datatype) .NE. NX_OK) stop
      if (NXgetattr (file_id, 't0_bin', NXM_histogram_t0_bin, ATTLEN, datatype) .NE. NX_OK) stop
      if (NXgetattr (file_id, 'first_good_bin', NXM_histogram_first_good_bin, ATTLEN, datatype) .NE. NX_OK) stop
      if (NXgetattr (file_id, 'last_good_bin', NXM_histogram_last_good_bin, ATTLEN, datatype) .NE. NX_OK) stop
      ATTLEN=12
      datatype = NX_CHAR
      if (NXgetcharattr (file_id, 'units', NXM_histogram_counts_units, ATTLEN, datatype) .NE. NX_OK) stop
      if (NXclosedata (file_id) .NE. NX_OK) stop

      if (debug) print *, 'Histogram number : ', NXM_histogram_number
      if (debug) print *, 'Histogram length : ', NXM_histogram_length
      if (debug) print *, 'Histogram t0_bin : ', NXM_histogram_t0_bin

```

```

if (debug) print *, 'Histogram first_good_bin : ', NXM_histogram_first_good_bin
if (debug) print *, 'Histogram last_good_bin : ', NXM_histogram_last_good_bin
if (debug) print *, 'Histogram units : ', NXM_histogram_counts_units

k=0
do i=1, NXM_histogram_number
  do j=1, NXM_histogram_length
    k=k+1
    NXM_histogram_counts (i,j) = temp_data (k)
  end do
end do

C read histogram_resolution
if (NXopendata (file_id, 'resolution') .NE. NX_OK) stop
if (NXgetdata (file_id, NXM_histogram_resolution) .NE. NX_OK) stop
ATTLEN=12
datatype = NX_CHAR
if (NXgetcharattr (file_id, 'units', NXM_histogram_resolution_units, ATTLEN, datatype) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Histogram_resolution : ', NXM_histogram_resolution
if (debug) print *, 'Histogram_resolution units : ', NXM_histogram_resolution_units

C read histogram time zero
if (NXopendata (file_id, 'time_zero') .NE. NX_OK) stop
if (debug) print *, 'Reading histogram time_zero...'
ATTLEN=1
datatype = NX_INT32
if (NXgetattr (file_id, 'available', NXM_histogram_t0_available, ATTLEN, datatype) .NE. NX_OK) stop
if (NXM_histogram_t0_available.GT.0) then

  if (NXgetdata (file_id, NXM_histogram_time_zero) .NE. NX_OK) stop
  if (debug) print *, "Histogram time_zero '", NXM_histogram_time_zero, "' value found"
  ATTLEN=12
  datatype = NX_CHAR
  if (NXgetcharattr (file_id, 'units', NXM_histogram_time_zero_units, ATTLEN, datatype) .NE. NX_OK) stop
  if (debug) print *, "Histogram time_zero units : ", NXM_histogram_time_zero_units
else
  if (debug) print *, 'Histogram time_zero not available'
end if
if (NXclosedata (file_id) .NE. NX_OK) stop

C read raw time
if (NXopendata (file_id, 'raw_time') .NE. NX_OK) stop
if (debug) print *, 'Reading raw time...'
if (NXgetdata (file_id, NXM_histogram_raw_time) .NE. NX_OK) stop
ATTLEN=14
datatype = NX_CHAR
if (NXgetcharattr (file_id, 'units', NXM_histogram_raw_time_units, ATTLEN, datatype) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Histogram raw time units : ', NXM_histogram_raw_time_units

C read corrected time
if (NXopendata (file_id, 'corrected_time') .NE. NX_OK) stop
if (debug) print *, 'Reading corrected time...'
if (NXgetdata (file_id, NXM_histogram_corrected_time) .NE. NX_OK) stop
ATTLEN=14
datatype = NX_CHAR
if (NXgetcharattr (file_id, 'units', NXM_histogram_ctime_units, ATTLEN, datatype) .NE. NX_OK) stop
if (NXclosedata (file_id) .NE. NX_OK) stop
if (debug) print *, 'Histogram corrected time units : ', NXM_histogram_ctime_units

C read grouping
if (NXopendata (file_id, 'grouping') .NE. NX_OK) stop
if (debug) print *, 'Reading histogram groupings...'
ATTLEN=1
datatype = NX_INT32
if (NXgetattr (file_id, 'available', NXM_histogram_grp_available, ATTLEN, datatype) .NE. NX_OK) stop
if (NXM_histogram_grp_available.GT.0) then
  if (debug) print *, 'NXM_histogram_grp_available, groups found'
  if (NXgetdata (file_id, NXM_histogram_grouping) .NE. NX_OK) stop
else
  if (debug) print *, 'No groupings available'
end if
if (NXclosedata (file_id) .NE. NX_OK) stop

C read alpha
if (NXopendata (file_id, 'alpha') .NE. NX_OK) stop
if (debug) print *, 'Reading alpha for grouped pairs...'
ATTLEN=1
datatype = NX_INT32
if (NXgetattr (file_id, 'available', NXM_histogram_alpha_available, ATTLEN, datatype) .NE. NX_OK) stop
if (NXM_histogram_alpha_available.GT.0) then

```

```

        if (debug) print *, NXM_histogram_alpha_available, ' alpha pairs found'
        if (NXgetdata (file_id, NXM_histogram_alpha) .NE. NX_OK) stop
    else
        if (debug) print *, 'No alpha pairs available'
    end if
    if (NXclosedata (file_id) .NE. NX_OK) stop

C      @@@@@@@@@@@@@ close NXdata group @@@@@@@@@@@@@@@@@@
    if (NXclosegroup (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Closing NXdata...'

C      @@@@@@@@@@@@@ open NXlog group @@@@@@@@@@@@@@@@@@
    if (NXopengroup (file_id, 'temperature_log_1', 'NXlog') .NE. NX_OK) stop
    if (debug) print *, ' '
    if (debug) print *, 'Opening NXlog...'

C      read name and available attribute
    if (NXopendata (file_id, 'name') .NE. NX_OK) stop
    if (NXgetchardata (file_id, NXM_temp_log_name) .NE. NX_OK) stop
    ATTLEN=1
    datatype = NX_INT32
    if (NXgetattr (file_id, 'available', NXM_temp_log_available, ATTLEN, datatype) .NE. NX_OK) stop
    if (NXclosedata (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Logged name : ', NXM_temp_log_name
    if (debug) print *, NXM_temp_log_available, ' values available'

    if (NXM_temp_log_available.gt.0) then

C      read temperature values and units
    if (NXopendata (file_id, 'values') .NE. NX_OK) stop
    if (debug) print *, 'reading temperature values...'
    if (NXgetdata (file_id, NXM_temp_log_values) .NE. NX_OK) stop
    ATTLEN=14
    datatype = NX_CHAR
    if (NXgetcharattr (file_id, 'units', NXM_temp_log_values_units, ATTLEN, datatype) .NE. NX_OK) stop
    if (debug) print *, 'First temperature value logged : ', NXM_temp_log_values (1)

C      read temperature time and units
    if (NXopendata (file_id, 'time') .NE. NX_OK) stop
    if (debug) print *, 'reading time values...'
    if (NXgetdata (file_id, NXM_temp_log_time) .NE. NX_OK) stop
    ATTLEN=14
    datatype = NX_CHAR
    if (NXgetcharattr (file_id, 'units', NXM_temp_log_time_units, ATTLEN, datatype) .NE. NX_OK) stop
    if (debug) print *, 'Logging started at time : ', NXM_temp_log_time (1)

    end if

C      @@@@@@@@@@@@@ close NXlog group @@@@@@@@@@@@@@@@@@
    if (NXclosegroup (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Closing NXlog...'

C      @@@@@@@@@@@@@ open NXlog group @@@@@@@@@@@@@@@@@@
    if (NXopengroup (file_id, 'event_log_1', 'NXlog') .NE. NX_OK) stop
    if (debug) print *, ' '
    if (debug) print *, 'Opening NXlog...'

C      read name and available attribute
    if (NXopendata (file_id, 'name') .NE. NX_OK) stop
    if (NXgetchardata (file_id, NXM_event_log_name) .NE. NX_OK) stop
    ATTLEN=1
    datatype = NX_INT32
    if (NXgetattr (file_id, 'available', NXM_event_log_available, ATTLEN, datatype) .NE. NX_OK) stop
    if (NXclosedata (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Logged name : ', NXM_event_log_name
    if (debug) print *, NXM_event_log_available, ' values available'

    if (NXM_event_log_available.gt.0) then

C      read event values and units
    if (NXopendata (file_id, 'values') .NE. NX_OK) stop
    if (debug) print *, 'reading event values...'
    if (NXgetdata (file_id, NXM_event_log_values) .NE. NX_OK) stop
    ATTLEN=14
    datatype = NX_CHAR
    if (NXgetcharattr (file_id, 'units', NXM_event_log_values_units, ATTLEN, datatype) .NE. NX_OK) stop
    if (debug) print *, 'First event value logged : ', NXM_event_log_values (1)

C      read event time and units
    if (NXopendata (file_id, 'time') .NE. NX_OK) stop
    if (debug) print *, 'reading time values...'
    if (NXgetdata (file_id, NXM_event_log_time) .NE. NX_OK) stop
    ATTLEN=14

```

```

        datatype = NX_CHAR
        if (NXgetcharattr (file_id, 'units', NXM_event_log_time_units, ATTLEN, datatype) .NE. NX_OK) stop
        if (debug) print *, 'Logging started at time : ', NXM_event_log_time (1)

    end if

C      @@@@@@@@@@@@@@@@@@ close NXlog group @@@@@@@@@@@@@@@@@@
    if (NXclosegroup (file_id) .NE. NX_OK) stop
    if (debug) print *, 'Closing NXlog...'

C      @@@@@@@@@@@@@@@@@@ close NXentry group @@@@@@@@@@@@@@@@@@
    if (NXclosegroup (file_id) .NE. NX_OK) stop
    if (debug) print *, ' '
    if (debug) print *, 'Closing NXentry...'

C      @@@@@@@@@@@@@@@@@@ close NeXus file @@@@@@@@@@@@@@@@@@
    if (NXclose (file_id) .NE. NX_OK) stop
    if (debug) print *, ' '
    if (debug) print *, 'Closing NeXus file...'

NXMread = 1
END

```

APPENDIX G: TEST_NEXUS_READER.FOR

```
program test_nexus_reader
  include 'muon_def.inc'
  character*60 fname
  integer status, i, j, k, histogram

  print *, ' '
  print *, ' '
  print *, 'Welcome, this is a test program to read ISIS Muon NeXus files'
  print *, 'If you have any questions or comments please contact.....'
  print *, 'Damian.Flannery@rl.ac.uk      Thank You.'
  print *, ' '
  print *, 'NB, When entering NeXus filename, please input complete '
  print *, 'filename in single quotes e.g. ''30000.nxs'' '
  print *, ' '

  WRITE (*,*) 'Please Enter Nexus File Name'
  READ (*,*) fname
  print *, ' '
  print *, ' '
  status = NXMread (fname)

  print *, 'NeXus file Contents...'
  print *, ' '

C   NXrun -----
  print *, 'Program name : ',NXM_run_program_name
  print *, 'Program version : ',NXM_run_program_version
  print *, 'Run number : ', NXM_run_number
  print *, 'Title : ', NXM_run_title
  print *, 'Notes : ', NXM_run_notes
  print *, 'Start time : ', NXM_run_start_time
  print *, 'Stop time : ', NXM_run_stop_time
  print *, 'Duration : ', NXM_run_duration
  print *, 'Switching states : ', NXM_run_switching_states

C   NXuser-----
  print *, 'User name : ', NXM_user_name
  print *, 'Experiment number : ', NXM_user_experiment_number

C   print user information
  do i=1, NXM_user_uif_array_length
    print *, 'user information : ', i
    print *, '  Data name : ', NXM_user_uif_element_name (i)
    print *, '  Data value : ', NXM_user_uif_element_value (i)
  end do

C   NXsample-----
  print *, 'Sample : ', NXM_sample_name
  print *, 'Temperature : ', NXM_sample_temperature,
&         ' ',NXM_sample_temperature_units
  print *, 'Magnetic field : ', NXM_sample_magnetic_field,
&         ' ',NXM_sample_magnetic_field_units
  print *, 'Shape : ', NXM_sample_shape
  print *, 'Magnetic field state : ', NXM_sample_magnetic_field_state
  if (NXM_sample_mfield_vec_available .NE. 0) then
    print *, 'Field vector : {', NXM_sample_mfield_vector(1),
&         ' ', ' ', NXM_sample_mfield_vector(2),
&         ' ', ' ', NXM_sample_mfield_vector(3),'}'
    print *, 'Coordinate system : ',NXM_sample_mfield_vector_coord
  else
    print *, 'Magnetic field vector not available'
  end if
  print *, 'Sample environment : ', NXM_sample_environment

C   NXinstrument-----
  print *, 'Instrument name : ', NXM_instrument_name
```

```

C   NXdetector-----
print *, 'Number of detectors : ', NXM_detector_number
print *, 'Detector orientation : ', NXM_detector_orientation

if (NXM_detector_angles_available.NE.0) then
  print *, 'Position of first detector : [',NXM_detector_angles(1,1),
&      ', ',NXM_detector_angles(2,1),
&      ', ',NXM_detector_angles(3,1), ']'
  print *, 'Coordinate System : ', NXM_detector_angles_coord
  print *, 'Solid angle of first detector : ', NXM_detector_angles(4,1)
else
  print *, 'Detector angles not available'
end if

if (NXM_deadtimes_available.NE.0) then
  print *, 'Deadtimes available in ', NXM_detector_deadtimes_units
  do i=1, NXM_detector_number
    print *, 'Deadtime of detector ', i, ', ', NXM_detector_deadtimes (i)
  end do
else
  print *, 'Deadtime information not available'
end if

C   NXcollimator-----
print *, 'Collimator type : ', NXM_collimator_type
print *, 'Collimator aperture : ', NXM_collimator_aperture

C   NXbeam-----
print *, 'Events : ', NXM_beam_total_counts
print *, 'Frames : ', NXM_beam_daereads
print *, 'DAEreads : ', NXM_beam_frames

C   NXdata-----
print *, 'Number of histograms : ', NXM_histogram_number
print *, 'Time histogram length : ', NXM_histogram_length
print *, 'Time zero bin : ', NXM_histogram_t0_bin
print *, 'First good bin : ', NXM_histogram_first_good_bin
print *, 'Last good bin : ', NXM_histogram_last_good_bin
print *, 'Histogram_resolution : ', NXM_histogram_resolution,
&      ' in ',NXM_histogram_resolution_units
print *, 'Histogram units : ', NXM_histogram_counts_units

C   print grouping information
if (NXM_histogram_grp_available.gt.0) then
  print *, 'Found ',NXM_histogram_grp_available,' histogram groups'
  do i=1,NXM_histogram_grp_available
    print *, 'group ', i, 'histograms : '
    do j=1, NXM_histogram_number
      if (NXM_histogram_grouping (j) .eq. i) print *, j
    end do
  end do
end if

C   print time_zero information
if (NXM_histogram_t0_available.gt.0) then
  print *, 'Histogram time_zero : ', NXM_histogram_time_zero
  print *, 'Histogram time_zero units : ', NXM_histogram_time_zero_units
end if

C   print individual histogram elements, choose histogram to manipulate (print every 50th value)
histogram = 1
print *, 'Histogram information - (raw_time, corrected_time, counts) (time in microseconds)'
do i=50, NXM_histogram_length, 50
  print *, '(',NXM_histogram_raw_time(i),
&      ', ',NXM_histogram_corrected_time(i),', ',NXM_histogram_counts(i, histogram),')'
end do

C   print alpha pairs
if (NXM_histogram_alpha_available.gt.0) then
  print *, 'Alpha defined for ',NXM_histogram_alpha_available,' pairs'
  do i=1, NXM_histogram_alpha_available
    print *, 'alpha for group pair ', i,
&      ' (groups ',NXM_histogram_alpha(1,i), ' and ',
&      NXM_histogram_alpha(2,i), ' : ',
&      NXM_histogram_alpha(3,i)
    end do
  end if

```

```

C    NXlog-----
    if (NXM_temp_log_available.gt.1) then
        print *, ' '
        print *, 'Temperature Log - time (seconds), temperature (kelvin)'
        do i=1, NXM_temp_log_available
            print *, '(', NXM_temp_log_time(i), ',', NXM_temp_log_values(i), ')'
        end do
    end if

C    NXlog-----
    if (NXM_event_log_available.gt.0) then
        print *, ' '
        print *, 'Event Log - time (seconds), events (counts)'
        do i=1, NXM_event_log_available
            print *, '(', NXM_event_log_time(i), ',', NXM_event_log_values(i), ')'
        end do
    end if

end

```